# **OPERATING MANUAL**

**Cole-Parmer®**

**GP Laboratory Mixer Set**

Model Number 55008-00

**GP Laboratory Mixer System**

Model Number 55008-01

## **ENGLISH**

Instruction Manual  
pages 3 - 15

## **FRANÇAIS**

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## **ESPAÑOL**

Manual de instrucciones  
paginas 35 - 48

## **ACCESSORY TABLE**

pages 50 - 51

# CONTENTS

|                                      |    |
|--------------------------------------|----|
| SAFETY PRECAUTIONS                   | 5  |
| GENERAL INSTRUCTIONS                 | 7  |
| INTRODUCTION                         | 7  |
| Packing List                         | 8  |
| Cord Sets                            | 8  |
| CONTROLLER AND LED FUNCTIONS         | 8  |
| INSTALLATION AND SETUP PROCEDURES    | 9  |
| Installation of Support Rod          | 9  |
| Installation Using a Mixer Stand     | 10 |
| Impeller Attachment                  | 10 |
| Chuck Removal and Attachment         | 10 |
| MIXER OPERATION                      | 11 |
| Powering On                          | 11 |
| SPEED Control                        | 11 |
| Overheat Indicator                   | 11 |
| Overload Indicator                   | 11 |
| SPECIFICATIONS                       | 12 |
| Torque Capabilities                  | 13 |
| Torque Curve                         | 13 |
| MAINTENANCE AND TROUBLESHOOTING      | 14 |
| Cleaning the Mixer/Controller System | 14 |
| Replacement Parts                    | 14 |
| Troubleshooting                      | 14 |
| WARRANTY                             | 15 |
| PRODUCT RETURN                       | 15 |
| TECHNICAL ASSISTANCE                 | 15 |

# SAFETY PRECAUTIONS

1. All operators must be familiar with the mixer and read this entire manual. If the mixer is used with accessories not provided or recommended by the manufacturer or used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.
2. The mixer must be securely fixed to a stable support. Mount to a Cole-Parmer stand using a Cole-Parmer double (rod) clamp. If other stands or clamps are used, the mixer must be secured such that it will not fall if tipped.
3. **CAUTION:** This is not an explosion proof mixer. Do not use with highly flammable or explosive materials. Do not operate the mixer in explosive atmospheres, with hazardous substances or under water.
4. Spinning impellers can cause severe personal injuries. Operators must use extreme care and good judgment when mixing at any speed.
5. All mixing paddles and impellers must be in good condition with straight shafts. If the mixer vibrates at high speeds check the paddle shaft for damage and repair or replace if needed.
6. Extreme care must be taken when mixing chemicals to ensure that no chemicals are splashed outside the mixing vessel. Care must be taken when changing to faster mixing speeds.
7. Ensure that the mixing impeller does not contact the containment vessel.
8. Do not operate while standing in water. Keep the mixer dry and do not immerse any part, except the mixing paddle, into any liquids. Protect from splashing.
9. Ensure that no loose clothing, jewelry, or hair can become entangled in any rotating parts. A fast spinning chuck can cause injury to operator.
10. To stop the mixer in an emergency, turn the SPEED control knob counter-clockwise to the "0" or Off position or disconnect the mains cord.
11. Ensure the SPEED control knob and mains power cord of the mixer can be accessed immediately, directly and without risk at any time, and above all, in case of fault.
12. The mixer must be used with the supplied power cord. The cord and grounded plug type are important safety elements. Always connect the mixer to a grounded outlet. Do not alter the plug.
13. The mixer must be supplied with the voltage specified on the serial label. See serial label.
14. Excluding procedures outlined under maintenance and troubleshooting, service must be done by an authorized service centre.
15. Be sure SPEED control knob is turned fully counter clockwise to the "0" or Off position before connecting power to the mixer.
16. Wear appropriate personal protective equipment such as eye protection and any clothing suitable for the application.
17. When running the mixer at full speed and load at higher ambient temperatures, the mixer motor may get very hot to the touch.

# GENERAL INSTRUCTIONS

**Read SAFETY PRECAUTIONS.** Unpack the product carefully. Ensure the product is not damaged and verify the contents against the packing list. If the product is damaged or the contents do not correspond with the packing list, notify your dealer immediately.

Read the manual carefully! Make sure that all users read this manual thoroughly before operating the mixer. For further questions or inquiries, please contact:

## **Cole-Parmer**

625 East Bunker Court  
Vernon Hills, IL 60061 United States

Tel: 847-549-7600

Toll Free: 1-800-323-4340

Website: [www.coleparmer.com](http://www.coleparmer.com)

E-mail: [sales@antylia.com](mailto:sales@antylia.com)

# INTRODUCTION

**GP Laboratory Mixer** covered in this manual is specifically designed for professional use in the laboratory and for use in industry to stir numerous liquids with various properties and applications.

**GP Laboratory Mixer** rotates at 40-4000 rpm with a maximum torque of 9.5 in-oz (6.7 Ncm).

The mixer uses a brushed DC motor. The speed is kept constant with continuous feedback to the motor. The mixer has integrated overload protection which ensures no damage can occur to the electronics or motor when an overload situation occurs, such as too great a torque demand or mixer overheating.

The mixer is designed for reliability and durability. There are no wear parts that require user service during the normal life of the mixer.

## Packing List

Compare the contents of the shipment to the items on the following lists to ensure all parts are received. Do not discard the packing materials until all parts are accounted for. It is recommended that the original packing materials be kept for transportation and storage.

**GP Laboratory Mixer Set Model 55008-00** covered in this manual consists of:

- Mixer motor with integrated controller and 6 ft (183 cm) cable
- Power supply (universal voltage 100-240 volt)
- North American (NA) cord set
- Stainless steel shaft for propeller, 12 in x 3/8 in (305 mm x 10 mm)
- Stainless steel propeller, 1.5 in (38 mm)
- Mounting rod including hex drive screw and lock washer, 6 in (152 mm)

**GP Laboratory Mixer System 55008-01** consists of the package above plus:

- Double (rod) clamp
- Mixer stand

## Cord Sets

North American cord set included. RoHS compliant. For international cord sets contact your dealer.

## CONTROLLER AND LED FUNCTIONS

The SPEED Control knob, Figure 1, is used to turn on and set the speed of the mixer motor and provides a 100:1 adjustment range.


Figure 1: Mixer Controller

The variable-speed controller ensures accurate, dependable motor speed (mixing) control, with “soft-start” and smooth, gradual adjustment to the speed setting. Constant speed with a 100:1 adjustment range is maintained electronically via automatic compensation for changes in load torque demands at both high and low speeds. Built-in motor current compensation minimizes speed drift with changes in load. Motor current limiting protects the system in case of a locked rotor condition. Output circuit protection protects against inadvertent output shorting.

LED indicates mixer is powered on when it is green. When LED flashes amber or red, it indicates a limit condition has been reached.

# INSTALLATION AND SETUP PROCEDURES

## Installation of Support Rod

The mixer is equipped with a 6 in (152 mm) mounting rod. The mounting assembly is designed to be attached to an optional mixer support stand.

Attach mixer motor to support rod, as shown in Figure 2, using hex drive screw and lock washer provided.

1. Remove the hex screw and lock washer from the support rod.
2. Insert the support rod into the back of the mixer so that the holes line up for the hex screw.
3. Secure with hex screw and lock washer and tighten with hex key.
4. Clamp the mixer to a recommended stand on a stable, level surface (see SAFETY PRECAUTIONS) or tank clamp (not supplied). Connect mixer to supplied power supply. Connect power supply to a grounded electrical outlet.


Figure 2: Mixer with Mounting Rod Attached

Recommended stands, tank clamps, and a double (rod) clamp are available through your dealer. Part numbers are listed in the ACCESSORIES section.

## Installation Using a Mixer Stand

1. Set up mixer stand in accordance with the manufacturer's instructions.
2. Attach a double (rod) clamp to the vertical rod on the mixer stand and tighten clamp screw.
3. Slide rod of mixer into clamp and secure with clamp screw.
4. Place mixer stand legs around vessel.
5. Position the height of the mixer motor by moving the clamp up or down on the mixer stand.


Figure 3: Typical Setup Using a Mixer Stand

Note: The power supply can be placed on any flat surface within 6 ft (183 cm) of the mixer.

6. Connect the 6 ft (183 cm) cable between the mixer motor and the power supply.
7. Check that the controller is in the "0" or Off position.
8. Connect the 6 ft (183 cm) line cord from the power supply to the appropriate 3-wire, grounded power source.

## Impeller Attachment

Select an impeller and shaft for the required operation and install as follows:

1. Place impeller on shaft and tighten set screw.
2. Open the chuck. Move the impeller shaft upwards into the open chuck until shaft is fully inserted.
3. Tighten the chuck.

## Chuck Removal and Attachment

1. Slide chuck onto motor shaft.
2. Tighten the set screw using the provided hex key.
3. To remove, loosen set screw and remove chuck from motor shaft.

# MIXER OPERATION

The mixer can operate continuously if the torque demand falls below the continuous torque operation line (see Torque Capabilities).

If the mixer is used for an extended period of time above the continuous torque operation line (see Torque Capabilities), then it may overheat and shut down. The mixer may be started after it has cooled down.

## Powering On

Make sure that the DC power cable at the motor is connected to the power supply and that the power supply is connected to the mains supply. A green LED illuminates on the power supply indicating it is energized with power.

SPEED Control in “0” position is Off. When SPEED Control is turned clockwise just past zero there is a click and the green LED illuminates indicating power is on.

## SPEED Control

The SPEED Control precisely sets the output speed of the shaft. The controller has a range of 100:1 depending on the shaft; speed can be as low as 40 rpm or as high as 4000 rpm.

Slowly adjust the SPEED Control from the “0” or Off position to a setting that provides the desired mixing speed. The controller will ramp up the mixing speed to match the setting of the control.

The table below outlines the typical speeds for each speed setting.

| Setting   | 3   | 4   | 5    | 6    | 7    | 8    | 9    | 10   |
|-----------|-----|-----|------|------|------|------|------|------|
| Speed rpm | 250 | 900 | 1600 | 2100 | 2700 | 3200 | 3800 | 4000 |

Always return the SPEED Control to the “0” or Off position when mixing is complete.

**Note:** Allow a 20 minute run-in time to warm-up mixer so speed drift is minimal.

If the SPEED Control knob is left On at one of the speed settings and remains plugged into AC power, the mixer will resume mixing after a power failure or when the power has been removed by control equipment (i.e., PLC) or when the AC power has been cut (i.e., the mixer unplugged). When power is restored, the mixer will automatically ramp up to the last speed. This function is inherent with the analog signal. It cannot be enabled or disabled by the end user; it is always enabled.

## Overheat Indicator

The green LED indicator will change to a flashing amber when the mixer control senses an over-temperature situation. An over-temperature fault is caused by long running time at high loads.

Allow the mixer to cool down.

The mixer will remain in a stopped position until it is turned off and back on again after it has time to cool down.

## Overload Indicator

The green LED indicator will change to a flashing red when the mixer control senses an over-current situation. Over-current faults can be caused by too much load (torque) on the mixer or by a locked motor condition.

Resolve the condition causing too much load or a locked motor.

The mixer will remain in a stopped position until it is turned off and back on again.

# SPECIFICATIONS

|   |   |
|---|---|
| <b>Direction of Rotation</b>              | Clockwise when viewed from mixer down the shaft   |
| <b>Speed Range</b>                        | 40 to 4000 rpm  |
| <b>Torque Output max</b>                  | 7.0 in-oz (4.9 Ncm) Continuous<br>9.5 in-oz (6.7 Ncm) Intermittent  |
| <b>Maximum Volume</b>                     | 10 gal (38 L) [1 cP]  |
| <b>Maximum Viscosity</b>                  | 100 cP [ $< 1$ gal (4 L)]   |
| <b>Keyless Chuck</b>                      | Three-jawed type chuck accepts up to 3/8 in (10 mm) diameter shafts   |
| <b>Motor Type</b>                         | Brushed DC Motor  |
| <b>Voltage</b>                            | 100-240 VAC, 50/60 Hz<br>+/- 10% of the nominal voltage   |
| <b>Maximum Current Draw</b>               | 2.0 Amps  |
| <b>Dimensions</b>                         |   |
| <b>Power Supply Cord Length</b>           | 6 ft (183 cm) from mixer to power supply - DC<br>6 ft (183 cm) from power supply to mains - AC  |
| <b>Weight</b>                             | 4 lb (1.8 kg)   |
| <b>Enclosure Rating (IP)</b>              | IP44 per IEC 60529  |
| <b>Environmental Operating Conditions</b> | Indoors <ol style="list-style-type: none"> <li>Altitudes up to 6500 ft (2000 m)</li> <li>Temperatures from 41°F to 104°F (5°C to 40°C)</li> <li>Maximum relative humidity 80% for temperatures up to 88°F (31°C) decreasing linearly to 50% relative humidity at 104°F (40°C)</li> <li>Pollution Degree 2 in accordance with IEC 664</li> </ol> |
| <b>Chemical Resistance</b>                | Exposed material is painted or anodized aluminum, plastic and vinyl   |
| <b>Compliance</b>                         | CAN/CSA C22.2 No. 61010-1:2019<br>CAN/CSA-C22.2 No. 61010-2-051:2015<br>UL 61010-1:2019<br>UL 61010-2-051:2015<br>IEC 61010-1:2010, IEC 61010-2-051:2015<br>EN 61010-1:2010, EN 61010-2-051:2015<br>FCC Part 15 Subpart B, Class A<br>ICES-003:2020<br>IEC 61326-1:2012<br>CISPR 11:2017<br>EN 61326-1:2013, EN 55011:2016                      |


- RoHS Restriction of Hazardous Substances (EU)
- WEEE Waste Electrical and Electronic Equipment (EU)
- CSA Canadian Standards Association
- CE Communauté Européenne (European Community)
- UKCA UK Conformity Assessed
- UL Underwriters Laboratories
- FCC Federal Communications Commission
- IEC The International Electrotechnical Commission
- ICES Interference Causing Equipment Standards
- EMC Electromagnetic Compatibility
- EMI Electromagnetic Interference

## Torque Capabilities

|                | Continuous          | Intermittent        |
|----------------|---------------------|---------------------|
| Model 55008-00 | 7.0 in-oz (4.9 Ncm) | 9.5 in-oz (6.7 Ncm) |

| Torque Conversion Table |                   |                           |
|-------------------------|-------------------|---------------------------|
| in-lb to in-oz          | Multiply by 16    | A in-lb = B in-oz = C Ncm |
| in-lb to Ncm            | Multiply by 11.36 |                           |
| Ncm to in-oz            | Multiply by 1.416 |                           |

## Torque Curve

# MAINTENANCE AND TROUBLESHOOTING

Maintenance is limited to cleaning the mixer. If a fault is detected when troubleshooting and cannot be corrected by the suggested remedy, contact your dealer for service.

## Cleaning the Mixer/Controller System

Keep the mixer housing and controller enclosure clean using mild detergents. Never immerse or use excessive fluid.

## Replacement Parts

There are no user-serviceable parts for this mixer. Contact your dealer for service.

The following parts are user-replaceable:

| Description                               | Part Number |
|---|-------------|
| Replacement Chuck - Cole-Parmer GP        | 04552-90    |
| Replacement Power Supply - Cole-Parmer GP | 55008-12    |

## Troubleshooting

| Symptom   | Cause  | Remedy  |
|---|--|---|
| <p>Mixer shaft does not rotate when SPEED Control is turned past "0" and click is experienced.</p> <p>Mixer LED POWER Indicator is not illuminated.</p> <p>Power supply LED POWER Indicator is not illuminate.</p>      | <ul style="list-style-type: none"> <li>• Disconnected AC line cord.</li> <li>• Defective line cord and/or power supply.</li> </ul> | <ol style="list-style-type: none"> <li>1. Check mixer is plugged into a live line.</li> <li>2. Check cord connections to the mains supply and the power supply to the mixer.</li> <li>3. Replace power supply. Contact dealer for service.</li> </ol> |
| <p>Mixer shaft does not rotate when SPEED Control is turned past "0" and click is experienced.</p> <p>Mixer LED POWER Indicator is not illuminated.</p> <p>Power supply LED POWER Indicator is illuminated green.</p>   | <ul style="list-style-type: none"> <li>• Defective DC connection.</li> <li>• Defective mixer motor or controller.</li> </ul>       | <ol style="list-style-type: none"> <li>1. Check cord connections from the power supply to the mixer.</li> <li>2. Contact dealer for service.</li> </ol>   |
| <p>Mixer shaft does not rotate when SPEED Control is turned past "0" and click is experienced.</p> <p>Mixer LED POWER Indicator is illuminated green.</p> <p>Power supply LED POWER Indicator is illuminated green.</p> | <ul style="list-style-type: none"> <li>• Defective mixer motor or controller.</li> </ul>   | <ol style="list-style-type: none"> <li>1. Contact dealer for service.</li> </ol>  |

|   |  |  |
|---|--|--|
| <p>Mixer shaft does not rotate when SPEED Control is set past "0".</p> <p>Mixer LED POWER Indicator flashed amber before mixer stopped.</p> | <ul style="list-style-type: none"> <li>Mixer overheat caused by long running time at high loads.</li> </ul>                                    | <ol style="list-style-type: none"> <li>Turn mixer off.</li> <li>Wait approximately 30 minutes to allow the mixer to cool down.</li> <li>Turn mixer on again to resume mixing.</li> </ol> |
| <p>Mixer shaft does not rotate when SPEED Control is turned past "0"</p> <p>Mixer LED POWER Indicator flashed red before mixer stopped.</p> | <ul style="list-style-type: none"> <li>Mixer overload caused by too much load (torque) on the mixer or by a locked motor condition.</li> </ul> | <ol style="list-style-type: none"> <li>Turn mixer off.</li> <li>Resolve the condition causing too much load or a locked motor.</li> <li>Turn mixer on again to resume mixing.</li> </ol> |

## WARRANTY

The manufacturer warrants this product to be free from significant deviations from published specifications. If repair or adjustment is necessary within the 2 year warranty period, the problem will be corrected at no charge if it is not due to misuse or abuse as determined by the manufacturer. Repair costs outside the warranty period, or those resulting from product misuse or abuse, may be invoiced. To register the warranty online, go to the Cole-Parmer website address below:

[www.coleparmer.com/product-warranty-registration](http://www.coleparmer.com/product-warranty-registration)

## PRODUCT RETURN

To limit charges and delays, contact your dealer for authorization and shipping instructions before returning the product, either within or outside of the warranty period. Clean the motor and controller so that it is free from any materials which may be harmful. Provide an SDS or other safety data sheets for all chemicals used in conjunction with this system where appropriate. When returning the product, please state the reason for the return. For your protection, pack the product carefully, add labels and follow shipping directions provided. Insure it against possible damage or loss. Any damages resulting from improper packaging are your responsibility.

## TECHNICAL ASSISTANCE

If you have any questions about the use of this product, contact your local dealer.

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**Cole-Parmer®**

## ACCESSORIES

The accessories listed below, used for mounting the mixer head and controller, are available from your dealer. Many other accessories such as shafts and impellers are also available. Contact your dealer for the full line of available accessories.

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|   | Part Number | Description   |
|---|-------------|---|
|    | 55008-06    | Cole-Parmer SP Controller<br>91 in-oz (64 Ncm) max torque<br>25 gal (95 L) [250 cP]<br>1000 cP [ $<$ 1 gal (4 L)]       |
|     | 55008-07    | Cole-Parmer SV Controller<br>113 in-oz (80 Ncm) max torque<br>31 gal (120 L) [250 cP]<br>2000 cP [ $<$ 1 gal (4 L)]     |
|    | 55008-08    | Cole-Parmer SP/SV Mixer Motor<br>50-2500 rpm<br>Keyless chuck<br>6 in (152 mm) mounting rod with tilt bracket           |
|     | 08041-22    | Double (rod) Clamp  |
|     | 04552-70    | Pivot Mount for Controller, includes clamp  |
|   | 04552-80    | Mixer Stand, 19-1/4 x 11 in (489 x 279 mm) base with 5/8 in dia x 28 in (16 x 711 mm) stainless steel rod               |
|  | 04552-65    | Tank Clamp  |
|  | 04552-25    | Shaft, Stainless Steel, 3/8 in dia x 12 in (10 x 305 mm)<br>• Accepts 5/16 in or 3/8 in (8 mm or 10 mm) bore propellers |
|   | 04552-35    | Shaft, Stainless Steel 3/8 in dia x 18 in (10 x 457 mm)<br>• Accepts 5/16 in or 3/8 in (8 mm or 10 mm) bore propellers  |
|   | 04552-00    | Shaft, Stainless Steel, 3/8 in dia x 24 in (10 x 610 mm)<br>• Accepts 5/16 in or 3/8 in (8 mm or 10 mm) bore propellers |
|   | 04552-05    | Shaft, Stainless Steel, 3/8 in dia x 30 in (10 x 762 mm)<br>• Accepts 5/16 in or 3/8 in (8 mm or 10 mm) bore propellers |

|   |          |   |
|---|----------|---|
|     | 04552-20 | Shaft, Stainless Steel 5/16 in dia x 12 in (8 x 305 mm)                               |
|   | 04552-30 | Shaft, Stainless Steel, 5/16 in dia x 18 in (8 x 457 mm)                              |
|    | 04552-40 | Three-blade Propeller, Stainless Steel, 1 in dia x 5/16 in bore dia (25 x 8 mm)       |
|   | 04552-45 | Three-blade Propeller, Stainless Steel, 1.5 in dia x 5/16 in bore dia (38 x 8 mm)     |
|   | 04552-50 | Three-blade Propeller, Stainless Steel, 2 in dia x 5/16 in bore dia (50 x 8 mm)       |
|   | 04552-60 | Three-blade Propeller, Stainless Steel, 3 in dia x 5/16 in bore dia (78 x 8 mm)       |
|    | 04552-55 | Three-blade Propeller, Stainless Steel, 2 in dia x 3/8 in bore dia (50 x 10 mm)       |
|    | 04552-10 | Turbine Propeller, Stainless Steel, 2 in dia x 5/16 in bore dia (50 x 8 mm)           |
|   | 04552-15 | Turbine Propeller, Stainless Steel, 2 in dia x 3/8 in bore dia (50 x 10 mm)           |
|    | 04552-61 | Straight Blade impeller only 1-9/16 in dia x 5/16 in bore dia (40 x 8 mm)             |
|   | 04552-62 | Dispersion Blade impeller only 1-7/8 in dia x 5/16 in bore dia (48 x 8 mm)            |
|  | 50801-53 | Square Blade impeller; 2-3/4 in (70 mm) dia with 5/16 x 16 in shaft (8 x 406 mm)      |
|  | 50801-56 | Cross Blade impeller; 2 in (50 mm) dia with 5/16 x 16 in shaft (8 x 406 mm)           |
|  | 50801-58 | Collapsible blade impeller; 3-1/2 in (90 mm) dia with 5/16 x 16 in shaft (8 x 406 mm) |
|  | 50801-60 | Straight blade impeller; 2 in (50 mm) dia with 5/16 x 16 in shaft (8 x 406 mm)        |
|  | 50801-62 | Pitched blade impeller; 2-1/2 in (64 mm) dia with 5/16 x 16 in shaft (64 x 406 mm)    |

# Cole-Parmer®

625 E. Bunker Court  
Vernon Hills, IL 60061 United States

847-549-7600

Toll Free: 1-800-323-4340 (USA)

[sales@antylia.com](mailto:sales@antylia.com)

[www.coleparmer.com](http://www.coleparmer.com)

Canada

1-514-355-6100

[info@antylia.ca](mailto:info@antylia.ca)

[www.coleparmer.ca](http://www.coleparmer.ca)

United Kingdom

+44 (0) 1480 272279

[uk.sales@antylia.com](mailto:uk.sales@antylia.com)

[www.coleparmer.co.uk](http://www.coleparmer.co.uk)

Germany

+49 (0) 9377 92030

[de.sales@antylia.com](mailto:de.sales@antylia.com)

[www.coleparmer.de](http://www.coleparmer.de)

France

+33 (0) 1486 37800

[fr.sales@antylia.com](mailto:fr.sales@antylia.com)

[www.coleparmer.fr](http://www.coleparmer.fr)

Italy

+39 (0) 2843 49215

[it.sales@antylia.com](mailto:it.sales@antylia.com)

[www.coleparmer.it](http://www.coleparmer.it)

India

+9122 61394444

[info@coleparmer.in](mailto:info@coleparmer.in)

[www.coleparmer.in](http://www.coleparmer.in)

China

+1 847-549-7600

[sales@antylia.com](mailto:sales@antylia.com)

[www.coleparmer.com](http://www.coleparmer.com)