**LANGE** ©  
UNITED FOR WATER QUALITY

# LCK 514

100–2000 mg/L

T1

NL

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EN

## LCK 514 COD Chemical Oxygen Demand

**Please check the "Edition Date"** (see data table) and read the "Note". Safety advice and expiry date on package.

### Principle

Oxidizable substances react with sulphuric acid – potassium dichromate solution in the presence of silver sulphate as a catalyst. Chloride is masked by mercury sulphate. The green coloration of Cr<sup>3+</sup> is evaluated.

### Range of Application

Waste water, process analysis

### Interferences

The method can be used for samples (or diluted samples) with chloride concentrations of up to 1500 mg/L.

The measurement results must be subjected to plausibility checks (dilute and/or spike the sample).

### Note

In contrast to the classic COD Cuvette Test (COD classic) the HT-COD is characterised by a higher digestion temperature and shorter digestion time.

Users are advised to carry out a comparison with the COD classic, in order to be sure that the results obtained from their own samples when using the HT-COD are comparable to the standard.

### Special note

For optimal stability until its expiry date, it is recommended the reagent LCK 514 is stored in a fridge.

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EN

1. Bring the sediment into suspension by inverting a few times.
2. **Carefully** pipette **2.0 mL** sample.
3. Close cuvette, thoroughly clean the outside.
4. Invert.
5. Heat in the thermostat.  
**HT 200 S:** in standard program **HT** for **15 min**  
**COD classic:** **2 h at 148 °C**
6. Remove the **hot** cuvette.
  - a. **HT 200 S:** After the lock opens, **carefully** invert **twice**.
  - b. **COD classic:** **Carefully invert twice**.
7. Allow to cool to room temperature.
  - a. **HT 200 S:** in the thermostat
  - b. **COD classic:** in a cooling rack
8. **HT 200 S:** Sediment must be completely settled before evaluation is carried out. Clean the outside of the cuvette and evaluate.  
**COD classic:** Clean the outside of the cuvette and evaluate.

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EN: If any of the following barcode instruments is used, an automatic evaluation is carried out after the sample cuvette is inserted:

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