
LANGE 

sensION™ + Portable Meter

02/2013, Edition 3

User Manual
Bedienungsanleitung
Manuale dell'utente
Manuel d'utilisation
Manual del usuario
Manual do utilizador
Uživatelská příručka
Brugsanvisning
Gebruikershandleiding
Instrukcja obsługi
Bruksanvisning
Käyttäjän käsikirja
Ръководство за потребителя
Használati útmutató
Manual de utilizare
Naudotojo vadovas
Руководство пользователя
Kullanım Kılavuzu
Návod na obsluhu
Navodila za uporabo
Korisnički priručnik
Εγχειρίδιο χρήστη
Kasutusjuhend
Priručnik za korisnika

English.....	3
Deutsch.....	18
Italiano.....	35
Français.....	51
Español.....	67
Português.....	83
Čeština.....	99
Dansk.....	114
Nederlands.....	130
Polski.....	147
Svenska.....	163
Suomi.....	178
български.....	193
Magyar.....	211
Română.....	227
lietuvių kalba.....	243
Русский.....	258
Türkçe.....	276
Slovenský jazyk.....	291
Slovenski.....	307
Hrvatski.....	322
Ελληνικά.....	337
eesti keel.....	355
Српски.....	370

Table of contents

- [Specifications](#) on page 3
- [General information](#) on page 4
- [Installation](#) on page 5
- [User interface and navigation](#) on page 6
- [Startup](#) on page 10
- [Operation](#) on page 10
- [Maintenance](#) on page 13
- [Troubleshooting](#) on page 13
- [Replacement parts and accessories](#) on page 15
- [Standard solutions](#) on page 15
- [Salinity values for dissolved oxygen meters](#) on page 17

Specifications

Specifications are subject to change without notice.

Specification	Details
Dimensions	18.6 x 7.3 x 3.8 cm (7.32 x 2.87 x 1.5 in.)
Weight	300 g (0.66 lb)
Meter enclosure	IP67
Power requirements (internal)	AA Alkaline or rechargeable Nickel Metal Hydride (NiMH) batteries (3); battery life: over 500 hours
Input connectors	MP-5 or MP-8 connector for sensION+ probes
Meter protection class	Class III
Storage temperature	-15 to +65 °C (5 to +149 °F)
Operating temperature	0 to 50 °C (32 to 122 °F)
Operating humidity	80% (non-condensing)

Specification	Details
Instrument range	pH: -2.00 to 19.99; ORP: ±1999 mV
	Conductivity: 0.01 to 500 mS/cm; TDS: 0 to 500 g/L; Salinity: 0.0 to 1999 mg/L, 2.0 to 50.0 g/L
	DO: 0.00 to 19.99 mg/L and 20.0 to 22.0 mg/L (25°C), 0.0 to 199.9% and 200 to 250% (25°C)
Resolution	pH: 0.01 pH; ORP: 1 mV (0.1 mV from -199.9 to 199.9 mV)
	DO: 0.1% (1% if ≥ 200%); 0.01 mg/L (0.1 mg/L if ≥ 20 mg/L)
	Conductivity: changes with range (auto-ranging)
	Temperature: 0.1 °C (0.1°F)
Measuring error (± 1 digit)	pH: ≤ 0.01 pH; ORP: ≤ 1 mV
	DO: ≤ 0.5% of the measured value
	Conductivity: ≤ 0.5% (0.01 μS/cm to 19.99 mS/cm); ≤ 1% (≥ 20.0 mS). Salinity/TDS: ≤ 0.5%
	Temperature: ≤ 0.2 °C (≤ 0.4 °F)
Reproducibility (± 1 digit)	pH: ± 0.01 pH; ORP ± 1 mV
	DO: ≤ 0.2 of the measured value
	Conductivity/Salinity/TDS: ± 0.1%
	Temperature ± 0.1 °C (± 0.1 °F)
Reference temperature (RT)	Conductivity: 20 or 25 °C (68 or 77 °F) (factory settings: 25 °C (77 °F))
Temperature coefficient (TC)	Conductivity: 0.00 to 5.00%/°C (factory settings: 2%/°C)
TDS conversion factor	Conductivity: 0.00 to 4.44 (factory settings: 0.64)
Energy management	Automatic power off after 5 minutes of inactivity
Certification	CE

General information

Revised editions are found on the manufacturer's website.

Safety information

NOTICE

The manufacturer is not responsible for any damages due to misapplication or misuse of this product including, without limitation, direct, incidental and consequential damages, and disclaims such damages to the full extent permitted under applicable law. The user is solely responsible to identify critical application risks and install appropriate mechanisms to protect processes during a possible equipment malfunction.

Please read this entire manual before unpacking, setting up or operating this equipment. Pay attention to all danger and caution statements. Failure to do so could result in serious injury to the operator or damage to the equipment.

Make sure that the protection provided by this equipment is not impaired. Do not use or install this equipment in any manner other than that specified in this manual.

Use of hazard information

▲ DANGER

Indicates a potentially or imminently hazardous situation which, if not avoided, will result in death or serious injury.

▲ WARNING

Indicates a potentially or imminently hazardous situation which, if not avoided, could result in death or serious injury.

▲ CAUTION

Indicates a potentially hazardous situation that may result in minor or moderate injury.

NOTICE

Indicates a situation which, if not avoided, may cause damage to the instrument. Information that requires special emphasis.

Precautionary labels

Read all labels and tags attached to the instrument. Personal injury or damage to the instrument could occur if not observed. A symbol on the instrument is referenced in the manual with a precautionary statement.

	This symbol, if noted on the instrument, references the instruction manual for operation and/or safety information.
	Electrical equipment marked with this symbol may not be disposed of in European public disposal systems after 12 August of 2005. In conformity with European local and national regulations (EU Directive 2002/96/EC), European electrical equipment users must now return old or end-of-life equipment to the Producer for disposal at no charge to the user. <i>Note: For return for recycling, please contact the equipment producer or supplier for instructions on how to return end-of-life equipment, producer-supplied electrical accessories, and all auxiliary items for proper disposal.</i>

Product overview

The sensION™+ series meters are used with sensION+ probes to measure various parameters in water. Primarily for field use, the portable meter operates on three AA batteries.

The sensION™+ series meters are available in 6 models:

- sensION™+ EC5—conductivity, TDS, salinity and temperature meter
- sensION™+ MM150—pH, ORP (Redox), conductivity, TDS and temperature meter
- sensION™+ MM110—pH, ORP (Redox) and temperature meter
- sensION™+ MM156—pH, conductivity, salinity, dissolved oxygen and temperature meter
- sensION™+ pH1—pH meter
- sensION™+ DO6—dissolved oxygen (concentration and % saturation) and temperature meter

Certification

Canadian Radio Interference-Causing Equipment Regulation, IECS-003, Class A:

Supporting test records reside with the manufacturer.

This Class A digital apparatus meets all requirements of the Canadian Interference-Causing Equipment Regulations.

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FCC Part 15, Class "A" Limits

Supporting test records reside with the manufacturer. The device complies with Part 15 of the FCC Rules. Operation is subject to the following conditions:

1. The equipment may not cause harmful interference.
2. The equipment must accept any interference received, including interference that may cause undesired operation.

Changes or modifications to this equipment not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference, in which case the user will be required to correct the interference at their expense. The following techniques can be used to reduce interference problems:

1. Disconnect the equipment from its power source to verify that it is or is not the source of the interference.
2. If the equipment is connected to the same outlet as the device experiencing interference, connect the equipment to a different outlet.
3. Move the equipment away from the device receiving the interference.
4. Reposition the receiving antenna for the device receiving the interference.
5. Try combinations of the above.

Installation

⚠ CAUTION

Multiple hazards. Only qualified personnel must conduct the tasks described in this section of the document.

Install the batteries

⚠ WARNING

Explosion hazard. Incorrect battery installation can cause the release of explosive gases. Be sure that the batteries are of the same approved chemical type and are inserted in the correct orientation. Do not mix new and used batteries.

NOTICE

The battery compartment is not waterproof. If the battery compartment becomes wet, remove and dry the batteries and dry the interior of the compartment. Check the battery contacts for corrosion and clean them if necessary.

NOTICE

When using nickel metal hydride (NiMH) batteries, the battery icon will not indicate a full charge after freshly charged batteries have been inserted (NiMH batteries are 1.2 V versus 1.5 V for alkaline batteries). Even though the icon does not indicate complete charge, 2300 mAh NiMH batteries will achieve 90% of instrument operation lifetime (before recharge) versus new alkaline batteries.

NOTICE

To avoid potential damage to the meter from battery leakage, remove the meter batteries prior to extended periods of non-use.


The meter can be powered with AA alkaline or rechargeable NiMH batteries. To conserve battery life, the meter will power off after 5 minutes of inactivity. This time can be changed in the Display Options menu.

For battery installation, refer to [Figure 1](#).

Note: A tool such as a flat screwdriver may be necessary the first time the cover is removed.

Make sure that the cover is tightly closed to maintain the IP67 enclosure rating.

Figure 1 Battery installation


Connect a probe

1. Plug the probe into the meter ([Figure 2](#)).
2. Push the connector so that the alignment pin goes into the groove on the female connector.

Note: Do not turn the connector.


Figure 2 Probe connection

User interface and navigation

User interface

Figure 3 Keypad description


1 CALIBRATION key	4 DOWN key: scroll to other options, change a value
2 UP key: scroll to other options, change a value	5 ON/OFF: turn on or turn off the meter, DISPLAY LIGHT: turn on and off the display light
3 MEASUREMENT key	

Display description

The display shows 3 types of screens:


- Measure—the most recent parameter measurement is shown with the temperature and stabilization time.
- Calibration—the calibration standard values and temperature are shown during a calibration.
- Standby—the screen is empty except for the measurement timer at the bottom.

Figure 4 pH1 screen display


1 Main measurement value (pH or ORP)	4 Measurement time (hh:mm:ss)
2 Main measurement units	5 Temperature
3 Battery indicator	

Figure 5 DO6 screen display


1 Main measurement value (DO)	4 Measurement time (hh:mm:ss)
2 Main measurement units	5 Temperature
3 Battery indicator	

Figure 6 MM110 screen display


1 Main measurement value (pH)	4 Main measurement time (hh:mm:ss)
2 Measurement units	5 Temperature
3 Battery indicator	6 Main measurement value (ORP)

Figure 7 MM150 screen display


1 Main measurement value (conductivity)	5 Measurement time (hh:mm:ss)
2 Main measurement value (pH)	6 Temperature
3 Main measurement units	7 Main measurement value (TDS, ORP)
4 Battery indicator	8 Reference temperature

Figure 8 MM156 screen display

1 Main measurement value (conductivity)	5 Measurement time (hh:mm:ss)
2 Main measurement value (pH)	6 Temperature
3 Main measurement units	7 Main measurement value (DO, salinity)
4 Battery indicator	8 Reference temperature

Figure 9 EC5 screen display

1 Reference temperature	4 Battery indicator
2 Main measurement units	5 Measurement time (hh:mm:ss)
3 Temperature	6 Main measurement value (conductivity, salinity, TDS)

Navigation

Use the calibrate key  to calibrate the probe. Use the measure key  to take a sample measurement. Use the arrow keys   to scroll to other options or to change a value.

Some options require that a key be pushed and held or that multiple keys be pushed at the same time. Be sure to look at the display during tasks as the screens change quickly. Refer to each task for specific instructions.

Startup

Turn the meter on and off

NOTICE

Make sure that the probe is connected to the meter before the meter is turned on.

Push and hold  to turn on or turn off the meter. If the meter does not turn on, make sure that the batteries are properly installed.

Note: The meter is set to automatically power off after 5 minutes of inactivity to maximize the battery life.

Operation

WARNING

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current material safety data sheets (MSDS) for safety protocols.

Setup for pH

The meter is ready for pH or ORP calibration without additional setup.

Setup for dissolved oxygen

Before calibration or measurement for dissolved oxygen, the probe must be polarized and the settings for atmospheric pressure and salinity must be entered.

Polarize the electrode

If the probe or batteries are disconnected, connect the probe or install the batteries and wait for polarization:

Disconnect time	Polarization time
< 5 minutes	10 minutes
5 to 15 minutes	45 minutes
> 15 minutes	6 hours

Change the settings

The settings for dissolved oxygen measurements can be changed when a dissolved oxygen probe is used. Push  to view the current settings.

1. Push . The first setting is shown.
2. Quickly use the arrow keys to change the value. Wait for the next setting to be shown and change any additional values:

Option	Description
	Atmospheric pressure in mbar (default: 1013 mbar)
SAL	Salinity in g/L NaCl (default: 0 g/L) Aut: automatic (only in model MM156)
SEn	Select the DO sensor model: 5120 or 5130

Note: The salinity value can be found from a conductivity measurement. Refer to [Table 3](#) on page 17.

Setup for conductivity

The settings for conductivity measurements can be changed when a conductivity probe is used. Push  to view the current settings.

1. Push . The first setting is shown.
2. Quickly use the arrow keys to change the value. Wait for the next setting to be shown and change any additional values:

Option	Description
CEL	Electrode type: platinum (default) or titanium. Note: This option is not available on all meters.
CAL	The number of standards for conductivity calibration. Select one standard (default). Select the standard nearest to the expected measurement range.
tC	Temperature compensation: 0 to 9.99%/°C (default: 2%/°C)
tr	Reference temperature: 20 or 25 °C (default: 25 °C)
F	TDS calculation factor: 0.01 to 4.44 (default: 0.64)

The new settings are automatically stored in the meter.

Calibration

▲ WARNING

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current material safety data sheets (MSDS) for safety protocols.

Calibration procedure

This procedure is for general use with liquid calibration solutions. Refer to the documents that are included with each probe for additional information.


1. Pour the buffers or calibration solution into the labeled calibration tubes.
2. Rinse the probe with deionized water and put the probe into the appropriate calibration tube. Make sure that there are no air bubbles under the probe tip.
3. Push . The parameter flashes.
4. Use the arrow keys to change the parameter, if applicable.
5. Push  to select the parameter.
6. Push the arrow keys to select the appropriate standard solution. For pH, the buffers are automatically recognized.
7. Push  to measure the calibration solution.
Note: For pH, the next calibration solution is shown. For conductivity, use one calibration point. Select the standard nearest to the expected measurement range.
8. Rinse the probe with deionized water and put the probe into the second calibration tube, if applicable. Make sure that there are no air bubbles under the probe tip.
9. Push  to measure the second calibration solution, if applicable. The next calibration solution is shown.
10. Rinse the probe with deionized water and put the probe into the third calibration tube. Make sure that there are no air bubbles under the probe tip.
11. Push  to measure the third calibration solution, if applicable.

When the calibration is good, the display briefly shows OK and then goes to the standby mode.

Note: To calibrate with only 1 or 2 standards when additional standards are supplied, push  after the first or second standard is measured.

View the calibration data

Data from the most recent calibration can be shown for pH, ORP and conductivity.

1. Push .
2. Use the arrow keys to change the parameter, if necessary, and push .
3. Push  and  at the same time. The calibration data is shown:
 - pH—the slope and offset values are shown alternating with the deviation (in %) and calibration temperature.
 - ORP—the measured mV value and calibration temperature are shown.
 - Conductivity—the cell constant and calibration temperature for each standard are shown.

Restore the factory calibration

The user calibration can be erased and the factory calibration restored for pH, ORP and conductivity.

1. Push . The parameter flashes.
2. Use the arrow keys to change the parameter, if necessary.
3. Push and hold  until **OFF** is shown.
The factory calibration for the selected parameter is restored.

Adjust the calibration

The instrument can be adjusted to read a specified solution value for pH, ORP and conductivity.


1. Rinse the probe. Put the probe in the solution.
2. Push  to read the value of the solution.
3. Push . The parameter flashes.

4. Use the arrow keys to select the parameter, if necessary.
5. Push and hold  and then push .
The measured value flashes.
6. Use the arrow keys to change the value.
7. Push . The display shows OK.

Adjust the temperature

The temperature measurement can be adjusted at 25 °C (77 °F) and/or 85 °C (185 °F) to increase accuracy for pH, ORP and conductivity.

Note: Refer to the probe documentation to find the temperature limits for the probe that is used.

1. Put the probe and a reference thermometer in a container of water at approximately 25 °C and allow the temperature to stabilize.
2. Compare the temperature read by the meter with that of the reference thermometer. The difference is the adjustment value for the meter.
Example: reference thermometer: 24.5 °C; meter: 24.3 °C.
Adjustment value: 0.2 °C.
3. Enter the adjustment value for the 25 °C reading:
 - a. Push . The parameter flashes.
 - b. Use the arrow keys to select **tEn**.
 - c. Push . 25 °C is shown.
 - d. Push  and then use the arrow keys to enter the adjustment value for 25 °C. Push . 85 °C is shown.
4. Put the probe and a reference thermometer in a container of water at approximately 85 °C and allow the temperature to stabilize.
5. Compare the temperature from the meter with that of the reference thermometer. The difference is the adjustment value for the meter.
6. Push  and then use the arrow keys to enter the adjustment value for the 85 °C reading. Push .

About sample measurements


Each probe has specific preparation steps and procedures for taking sample measurements. For step-by-step instructions, refer to the documents that are included with the probe.

Push  to take a sample measurement. During measurements, the parameter flashes and the timer shows the stabilization time. To change the measurement parameter (if applicable), push and hold .

To measure the sample continuously, push  2 times. The parameter flashes to indicate the continuous measurement mode.

Change the temperature units

The temperature units can be changed when the measurement screen is shown.

1. Make sure that a stable measurement reading is shown.
2. Push  and  at the same time. The temperature units change to °C or °F.

Disable the auto-shutoff option

The meter is set to automatically power off after 5 minutes of inactivity to maximize the battery life. This option can be temporarily disabled.


1. Make sure that the meter is powered off.
2. Push and hold  until the display shows **OFF** and **no**. The power remains on indefinitely.
3. To power off the meter, push and hold  until the display shows **OFF**.

Note: The auto-shutoff option is enabled again when the meter is powered on normally.

Set the backlight time

The display is illuminated when  is pushed. The user can set the length of time that the backlight stays on.

Note: The battery life is decreased when the backlight time is increased.

1. Make sure that the instrument is powered on.
2. Push  and immediately push  2 times. The backlight time flashes.
3. Push  or  to change the backlight time (range: 3 seconds to 2 minutes).

Maintenance

⚠ WARNING

Multiple hazards. Do not disassemble the instrument for maintenance or service. If the internal components must be cleaned or repaired, contact the manufacturer.

⚠ CAUTION

Personal injury hazard. Only qualified personnel should conduct the tasks described in this section of the manual.

Clean the instrument

NOTICE

Never use cleaning agents such as turpentine, acetone or similar products to clean the instrument including the display and accessories.

Clean the exterior of the instrument with a moist cloth and a mild soap solution.

Store the meter

NOTICE

To avoid potential damage to the meter from battery leakage, remove the meter batteries prior to extended periods of non-use.


Replace the batteries

⚠ WARNING

Explosion hazard. Incorrect battery installation can cause the release of explosive gases. Be sure that the batteries are of the same approved chemical type and are inserted in the correct orientation. Do not mix new and used batteries.

For battery replacement, refer to [Figure 10](#).

Figure 10 Battery replacement

Troubleshooting

Refer to the following table for common problem messages or symptoms, possible causes and corrective actions.

Error/Warning	Description	Solution
 	Measurement out of range.	<ul style="list-style-type: none"> • Insert the probe into the appropriate standard and read again. • Clean the probe and the membrane. • Make sure that there are no air bubbles in the membrane. • Disconnect the probe and then connect the probe again. • Connect a different probe to verify if problem is with probe or meter.
	Battery is low.	Insert new batteries.
E1	Unstable reading during a stability measurement or during calibration.	Make sure that the probe is properly immersed in the sample.
E2	Probe current: exceeded upper limit during measurement ≥ 250 nA.	<ul style="list-style-type: none"> • Clean the probe and the membrane. • Make sure that there are no air bubbles in the membrane. • Disconnect the probe and then connect the probe again. • Connect a different probe to verify if problem is with probe or meter.
E2	pH probes sensitivity out of range (accepted values 70–105%).	<ul style="list-style-type: none"> • Make sure that the correct probe is connected to the meter. • Connect a new probe.

Error/Warning	Description	Solution
E2	During calibration, the difference between the constants gets $> 30\%$.	<ul style="list-style-type: none"> • Insert the probe into the appropriate standard and read again. • Clean the probe and the membrane. • Make sure that there are no air bubbles in the membrane. • Disconnect the probe and then connect the probe again. • Connect a different probe to verify if problem is with probe or meter. • Connect a new probe.
E3	During calibration, the cell constant gets $< 0.05 \text{ cm}^{-1}$.	<ul style="list-style-type: none"> • Insert the probe into the appropriate standard and read again. • Clean the probe and the membrane. • Make sure that there are no air bubbles in the membrane. • Disconnect the probe and then connect the probe again. • Connect a different probe to verify if problem is with probe or meter.
E3	Slope out of range (accepted values $\pm 58 \text{ mV}$).	Connect a new probe.
E4	Buffer solution not recognized.	<ul style="list-style-type: none"> • Clean the probe and the membrane. • Make sure that there are no air bubbles in the membrane. • Disconnect the probe and then connect the probe again. • Connect a different probe to verify if problem is with probe or meter. • Make sure that the buffer used matches the buffer specified in configuration. • Replace the buffer solution. • Make sure of the temperature specification in configuration.

Error/Warning	Description	Solution
E4	Salinity cannot be calculated TC=0	Modify TC
E5	Same buffer solutions.	<ul style="list-style-type: none"> • Clean the probe and the membrane. • Make sure that there are no air bubbles in the membrane. • Disconnect the probe and then connect the probe again. • Connect a different probe to verify if problem is with probe or meter. • Make sure that the buffer used matches the buffer specified in configuration. • Replace the buffer solution. • Make sure of the temperature specification in configuration.
E6	The calibration solutions have different temperatures.	Make sure that the calibration solutions have the same temperature.
ErA (only MM156)	Salinity is set to automatic and measured value > 50 g/L.	Modify TC

Replacement parts and accessories

Note: Product and Article numbers may vary for some selling regions. Contact the appropriate distributor or refer to the company website for contact information.

Replacement parts

Description	Item no.
Batteries, Alkaline AA	1938004
Carrying case	LZW9990.99

Accessories

Description	Item no.
Electrolyte for DO probe 51 20, 25 mL	LZW9811.99
Electrolyte for DO probe 51 30, 50 ml	2759126
Tube for DO probe calibration	LZW5123.99
Service kit for 51 30 DO probe (contains 2 membrane modules and DO filling solution)	5196800
Replacement DO membrane for DO probe 51 20	LZW5125.99
147 µS/cm standard solution (at 25 °C, 77 °F), 125 mL bottle	LZW9701.99
1413 µS/cm standard solution (at 25 °C, 77 °F), 125 mL bottle	LZW9711.99
12.88 mS/cm standard solution (at 25 °C, 77 °F), 125 mL bottle	LZW9721.99
Set of calibration tubes, 147 µS/cm, 1413 µS/cm and 12.88 mS/cm	LZW9138.99
220 mV standard solution, 125 mL	LZW9402.99
Buffer solution pH 4.01, 125 mL	LZW9460.99
Buffer solution pH 7.00, 125 mL	LZW9461.97
Buffer solution pH 10.01, 125 mL	LZW9470.99
Electrolyte solution, KCl 3M, 50 mL	LZW9509.99
Calibration tubes pH 4.01, pH 7.00 and pH 10.01	LZW9137.97
Tube for ORP standard calibration 220 mV	LZW9136.99
Work protector for electrodes 50 51 T and 50 52 T	LZW9162.99

Standard solutions

Technical buffer solutions (DIN 19267)

Refer to [Table 1](#) pH and ORP (mV) values of specific buffer sets at varying temperatures.

Table 1 pH, ORP (mV) and temperature values

Temperature		pH					mV
°C	°F						
0	32	2.01	4.01	7.12	9.52	10.30	—
10	50	2.01	4.00	7.06	9.38	10.17	245
20	68	2.00	4.00	7.02	9.26	10.06	228
25	77	2.00	4.01	7.00	9.21	10.01	220
30	86	2.00	4.01	6.99	9.16	9.96	212
40	104	2.00	4.03	6.97	9.06	9.88	195
50	122	2.00	4.06	6.97	8.99	9.82	178
60	140	2.00	4.10	6.98	8.93	9.76	160
70	158	2.01	4.16	7.00	8.88	—	—
80	176	2.01	4.22	7.04	8.83	—	—
90	194	2.01	4.30	7.09	8.79	—	—

Conductivity standard solutions

Refer to [Table 2](#) for conductivity values of standard solutions at varying temperatures.

Table 2 Conductivity and temperature values

Temperature		Conductivity (EC)			
°C	°F	µS/cm	µS/cm	mS/cm	mS/cm
15.0	59	119	1147	10.48	92.5
16.0	60.8	122	1173	10.72	94.4
17.0	62.6	125	1199	10.95	96.3
18.0	64.4	127	1225	11.19	98.2
19.0	66.2	130	1251	11.43	100.1

Table 2 Conductivity and temperature values (continued)

Temperature		Conductivity (EC)			
°C	°F	µS/cm	µS/cm	mS/cm	mS/cm
20.0	68	133	1278	11.67	102.1
21.0	69.8	136	1305	11.91	104.0
22.0	71.6	139	1332	12.15	105.4
23.0	73.4	142	1359	12.39	107.9
24.0	75.2	145	1386	12.64	109.8
25.0	77	147	1413	12.88	111.8
26.0	78.8	150	1440	13.13	113.8
27.0	80.6	153	1467	13.37	115.7
28.0	82.4	156	1494	13.62	—
29.0	84.2	159	1522	13.87	—
30.0	86	162	1549	14.12	—
31.0	87.8	165	1581	14.37	—
32.0	89.6	168	1609	14.62	—
33.0	91.4	171	1638	14.88	—
34.0	93.2	174	1667	15.13	—
35.0	95	177	1696	15.39	—

Salinity values for dissolved oxygen meters

Table 3 can be used to convert conductivity values to salinity for dissolved oxygen meters.

Table 3 Conversion from conductivity to salinity

Conductivity (mS/cm)	Salinity (g/L NaCl)
1.9	1.0
3.7	2.0
7.3	4.0
10.9	6.0
17.8	10.1
25.8	15.1
33.6	20.2
41.2	25.4
48.9	30.6
56.3	35.8
62.2	40.0
69.4	45.3
75.7	50.0

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P.O. Box 389, Loveland, CO 80539-0389 U.S.A.
(800) 227-4224 (U.S.A. only)