20000005226

# IKA®

Dry Block Heater\_022016

## IKA® Dry Block Heater 1 IKA® Dry Block Heater 2 IKA® Dry Block Heater 3 IKA® Dry Block Heater 4

Indicaciones de seguridad	ES	31
Veiligheidsinstructies	NL	32
Norme di sicurezza	IT	33
Säkerhetsanvisningar	SV	34
Sikkerhedshenvisninger	DA	35
Sikkerhedsanvisninger	NO	36
Turvallisuusohjeet	FI	37
Normas de segurança	PT	38
Υποδείξεις ασφαλείας	EL	39

Betriebsanleitung	DE	3
Ursprungssprache		
Operating instructions	EN	10
Mode d'emploi	FR	17
Инструкция по эксплуатации	RU	24
Wskazówki bezpieczeństwa	PL	40
Bezpečnostní upozornění	CS	41
Biztonsági utasítások	HU	42
Varnostna navodila	SL	43
Bezpečnostné pokyny	SK	44
Ohutusjuhised	ET	45
Drošības Norādes	LV	46
Saugos Reikalavimai	LT	47
Инструкции за безопасност	BG	48
Indicații de siguranță	RO	49

IKA®-Werke, Germany  
Reg. No. 004343

**IKA® Dry Block Heater 1**  
**IKA® Dry Block Heater 2**  
**IKA® Dry Block Heater 3**  
**IKA® Dry Block Heater 4**


Fig. 1

Pos.	Bezeichnung	Item	Designation	Pos.	Désignation
A	Netzbuchse	A	Power socket	A	Prise secteur
B	KT-Buchse	B	KT socket	B	Prise KT
C	Stativgewindebohrung	C	Threaded hole for stand	C	Alésage fileté du statif
D	Kontaktstecker	D	Contact plug	D	Fiche de contact
E	Geräteschalter	E	Main switch	E	Commutateur
F	Taste "Start/Stop"	F	"Start/Stop" button	F	Touche "Start/Stop"
G	Taste Timer "+"	G	Button timer "+"	G	Touche Timer "+"
H	Taste Timer "-"	H	Button timer "-"	H	Touche Timer "-"
I	Taste Timer/ Counter	I	Button Timer/ Counter	I	Touche Timer/ Counter
J	Taste Heizung "+"	J	Button heating "+"	J	Touche chauffage "+"
K	Taste Heizung "-"	K	Button heating "-"	K	Touche chauffage "-"
L	LED Display Heizung	L	LED Display heating	L	Affichage DEL chauffage
M	LED (Set=Sollwert)	M	LED (set=set point value)	M	DEL (set=valeur théorique)
N	LED externer Temperatursensor	N	LED external temperature sensor	N	DEL capteur de température externe
O	LED Heizplatte	O	LED heating plate	O	DEL plaque chauffante
P	LED Display Timer	P	LED display timer	P	Affichage DEL minuteur

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## Contents

	Page
Declaration of conformity	10
Explication of warning symbols	10
Safety instructions	11
Unpacking	12
Correct use	12
Assembling the stand	12
Operation	13
Setting the operating mode	13
Controlling the temperature via contact thermometer	14
Maintenance and cleaning	14
Error codes	15
Technical data	16
Warranty	16

## Declaration of conformity

We declare under our sole responsibility that this product corresponds to the regulations 2006/95/EEC and 2004/108/EC and conforms with the standards or standardized documents EN 61010-1, -2-010 and EN 61326-1.

## Explication of warning symbols

General hazard

This symbol identifies information **that is of vital importance for safeguarding your health and safety**. Disregarding this information can lead to health impairment and injuries.

This symbol identifies information **that is of importance for the technically correct functioning of the system**. Disregarding this information can result in damage to the appliance or to system components.

This symbol indicates information **which is important for proper use of the appliance and/or ensuring that the appliance functions correctly**. Failure to observe this information can lead to inaccurate results.

Danger of being burnt! - Note on hazards arising from a hot surface.

# Safety instructions

## General safety instructions

- **Read the operating instructions in full before starting up and follow the safety instructions.**
- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations
- Socket must be earthed (protective ground contact).
- Before starting the device for the first time, ensure that the connector cable is suitable for the type of socket used and that a safe protective earth connection is assured.

**⚠ DANGER** **Risk of burns!** Exercise caution when touching the housing parts and the heating plate. The heating plate can reach temperatures in excess of 120 °C. Pay attention to the residual heat after switching off. **The device may not be transported unless it is cold!**

## Safety instructions - device arrangement

**⚠ DANGER** **Do not** operate the appliance in explosive atmospheres, with hazardous substances or under water.

- Set up the appliance in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.
- The feet of the appliance must be clean and undamaged.
- The voltage stated on the type plate must correspond to the mains voltage.
- The socket for the mains cord must be easily accessible.
- The power supply cable and cables to the external sensors must not be allowed to come into contact with the hot mounting plate.
- Check the appliance and accessories beforehand for damage each time you use them. Do not use damaged components.

## Safety instructions - permissible media

**⚠ WARNING** **Caution!** Only media whose flashpoint lies above the safety temperature limit of 130 °C may be processed or heated with this device.

**⚠ WARNING** Beware of hazards due to:

- flammable materials
- combustible media with a low boiling temperature
- glass breakage
- incorrect container size
- overfilling of media (in a beaker)
- unsafe condition of container.

- Process pathogenic materials only in closed vessels under a suitable extractor hood. Please contact **IKA®** if you have any questions.

**⚠ DANGER** Only process media that will not react dangerously to the extra energy produced through processing. This also applies to any extra energy produced in other ways, e.g. through light irradiation.

- Bear in mind the possibility of contamination that might lead to undesirable chemical reactions.

## Safety instructions - performance of trials

**⚠ WARNING** Wear your personal protective equipment in accordance with the hazard category of the media to be processed. Otherwise

there is a risk from:

- splashing and evaporation of liquids
- ejection of parts
- release of toxic or combustible gases.

- Immerse external temperature sensors (PT 1000 ...) at least 20 mm deep into the medium or insert the sensor into the hole provided for the purpose in the temperature control block.
- The PT 1000 external temperature sensor must always be inserted in the block / media when connected.
- Do not heat liquid media except with the block. Never pour liquid media directly on to the heating plate!
- Safe operation is only guaranteed with the accessories described in the "Accessories" chapter.
- Always disconnect the plug before fitting accessories.
- Accessories must be securely attached to the device and cannot come off by themselves. The centre of gravity of the assembly must lie within the surface on which it is set up.
- Refer to the operating instructions for the accessories.

## Safety instructions - power supply/ switching off the device

**⚠ WARNING** The appliance can only be disconnected from the mains supply by pulling out the mains plug or the connector plug.

- The device will automatically restart in mode B following any interruption to the power supply.

## For protection of the equipment

- The appliance may only be opened by experts.
- Do not cover the device, even partially e.g. with metallic plates or film. This results in overheating.
- Ensure that the base plate is kept clean.
- Protect the appliance and accessories from bumps and impacts.
- Observe the minimum distances between the devices, between the device and the wall and above the assembly (min. 800 mm), see **Fig. 2**.


Fig. 2

## Unpacking

### • Unpacking

- Please unpack the device carefully
- In the case of any damage a fact report must be sent immediately (post, rail or forwarder)

### • Delivery scope

- IKA® Dry Block Heater,  
Depending on the IKA® Dry Block Heater model 1.2, 3, or 4
- Mains cable
- Temperature sensor PT 1000.60
- Operating instructions

## Correct use

### • Use

- For heating media in block heaters

### • Range of use (interior areas only)

- Laboratories                      - Schools
- Pharmacies                        - Universities

The safety of the user cannot be guaranteed if the appliance is operated with accessories that are not supplied or recommended by the manufacturer or if the appliance is operated improperly contrary to the manufacturer's specifications or if the appliance or the printed circuit board are modified by third parties.

## Assembling the stand

- ☞ Remove the protective cap from the support rod.
- ☞ Put the washer between housing and nut.
- ☞ Screw the support rod onto the device by hand until the end stop is reached.

- ☞ Using an open-jaw spanner 17 AF **gently** tighten the nut (**tightening torque max. 0.5 Nm**).
- ☞ Accessories should be attached using cross sleeves.


Fig. 3

## Operation

	IKA® Dry Block Heater 1	IKA® Dry Block Heater 2	IKA® Dry Block Heater 3	IKA® Dry Block Heater 4
<b>Commissioning</b>	<ul style="list-style-type: none"> <li>☞ Move device switch <b>(E)</b> to the OFF position</li> <li>☞ Insert the mains power cable into the power socket <b>(A)</b></li> <li>☞ Move device switch <b>(E)</b> to the ON position</li> <li>⇒ The unit will be set to factory setting mode A (see "Operating modes")</li> </ul>			
<b>Heating</b>	<ul style="list-style-type: none"> <li>☞ Set the target temperature using the heater buttons +/- <b>(J)</b> and <b>(K)</b></li> <li>⇒ The value selected will be shown on the display <b>(L)</b></li> <li>☞ Start the heating function by pressing the "Start / Stop" button <b>(F)</b></li> </ul>			
	 	<p>The set-point and actual temperatures will be shown alternately on the display (<b>(L)</b> is blinking):</p> <ul style="list-style-type: none"> <li>⇒ When the heating is on, the LED <b>(O)</b> lights up.</li> <li>⇒ During standby operation, the display <b>(L)</b> will show <span style="border: 1px solid black; padding: 0 2px;">HOT</span> if the base plate temperature is above 50 °C</li> </ul>		
<b>Connecting external thermometers</b> <small>(direct temperature control in the media)</small>	<ul style="list-style-type: none"> <li>☞ Move device switch <b>(E)</b> to the OFF position</li> <li>☞ Detach contact plug <b>(D)</b></li> <li>☞ Attach a DIN 12878 (Class 2) compliant safety contact thermometer or a PT 1000.60 temperature sensor to the socket <b>(B)</b></li> <li>☞ Move device switch <b>(E)</b> to the ON position</li> </ul>			
		<p>Temperature sensor PT 1000</p> <ul style="list-style-type: none"> <li>⇒ The actual temperature for the temperature sensor shown on display <b>(L)</b> will correspond to the temperature of the media.</li> <li>The LED "external temperature sensor" PT 1000 <b>(N)</b> will be lit.</li> </ul>		

## Setting operating mode

### Operating the device in mode A or B

#### Mode A

All settings will be stored if the device is switched off or disconnected from the power supply. The heating functions will be set to OFF when the device is powered on.

#### Mode B

All settings will be stored if the device is switched off or disconnected from the power supply. The heating function will be set to ON or OFF when the device is powered on, depending on the previous status of the device..

*Factory setting: mode A*

The mode selected will be shown on the display when the device is started up.

### Change the mode

- ☞ Move device switch **(E)** to the OFF position
- ☞ Keep the selection button **(F)** depressed
- ☞ Move device switch **(E)** to the ON position
- ☞ Release the selection button **(F)**

 ⇒ The set mode is indicated on the Display **(P)**   b

## Timer/Counter

After the device has been switched on, the "Timer/ Counter" button **(I)** can be used to toggle between timer mode and counter mode.

### Timer

Pressing the buttons "+" **(G)** and "-" **(H)** allows selection of a time period of up to 99 hours and 59 minutes. Once the heating function **(F)** is activated, the preset time starts to run. When the total time has run down to 00 hours and 00 minutes, a warning buzzer sounds for 30 seconds and in addition the display flashes. The heating function is then deactivated.

### Counter

When the heating function is started by pressing the button **(F)** the counter is activated. In counter mode the time is counted upwards from 00 hours 00 minutes until the heating function is stopped. The counter always starts from 00:00 hh:mm.

### Selection of Timer/Counter

Pressing the button **(I)** toggles between timer mode and counter mode. The mode can be activated only when the heating function is off.

## Controlling the temperature via contact thermometer

The preferable method for controlling the average temperature is with the contact thermometer. After the set point temperature has been adjusted, this results in a short heating-up time, practically no temperature drift and only minor fluctuation in temperature.

A 6-pin jack is located on the rear side of the instrument for connecting the contact thermometer or the contact plug. The electronics of the instruments returns a test current that must flow via connector pins 3 and 5 for the heating plate to heat up.

### Settings

For the detailed instructions for settings and limit values, refer to the operating instructions for the instrument you are connecting.

**The target temperature of the medium is set at the contact thermometer. The heater buttons +/- **(J)** and **(K)** on the device can be used to preselect the required surface temperature of the heating plate.**

### Safety contact thermometers

acc. to DIN 12 878 class 2 or acc. Gerstel are connected with a three-wire cable, the test current flows through the contact thermometer.

### Safety function

If the test current is interrupted because of e.g. breakage of contact thermometer or falling out of the cable plug, the heating cuts off. The required plug allocation is indicated on a reference plate above the diode bushing.


### Contact thermometer without safety circuit

acc. to DIN 12 878 class 0. The instrument only heats if the test current circuit is closed by an electrical connection of the plug pins 3 and 5.

## Maintenance and cleaning

The device is maintenance-free.

### Cleaning

For cleaning disconnect the main plug.

Only use cleansing agents which have been recommended by **IKA®**.

Dyes	Isopropyl alcohol
Construction materials	Water containing surfactant, isopropyl alcohol,
Cosmetics	Water containing surfactant, isopropyl alcohol
Foodstuffs	Water containing surfactant
Fuels	Water containing surfactant

- Do not allow moisture to get into the appliance when cleaning.
- Wear protective gloves when cleaning the devices.
- Please consult with **IKA®** before using any cleaning or decontamination methods, other than those recommended here.

### Ordering spare parts

When ordering spare parts, please give:

- Machine type
- Manufacturing number, see type plate
- Item number and designation of the spare part, see **www.ika.com**
- Software version

### Repair

In case of repair the device has to be cleaned and free from any materials which may constitute a health hazard. Please request the form to "clearance certificate" from IKA. Or, use the download-print the form to the **IKA®** website **www.ika.com**.

If you require servicing, return the appliance in its original packaging. Storage packaging is not sufficient. Please also use suitable transport packaging.

## Error codes

Any malfunctions during operation will be identified by an error message on the display (LED **(P)**).

Proceed as follows in such cases:

- ☞ Switch off device using the main switch **(E)**
- ☞ Carry out corrective measures
- ☞ Restart device

Error code	Cause	Effect	Solution
E 3	Temperature inside device too high	Heating off	Switch off device and allow to cool down
E 6	Break in safety circuit	Heating off	<ul style="list-style-type: none"> <li>- Plug in contact plug <b>(D)</b></li> <li>- Plug in PT 1000 contact thermometer / temperature sensor</li> <li>- Replace faulty connecting cable, plug, or contact thermometer</li> </ul>
E 25	Heating and switching element monitoring	Heating off	<ul style="list-style-type: none"> <li>- Switch off device</li> <li>- <b>Warning!</b> Only to be carried out by authorised service personnel: Carry out an internal test on the device to check the plug-in connector for the heating element</li> </ul>
E 26	Difference between temperature of safety sensor and temperature of control sensor <b>control temperature</b> > (safety temperature + 40 K)	Heating off	<ul style="list-style-type: none"> <li>- Switch off device</li> <li>- <b>Warning!</b> Only to be carried out by authorised service personnel: Carry out an internal test on the device to check the temperature sensor</li> </ul>

If the actions described fail to resolve the fault or another error code is displayed then take one of the following steps:

- contact the **IKA**® service department,
- send the device for repair, including a short description of the fault.

## Technical Data

		IKA® Dry Block Heater 1	IKA® Dry Block Heater 2	IKA® Dry Block Heater 3	IKA® Dry Block Heater 4
Operating voltage range	Vac	220 - 230 ± 10%			
Rated voltage	Vac	115 ± 10%			
	Vac	100 ± 10%			
	Vac	230/ 50 Hz			
Frequency	Vac	115/ 60 Hz			
	Hz	100/ 60 Hz			
Display		50 / 60			
Temperature range	°C	digital			
Timer	min/h	ambient temperature +5 - 120			
Perm. duration of operation	%	1 min to 99 h 59 min			
Perm. ambient temperature	°C	100			
Perm. humidity	%	+5 to +40			
Protection acc. to DIN EN 60529		80			
Protection class		IP 21			
Overvoltage category		I			
Contamination level		II			
Operation at a terrestrial altitude	m	2			
Dimension (W x D x H)	mm	151 x 228 x 73	151 x 304 x 73	151 x 380 x 73	151 x 456 x 73
Weight	kg	1	1.3	1.6	2
<b>Heating</b>					
Power consumption (+10%) max. at 230 Vac	W	165	250	330	412
	115 Vac	165	250	330	412
	100 Vac	165	250	330	412
Temperature sensor PT 1000 variation DIN EN 60751 Kl. A	K	≤ ± 0.1*			
Temperature stability within the blocks at 37 °C	°C	± 0.2	± 0.2	± 0.2	± 0.2
Temperature stability within the blocks at 60 °C	°C	± 0.2	± 0.2	± 0.2	± 0.4
Temperature homogeneity within the the blocks at 37 °C	°C	± 0.2	± 0.2	± 0.3	± 0.3
Temperature homogeneity within the the blocks at 60 °C	°C	± 0.4	± 0.4	± 0.6	± 0.6
Heating rate (PT 1000 in solid block)	K/min	5.0	4.5	4.0	4.0
<b>Fixed safety circuit</b>					
Safety temperature limit	°C	150			

\* Test parameter: with test block and glass vessel filled with M50 oil

Subject to technical changes!

## Warranty

In accordance with IKA® warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our factory, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover worn out parts, nor does it apply to faults resulting from improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating manual.

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**IKA®-Werke GmbH & Co.KG**