



IRIS OFFERS YOU THE BEST TECHNICAL SPECIFICATIONS, MAKING IT A COMPELLING CHOICE FOR YOUR RESEARCH.

BORN FROM A PASSION TO REVEAL THE UNSEEN, IRIS IS AN INNOVATIVE INSTRUMENT FOR YOUR NEXT SCIENTIFIC EXPLORATIONS.

GAIN CONFIDENCE IN YOUR RESULTS, PUSH THE BOUNDARIES OF YOUR PROTEIN DETECTION,



Uvitec Cambridge

WWW.UVITEC.CO.UK UVI@UVITEC.CO.UK



Ask for a demo





"Nothing in life is to be feared, it is only to be understood."

Marie Curie (1867-1934)

In honor of Rosalind Franklin and the many women in science.





In 1953, at the heart of Cambridge University, James Watson and Francis Crick discovered DNA double helix structure based on Rosalind Franklin famous "Photograph 51".

This breakthrough not only unveiled the secrets of genetic information; it also brought certainty on the very essence of life itself.

Today, following in the footsteps of these great scientists, we are committed to provide you certainty by guaranteeing quality images with IRIS.

The new Alliance IRIS is the latest generation top end imaging system on the market for chemiluminescence and fluorescence Western blots. IRIS not only presents the highest optic sensitivity with an aperture of f/0.75 but also strong innovative features for unmatched precision and detection.



UVITEC was born in Cambridge and inspired by the future. With years lab experience, our product design team are very familiar with researchers daily workflows. Today, practical and user-oriented design are always at the core of our developments.

Red is the new black

The standalone design of IRIS is a testament to its practicality and user-friendliness. It is easy to clean and compact to fit seamlessly into your laboratory environment. Inspired by the red square in our logo and the distinctive shape of the DNA molecule, its design connects IRIS to the very essence of molecular research, symbolizing the pursuit of knowledge and discovery.

Get the Cambridge Touch

Enjoy the comfort of the widest touch screen. With its impressive size of over 15.6 inches (40 centimeters), our touchscreen offers an immersive and smooth viewing experience. The screen is also adjustable, with a remarkable resolution to ensure that every detail of your molecular imaging experiment is vividly displayed.

Avoid time consuming manipulation

Our mobile sample tray has been designed to avoid timeconsuming manipulations ensuring that you can handle your samples with ease. Additionally, our Alliance IRIS imager features a full door aperture and slide-out tables, specifically engineered to alleviate the discomfort usually associated with draw-out transilluminators.



Don't miss the next UV-LED Tables generation

In the spirit of a sustainable environment, UVITEC Cambridge is the first to introduce a UV-LED technology to ensure a smooth transition from our traditional UV tables. UV-LEDs have a longer lifespan: they can operate for thousands of hours without significantly losing light intensity and efficiency. In addition, these LED conditions prevent damage to your most sensitive samples.

IRIS IS INNOVATIVE

UV-LED Table









IRIS IS RELIABLE





Epifluorescence users seek consistent and uniform illumination in their applications.

A brief history of excitation light sources

Traditional imagers that you will find on the market often rely on the following illumination modes.

White light:

The very first epifluorescence imaging systems started by using white light-based technology converting white light into red, green and blue via filters. However, multiplexing is impossible because the quantity of light is divided by three.

Spot LED:

The next generation and most commonly found epifluorescence imaging systems use Spot LED technology. However, NIR/IR excitation channels originate from the same excitation source, decreasing significantly the light intensity in these channels. Depending on the manufacturer, filter quality varies, influencing crosstalk and multiplexing capacities.

Scanner:

The introduction of laser-based technology raised the interest towards near-infrared and infrared applications.

However, this technology is expensive and requires long acquisition time, and is limited to 2 channels of excitation.

Find the right light source

Our well-known illumination system overcomes these challenges, providing strong light and multiplexing possibilities. Our unique Chromascan concept goes beyond, offering the best of our technologies.

Discover the next generation illumination

Through our new technology, IRIS offers more accurate quantification, and fewer crosstalk issues. The acquisition time is rapid. IRIS leaves no room for uncertainties by guaranteeing that your samples are scanned in the same way, everywhere, with a choice of up to 12 Chromascan modules packs.

Choose a trusted companion

IRIS is made to last. In line with our top-end commitment and quality-oriented values, our system is entirely made of stainless steel and recycled components.



IRIS IS INTUITIVE

A revolutionary molecular imager

Experience effortless data analysis. IRIS integrates an unprecedented technology in the world of molecular imaging. IRIS sets up your data as never seen before and allows you to go straight to analysis. The integrity of your raw data is maintained. IRIS serves as a faithful assistant, automating routine tasks.

Get your images in one click

Auto is our motto. Enjoy our 1-click acquisition process with fully automatic exposure, lighting and focus mode. Our software allows you to get outstanding quantifiable images and perfect band saturation level. IRIS comes with a complete license free software analysis pack (molecular weight, protein quantification or distance calculations) that you can install on as many computers as necessary.

Upgrade your system at anytime

IRIS is fully customizable, tailored to your workflow, and upgradable. At UVITEC, to answer the various numbers of epifluorescent dyes that may be used in a laboratory, we have designed 12 different packs including excitation and emission modules for you to easily insert into your system at any time you decide to start working on epifluorescence.









IRIS IS SENSITIVE

Just as the human iris opens a window to the world, our IRIS imager opens a portal to a molecular universe where your protein of interest is revealed in unprecedented detail. Combining all the parameters presented below, IRIS guarantees sensitivity.

Collect more light

IRIS has its eves wide open: its lens features an impressive f/0.75 aperture, which allows it to capture an incredible amount of light. As you may know, in optics, the smaller the focal number, the more light your device collects.

Get the best quantifiable pictures

IRIS provides a remarkable resolution of up to 30 megapixels To maintain its sensitivity, IRIS relies on a sophisticated cooling (9.2MP native), which ensures that you can discern even the finest system, using a three-stage Peltier normalized camera. This details in your imaging experiments. In addition, our system cooling system allows you to minimize background noise and provides pictures with high density of grayscale within the camera, detect your lowest signals. enabling you to precisely analyze your samples. The closer you



get to 65,535 grev levels, the more quantifiable information you have access to.

Detect even the weakest signals

IRIS' optical system is further enhanced by its exceptional dynamic range, with an optical density (OD) of 4.8. This means that IRIS can capture a wide range of signal intensities without losing crucial information. Whether you are working with faint signals or intense ones, IRIS can accommodate your needs with ease.

Keep it cool



Practicality

1-click to image > effortless acquisition, in no time **Full automation** > hands-off, automated routines **Interchangeable Tables** > hassle-free sample positioning **Innovative software** > revolutionized analysis

• Detection capabilities

Chromascan concept > up to 8 fluorescent excitation light **Uvipure technology** > enhanced UV for EtBr and all safe stains **Confocal discs** > boosted detection of fluorescent signals **Multiplexing** > imaging of several proteins

THEY TRUST US

Imaging

Up to 30 megapixels> 9.2MP native camera resolution f/0.75 custom lens > unrivalled camera sensitivity **3 stages Peltier Cooling** > reduction of background noise **OD 4.8 dynamic range** > outstanding detection

Design

Standalone system > with integrated PC Comfort > widest 15.6 inches touchscreen Sustainable UV > LED 312nm - UV tubes available Robustness > stainless steel and epoxy paint

Dims > Height: 672 mm - Width: 442 mm - Depth 534 mm - **Weight** > 70 Kg - **FOV** > 21 x 26 cm

vvestern blotting	Chemiluminescence
DNA and RNA gels with fluorescent stains	Fluorescence
Optional selection of NIR / RGB modules from 12 combinations available	Epi Fluorescence
In-vivo Luciferase and Fluorescence	Bioluminescence
Optional Colorimetry and Protein gels	Visible Imaging















More than 10,000 users worldwide