## Kunststoff - Dispergierwerkzeuge Plastic dispersion tools Outils dispersants

S18D - 14G - KS







S18D - 10G - KS

S25D - 14G - KS


S25D - 10G - KS

**IKA®-Werke GmbH & Co. KG**  
Janke & Kunkel-Str. 10  
D-79219 Staufen  
Tel. +49 7633 831-0  
Fax +49 7633 831-98  
sales@ika.de

[www.ika.com](http://www.ika.com)

34 51 6005

Technische Information  
Technical Information  
Informations techniques

DE  
EN  
FR

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**Safety instructions**

- The dispersion tool may heat up during operation. In particular in the processing of hot medium, the tool must be cooled before it can be removed from the drive.
- When in operation the dispersion tool can become hot. If the dispersion tool is not inserted into the drive flange correctly or the medium is hot, it can become extremely hot.
- Don't touch rotating parts during operation.
- No liquid may get into the drive unit of the dispenser.
- Abrasion of the dispersion equipment or the rotating accessories can get into the medium you are working on.
- While working with the dispenser, the user must wear his personal protective equipment, in particular safety goggles.
- Follow the health and safety regulations accident prevention regulations.
- Do not use any damaged dispersion tools e.g. with hairline cracks or splits.
- Note the operating instructions of the dispenser unit.
- When using disposable tube **518/25-ET150**, ensure that the dispersing element has fully engaged in the cover.
- The maximum capacity is 30 ml.
- Only operate in a vertical working position.
- The distance between the stand and the dispersing instrument must not exceed 80 mm.

**Product information**

- The dispersion tool comprises shaft tube, rotor shaft and two bearing shells.
- Dispersion tools made of plastics are not sterile and not free of DNA or RNA.
- The tools are ideal for applications, with which no cross contamination is permissible.
- The tools are ideally appropriate for a single-use. However they can be used, depending upon the application, several times too.
- In this case please consider the cleaning references.
- Use disposable tube **518/25-ET150** for working with the plastic dispersing element in a closed system.
- Attention! Never run the dispersion tool dry. The bearing will be destroyed without cooling by the medium.
- Always ensure that the dispersing element does not run dry in some cases, should spouts develop.
- The fill level of the medium must be between the minimum and maximum marks.
- Medium may leak out of the overflow hole around the maximum mark.
- Dispersion tools are not suitable for continuous operation.
- Dispensing causes the medium to heat up.
- The dispersing element must not be operated at the maximum permitted speed of 25000 rpm for more than one minute.
- Store the dispersing elements in a dry place, in their packaging and at room temperature.
- The solid bodies in samples must not be frozen or too hard.
- Likeover the optimal dispersion duration and rotating frequency must be determined by attempts. Usually a few seconds are sufficient, maximum duration is one minute.
- Longer application times bring no improvements, increase however the sample temperature substantially.
- Dispose of the tools after use according to the regulations. However the laboratory regulations.

**Application instruction**

- The shaft tube are made of polycarbonate (PC), the rotor shafts and bearing shells are made of polyethylenetherton (PEEK). The plastic materials used conform to the FDA (Food and Drug Administration Approval) - rules.
- The parts have a good stability against weak acids, chlorides, hypochlorides and many other chemicals.
- Cleaning, disinfection and sterilization
- The dispersion tools are not sterile, but they can be sterilised before use.
- Method of sterilisation
- Autoclave up to 125 °C with 2 bar  
Ethylene oxides  
Formaldehyde (2%)  
Gamma radiation up to 5 MrAD
- UV-irradiation  
Hot-air sterilisation 160-190 °C
- Biological disinfection
- Alcohol (70%)  
Hypochlorides (5%)  
Formaldehyde (40%)
- Penicill penazone  
Expanded uses of guanidine
- Radioactive decontamination
- Radioazsaun  
Count- Off  
Water / ethyl alcohol / SDS

**Materials and chemical resistance**

The shaft tube are made of polycarbonate (PC), the rotor shafts and bearing shells are made of polyethylenetherton (PEEK). The plastic materials used conform to the FDA (Food and Drug Administration Approval) - rules. The parts have a good stability against weak acids, chlorides, hypochlorides and many other chemicals.

**Cleaning, disinfection and sterilization**

The dispersion tools are not sterile, but they can be sterilised before use.

	<b>Allowed procedures</b>	<b>Not allowed procedures</b>		
Method of sterilisation	Auoclave up to 125 °C with 2 bar Ethylene oxides Formaldehyde (2%) Gamma radiation up to 5 MrAD	UV-irradiation Hot-air sterilisation 160-190 °C		
Biological disinfection	Alcohol (70%) Hypochlorides (5%) Formaldehyde (40%)	Penicill penazone Expanded uses of guanidine		
Radioactive decontamination	Radioazsaun Count- Off Water / ethyl alcohol / SDS			
<b>Technical data</b>				
	<b>518D-10G-KS</b>	<b>518D-14G-KS</b>	<b>525D-10G-KS</b>	<b>525D-14G-KS</b>
Suitable for dispenser	T 18 digital	T 18 digital	T 25 digital	T 25 digital
Working range	10-100	10-500	10-100	10-500
Stator /rotor ØD)	10 / 6,8	14 / 9,8	10 / 6,75	14 / 9,5
Max. allowable speed	rpm	25000		
Min. / max. circumferential speed	mm/s	8,9	12,8	8,5
Min. / max. immersion depth	mm		15 / 85	12
Materials in contact with medium				
Max. working temperature	°C		100	
Sterilization methods		autoclavable		

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