Data File 18-1107-22 AB

Gel filtration

Sephadex LH-20

Sephadex[™] LH-20 is a liquid chromatography medium designed for molecular sizing of natural products such as steroids, terpenoids, lipids, and low molecular weight peptides (up to 35 amino acid residues). Depending on the chosen solvents, this medium can also separate sample components by partition between the stationary and mobile phases.

Sephadex LH-20 is useful for both analytical and industrial scales for the preparation of closely related molecular species. Due to the unique physico-chemical properties of this medium, it can be used either during initial purification prior to polishing by high performance ion exchange or reversed phase chromatography, or as the final polishing step, for example, during the preparation of diastereomers.

Sephadex LH-20 is characterized by:

- Unique chromatographic selectivity due to dual hydrophilic and lipophilic nature of the matrix
- Easily predicted elution behavior based on the chemical structure of the sample (1–3, 14)
- Chemical and physical robustness, see Table 1
- Excellent batch-to-batch reproducibility

Physical and chemical characteristics

Sephadex LH-20 is beaded, cross-linked dextran that has been hydroxypropylated to yield a chromatography medium with both hydrophilic and lipophilic character. Due to its dual character, Sephadex LH-20 swells in water and a number of organic solvents. The chemical structure of the medium is shown in Figure 1. Both the wet particle size and exclusion limit for the medium vary depending on the solvent used for swelling, see Table 2. Table 1 summarizes the general physico-chemical and chromatographic properties of Sephadex LH-20.




Fig 1. Partial structure of Sephadex LH-20. The medium is based on hydroxypropylated dextran that has been cross-linked to yield a polysaccharide network.





Table 1. General physico-chemical properties and chromatographic performance characteristics of Sephadex LH-20

Matrix	Hydroxypropylated, cross-linked dextran	
	5 51 15	
Bead form	Spherical, porous	
Average particle size	70	
Dry In methanol	70 μm 107 μm	
	103 µm	
pH stability	0.17	
Working Cleaning	2-13 2-13	
5		
Chemical stability	Stable in most aqueous and organic eluent systems. Not stable below pH 2 nor to strong oxidizing agents	
Autoclavable	20 minutes at 121°C	
Maximum linear flow rate	700 cm/h	
Recommended linear flow rate	60 cm/h	
Operating temperature	4°C to 40°C	
Sample loading volumes		
Adsorption mode	Depends on resolution required	
Molecular sizing	< 2% total medium volume	
Partition mode	< 1% total medium volume	
Exclusion limit	4000–5000 (depends on solvent)	

Isolation of an HIV-1 reverse transcriptase (HIV-1 RT) inhibitor from *Phyllanthus niruri*

Phyllanthus niruri (P. niruri) has been used for many years as a natural medicine against edema and jaundice (13). Early studies showed isolations of a variety of compounds from this plant but none of the isolates were clearly correlated to any antiviral activities. Then in separate studies (5, 6), extracts of *P. niruri* showed inhibitory effects on an endogenous DNA polymerase from hepatitis B virus (HBV pol), on a reverse transcriptase from avian myeloblastosis virus, and on other DNA dependent polymerases. During a recent study (4), an active compound, identified as repandusinic acid A monosodium salt (RA), was isolated from an aqueous extract which showed inhibition of HIV-1 RT.

RA was found to exert significant inhibitory effects on the production of HIV-1 capsid protein in an *in vitro* culture system. These results and others (7) suggest that the tannin-like structure of RA might fit the binding site of template-primer on the enzyme. Table 3 shows the purification results of RA from *P. niruri*. The structure of the free acid is illustrated in Figure 2.




Fig 2. Structure of free acid form of repandusinic acid A.

 Table 2. Approximate packed medium volumes of Sephadex LH-20 as

 swollen in different solvents

Solvent	Approx. medium volume (ml/g dry powder)
Dimethyl sulfoxide	4.4-4.6
Pyridine	4.2-4.4
Water	4.0-4.4
Dimethylformamide	4.0-4.4
Methanol	3.9-4.1
Saline	3.8-4.2
Ethylene dichloride	3.8-4.1
Chloroform ¹	3.8-4.1
Propanol	3.7-4.0
Ethanol ²	3.6-3.9
Isobutanol	3.6-3.9
Formamide	3.6-3.9
Methylene dichloride	3.6-3.9
Butanol	3.5–3.8
Isopropanol	3.3-3.6
Tetrahydrofuran	3.3-3.6
Dioxane	3.2-3.5
Acetone	2.4-2.6
Acetonitrile ³	2.2-2.4
Carbon tetrachloride ³	1.8-2.2
Benzene ³	1.6-2.0
Ethyl acetate ³	1.6-1.8
Toluene ³	1.5-1.6

¹ Containing 1% ethanol.

² Containing 1% benzene.

³ Solvents giving a medium volume of less than about 2.5 ml/g dry powder are generally not useful.

Table 3. Summary of data for the purification of repandusinic acid A from P. niruri

Purification step	Yield (mg)	ID50 ¹ (µg/ml)	Specific activity (× 10² IU/mg)²	Total activity (× 10 ³ IU) ²
H ₂ O extract	6600	50	4	2640
Methanol insoluble	2500	20	10	2500
Sephadex LH-20, fr. 4–11 ³	247	3.0-3.6	56-67	1616
Cellulose				
Fr. 1	189	7.8	26	484
Fr. 2	24	5.0	40	96
Fr. 3	18	2.4	83	150
Fr. 4	9	3.4	58	52
Fr. 5	14	1.8	111	156
RA (pure substance)	5.9	0.3	668	394

¹ ID50 indicates the effectiveness of inhibitors expressed as concentrations that cause 50% inhibition of HIV-1 RT. Crude HIV-1 RT was used in this experiment.

² IU are arbitrary inhibitory activity units obtained by dividing the total weight of the fraction at each step by the weight of each fraction required to achieve 50% inhibition of [3H]dTTP incorporation into the polymer in the HIV-1 RT assay.

 3 Fractions 4–10 and fraction 11 were combined as both fractions had the inhibitory activity.

Preparative scale separation of 2-acetamidobenzoic acid from 4-acetamidobenzoic acid (1)

The dual character of Sephadex LH-20, hydrophilicity and lipophilicity, provides unique chromatographic selectivity and extremely high resolution of closely related molecular species. Figure 3 shows a preparative scale separation of a mixture containing both 2-acetamidobenzoic acid and 4-acetamidobenzoic acid. The two molecules differ only by the position of the acetamide function on the benzene ring. The separation was possible due to the unique selectivity of the Sephadex LH-20.




Fig 3. Separation of 2- and 4-acetamidobenzoic acid on Sephadex LH-20.

Isolation of BE-23372M, a novel protein tyrosine kinase inhibitor

BE-23372M was isolated (8) from the culture broth of a fungus and the producing strain, F23372, was identified as *Rhizoctonia solani*. The protein kinase inhibitory activity was purified by a sequence of solvent extractions and then by silica and Sephadex LH-20 column chromatographies.

The particular compound, BE-23372M, was discovered while screening for EGF protein tyrosine kinase inhibitors, which could prove useful as selective drugs for treating cancers. BE-23372M showed strongly specific inhibitory activity to EGF receptor kinase with IC_{s0} values of 0.02 and 0.03 μ M on two substrates, see Table 4.

The compound also inhibited the growth of A431 human epidermoid carcinoma and MKN-7 human stomach cancer cell lines.

Sephadex LH-20 (2.5 \times 50 cm column) was used to prepare 5.7 mg of BE-23372M in the form of a reddish orange solid substance. The physico-chemical properties, structure elucidation, and synthetic studies of this compound are reported in references 9 and 10. The specificity of inhibition of BE-23372M is illustrated through the data in Table 4.

Table 4. The inhibitory effect of BE-23372M on protein kinases

Protein kinase	Substrate	IC ₅₀ (μΜ)
EGF receptor kinase	Poly(Glu: Tyr)	0.02
EGF receptor kinase	RR-SRC	0.03
Protein kinase C	Histone (Lys rich)	4.5
cAMP-dependent protein kinase	Histone H2B	> 20

In vitro test to measure physical and chemical changes in doxorubicin (DXR) in liposomal formulations

Two different liposomal doxorubicin (DXR) preparations were solubilized in ethanol and then run chromatographically on a column packed with Sephadex LH-20 in ethanol. The DXR and phospholipid were determined by A₄₈₀ and modified Bartlett assay, respectively (11). Figure 4A shows the elution profile of L-DXR that had been stored at 4°C. All the phospholipids eluted in the void volume probably as micelles. The volume of elution for DXR suggested a molecular weight of 550, which is very similar to its actual molecular weight of 580. Figure 4B shows a comparison of the elution profile of free DXR and L-DXR stored in lyophilized form at 60°C.

The earlier elution of the DXR, which had been stored at elevated temperature (60°C) and in liposomal form, suggests a molecular weight of 1000 to 1100. The appearance of the higher molecular weight of DXR could be the result of either physical aggregation or chemical oligomerization. It is expected that the molecular weight change of DXR might restrict its biological availability. Therefore, this assay might be useful to assess the therapeutic efficacy of liposomal DXR by *in vitro* testing.

The resolving power of Sephadex LH-20 was further demonstrated by its ability to separate DXR and its polar degradation product, which was obtained after 12 hours of incubation in plasma at 37°C (not shown). The elution volume of this compound indicated a molecular weight of approximately 450, which together with spectral data suggested 7-deoxydemethyl DXR or its aglycone.




Fig 4A. Sephadex LH-20 (1 × 20 cm) elution profiles of L-DXR clinical batch LT4C stored at 4°C and phospholipid by modified Bartlett assay (•...•...•) and $A_{\tiny ueo}$ (o...o...o).




Fig 4B. Sephadex LH-20 (1 \times 20 cm) elution profile of L-DXR clinical batch LT60C stored at 60°C (... \square ... \square ...) and free DXR (... \blacksquare ...). Results of the Bartlett assay are not shown.

Preparative separation of an epimeric mixture of budesonide

Attempts to resolve epimeric mixtures of steroids similar to budesonide by fractional crystallization have met with limited success. The development of a chromatographic separation method on Sephadex LH-20 proved to provide a more effective preparative method (12). As synthesized, the glucocortoid, budesonide, is a mixture of 1:1 of the C-22 epimer, see Figure 5. During the evaluation of budesonide as an inflammatory, it was desirable to evaluate the chemical and pharmacological properties of the two epimers. Sephadex LH-20 was packed into a preparative column (6.3 × 75 cm) and different loadings of the epimeric mixture of budesonide were run to see the trade-off between resolution and loading, see Figure 6. When restricting the chromatography to a single cycle, sample loading at 1500 mg proved to be a practical preparative load.




Fig 5. Structure of budesonide C-22 epimer.

Medium:	Sephadex LH-20
Mobile phase:	Heptane-chloroform-ethanol, 20:20:1
Solute:	Budesonide 21-acetate (XVII)
Sample volumes:	A) 500 mg, B) 1500 mg, C) 3000 mg, and D) 5000 mg





Fig 6. Effect of loading with high sample mass in preparative gel chromatography. Amount of sample/fraction plotted against elution volume.

Isolation of the cyclic peptide, RP71955, an inhibitor of HIV-1 aspartyl protease

RP71955 (15), isolated from *Streptomyces* strain SP9440, has been shown by Helynck *et al.* (1993) to be an effective inhibitor of the HIV-1 protease. After centrifugation of 440 liters from an 800-liter fermentation, the mycelium was extracted with acetone and again centrifuged. The supernatant was then further reduced to 60 liters. Some of the feedstock was then purified through a number of chromatography steps including Sephadex LH-20.

The basic purification scheme is as follows:



References

- Henke, H. Präparative Gelchromatographie an Sephadex LH-20, Salz und Druck: Alpha Druck GmbH, 63773 Goldbach, p. 612 (1994). (For information on an English version contact your local GE Healthcare representative.)
- Sund, R. B. Chromatographic separation of danthron metabolites, Acta Pharm. Nordica 2, 89–98 (1989).
- Sephadex LH-20 chromatography in organic solvents, Pharmacia Fine Chemicals AB, p. 10 (1978).
- Ogata, T. et al. HIV-1 reverse transcriptase inhibitor from Phyllanthus niruri. AIDS Res Hum Retroviruses 8, 1937–1944 (1992).
- Venkateswaren, P. S. et al. Effects of an extract from Phyllanthus niruri on hepatitis B and woodchuck hepatitis viruses: In vitro and in vivo studies, Proc Natl Acad Sci USA 84, 274–278 (1987).
- Yanagi, M. et al. Proceedings of the Japanese Cancer Association, 48th Annual Meeting, p. 99 (1989).
- Kakiuchi, N. et al. Inhibitory effect of tannins on reverse transcriptase from RNA tumor virus. J. Nat. Prod. 48, 614–621 (1985).
- Okabe, T. *et al.* BE-23372M, a novel protein tyrosine kinase inhibitor. I. Producing organism, fermentation, isolation and biological activities. J. Antibiot. (Tokyo) 47, 289–293 (1994).
- Tanaka, S. et al. BE-23372M, a novel protein tyrosine kinase inhibitor. II. Physicochemical properties and structure elucidation. J. Antibiotics 47, 294–296 (1994).
- BE-23372M, a novel protein tyrosine kinase inhibitor. III. Synthesis. J. Antibiotics 47, 297–300 (1994).
- 11. Amselem, S. *et al. In vitro* tests to predict *in vivo* performance of liposomal dosage forms. *Chem. Phys. Lipids* **64**, 219–237 (1993).
- Thalén, A. and Nylander, B. Epimers of budesonide and related corticosteroids: Preparative resolution by chromatography on Sephadex LH-20. Acta Pharm Suec 19, 247–266 (1982).
- 13. Chopra, R. N. *et al. Glossary of Indian Medicinal Plants*, Council of Scientific and Industrial Research, New Delhi, p. 191 (1956).
- Henke, H. Selektive präparative gelchromatographische Trennung niedermolekularer Verbindungen an Sephadex LH-20. J. Chromatogr. 254, 296–308 (1983).
- Helynck, G. *et al.* Isolation of RP71955, a new anti-HIV-1 peptide secondary metabolite. *J. Antibiotics* 46, 1756–1757 (1993).

Additional applications of Sephadex LH-20

The following references are a few more publications where Sephadex LH-20 has been employed in the purification strategy. The articles are grouped according to the type of molecule and/or area of research.

Peptides and non-peptidic antibiotics

Liu, S. Y. *et al.* Studies on the isolation of endogenous sleep factors from Tupaia Belangeri Chinesis (TBC) after sleep deprivation. *Sheng Li Hsueh Pao* **46**, 83–89 (1994).

Asano, T. et al. Novel retrovirus protease inhibitors, RPI-856 A, B, C, and D, produced by Streptomyces sp. AL-322. J. Antibiot. (Tokyo) 47, 557–565 (1994).

Ueno, M. *et al.* IC 101, extracellular matrix antagonist produced by *Streptomyces* sp. MJ 202-72F3. Production, isolation, structure determination and biological activity. *J. Antibiot. (Tokyo)* **46**, 1658–1665 (1993).

Ueno, M. *et al.* Del aminomycins, novel nonpeptide extracellular matrix receptor antagonist and a new class of immunomodulator. I. Taxonomy, fermentation, isolation and biological activity. *J. Antibiot. (Tokyo)* **46**, 719–727 (1993).

Tabata, N. *et al.* Diolmycins, new anticoccidial agents produced by *Streptomyces* sp. I. Production, isolation and physico-chemical and biological properties. *J. Antibiot.* (*Tokyo*) **46**, 756–761 (1993).

Lipids/surfactants

Miller, P. H. *et al.* New method for the isolation of polyglycosyleceramides from human erythrocyt membranes. *Biochim. Biophys. Acta* **1168**, 330–339 (1993).

Rider, E. D. *et al.* Different ventilation strategies alter surfactant responses in preterm rabbits. *J. Appl. Physiol.* **73**, 2089–2096 (1992).

Fidelio, G. D. et al. Molecular parameters of gangliosides in monolayers: comparative evaluation of suitable purification procedures. J. Biochem. (Tokyo) **110**, 12–16 (1991).

Crick, D. C. *et al.* Characterization and localization of a long-chain isoprenyltransferase activity in porcine brain: proposed role in the biosynthesis of dolichyl phosphate. *J. Neurochem.* **57**, 1354–1362 (1991).

Lam, C. *et al.* Immunostimulatory, but not antiendotoxin, activity of lipid X is due to small amounts of contaminating N, O-acylated disaccharide-1-phosphate: *in vitro* and *in vivo* reevalutation of the biological activity of synthetic lipid X. *Infect. Immun.* **59**, 2351–2358 (1991).

Tannins/plant extracts

Inaoka, Y. *et al.* Studies on active substances in herbs used for hair treatment. I. Effects of herb extracts on hair growth and isolation of an active substance from Polyporus umbellatus F. *Chem. Pharm. Bull. (Tokyo)* **42**, 530–533 (1994).

Tanaka, T. *et al.* Sensitive determination of zearalenone and alphazearalenol in barley and Job's tears by liquid chromatography with fluorescence detection. *J. AOAC Int.* **76**, 1006–1009 (1993).

Nakahara, K. *et al.* Inhibitory effect of oolong tea polyphenols on glycosyltransferases of mutans streptococci. *Appl. Environ. Microbiol.* **59**, 968–973 (1993).

Tang, H. R. *et al.* Study on the composition and structure of commercial chestnut tanning agent. *Basic Life Sci.* **59**, 221–243 (1992).

Kumar, R. Proposis cineraria leaf tannins: their inhibitory effect upon ruminal cellulase and the recovery inhibition by polyethylene glycol-4000. *Basic Life Sci.* **59**, 699–704 (1992).

Steroids

Nahoul, K. Plasma 17-hydroxyprogesterone determination with two commercial immunoassays. J. Steroid. Biochem. Mol. Biol. **50**, 197–203 (1994).

Remer, T. *et al.* Measurement of urinary androgen sulfates without previous hydrolysis: a tool to investigate adrenarche. Determination of total 17-ketosteroid sulfates. *Steroids* **59**, 16–21 (1994).

Renoir, J. M. et al. Effects of immunosuppressants FK506 and rapamycin on the heteroligomeric form of the progesterone receptor. J. Steroid. Biochem. Mol. Biol. **48**, 101–110 (1994).

Minoto, K. *et al.* Biliary metabolites of levonorgestrel in rats. *Yakugaku Zasshi* **113**, 781–791 (1993).

Dorr, H. G. and Sippel, W. G. Prenatal dexamethasone treatment in pregnancies at risk for congential adrenal hyperplasia due to 21-hydroxylase deficiency: effect on midgestrational amniotic fluid steroid levels. *J. Clin. Endocrinol. Metab.* **76**, 117–120 (1993).

Toxins/pollutants

Lee, H. *et al.* Mutagenicity of airborne particulates from combustion of electric cables in a waste metal retrieval area. *Mutat. Res.* **324**, 77–84 (1994).

Lee, H. *et al.* Correlation between meterological conditions and mutagenicity of airborne particulate samples in a tropical monsoon climate area from Kaohsiung City, Taiwan. *Environ. Mol. Mutagen.* **23**, 200–207 (1994).

Craig, M. *et al.* Identification and characeterization of hydrophobic microcystins in Canadian freshwater cyanobacteria. *Toxicon* **31**, 1541–1549 (1993).

Singh, Y. N. *et al.* Studies on the muscle-paralyzing components of the juice of the banana plant. *Arch. Int. Pharmacodyn. Ther.* **324**, 105–113 (1993).

Sayato, Y. *et al.* Identification of polycyclic aromatic hydrocarbons in mutagenic adsorbates to a copper-phthalocyanine derivative recovered from municipal river water. *Mutat. Res.* **300**, 207–213 (1993).

Ordering information

Products	Quantity	Code no.
Sephadex LH-20	25 g	17-0090-10
Sephadex LH-20	100 g	17-0090-01
Sephadex LH-20	500 g	17-0090-02
Sephadex LH-20	5 kg	17-0090-03

Related products

SR 25/45 column	1	19-0879-01
SR 25/100 column	1	19-0880-01

www.gelifesciences.com/protein-purification

GE Healthcare Bio-Sciences AB Björkgatan 30 751 84 Uppsala Sweden GE, imagination at work, and GE Monogram are trademarks of General Electric Company.

Drop design and Sephadex are trademarks of GE Healthcare companies.

All third party trademarks are the property of their respective owners. © 1995-2007 General Electric Company – All rights reserved. First published March 1995.

All goods and services are sold subject to the terms and conditions of sole of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

GE Healthcare Europe GmbH Munzinger Strasse 5, D-79111 Freiburg, Germany

GE Healthcare UK Limited Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK

GE Healthcare Bio-Sciences Corp. 800 Centennial Avenue, P.O. Box 1327, Piscataway, NJ 08855-1327, USA

GE Healthcare Bio-Sciences KK Sanken Bldg. 3-25-1, Hyakunincho Shinjuku-ku Tokyo 169-0073, Japan

Asia Pacific Tel: +85 65 62751830 Fax: +85 65 62751829 • Australasia Tel: +61 2 8820 8299 Fax: +61 2 8820 8209 • Austria Tel: 01/57606 1613 Fax: 01/57606 1614 • Belgium Tel: 0800 73 890 Fax: 02 416 8206 • Canada Tel: 1800 463 5800 Fax: 1800 567 1008 • Central, East, & South East Europe Tel: +43 1 972 720 Fax: +43 1 972 722 750 • Denmark Tel: +45 70 25 24 50 Fax: +45 45 16 2424 • Eire Tel: 1800 709992 Fax: +44 1494 542010 • Finland & Baltics Tel: +358 9 512 3940 Fax: +35 19 512 39439 • France Tel: 01 69 35 567 00 Fax: 01 69 41 98 77 • Germany Tel: 0800 9080 712 • Greater China Tel: +852 2100 6300 Fax: +852 2100 6338 • Italy Tel: 02 26001 3207 Fax: 02 26001 399 • Japan Tel: 813 5331 9370 • Korea Tel: 82 2 6201 3700 Fax: +25 11 3933 7300 Fax: +55 11 3933 7304 • Middle East & Africa Tel: +30 0 687 Fax: +30 210 96 00 687 ex: +43 196 77 Fax: +47 815 65 777 Fax

