

easySpiral®

Automatic plater

Ref 412 000



Standards:



Manufactured under:



easySpiral® is an automatic plater for Petri dishes plating a liquid sample from 300 to 1.3×10^5 CFU/mL on one Petri dish without prior dilution.

TECHNICAL SPECIFICATIONS

- All stainless steel
- Controlled by microprocessor
- Fast rotating arm (Patented)
- Multi-use, autoclavable stylus
- Plating modes: exponential and circle
- Petri dish diameter: 90 mm
- Test key for checking plating volumes (volumetric volumes)
- Counting range: from 300 to 1.3×10^5 CFU/mL
- Preset dispensed volumes: 50 or 100 μ L
- Volumes accuracy: $\pm 0,5$ % (in normal mode)
- Syringe capacity: 1000 μ L
- Plating time: 8 seconds
- Plating cycle: 1 disinfection + 1 sample intake + 1 plating: 25 seconds
- Exponential mode: 4 dilutions logs on $\varnothing$ 90 mm Petri dish
- Successive plating capacity with the same sample: up to 20 Petri dishes (50 μ L)
- Sample intake possible in Eppendorf™ 0.5/1.5/ 2 mL beaker or Dilucup™
- Cleaning and disinfection of the stylus by "Overflow" technology. Pression: up to 8 bars (Patented)
- Disinfection process
- Approved disinfectants: bleach 1 % active chlorine or ethanol at 70 % or a solution of hydrogen peroxide H_2O_2 (1 to 5 %) and peracetic acid PAA (0.05 to 1 %)
- Security detection in case of lack of disinfection liquid
- Autonomy in disinfectant (in normal mode): 1000 cycles (2 L bottles)
- Autonomy in distilled water: 600 cycles (2 L bottles)
- Voltage - Frequency: 100-240V ~ 50-60Hz
- Fuses: T1.6A L 250V 5x20 mm
- Power: 45 W
- International machine marking : IP10
- Dimensions (w x d x h): 39,4 x 42,1 x 28,9 cm, weight: 16.3 kg
- Box (w x d x h): 60 x 50 x 55 cm, weight: 24.6 kg

- Warranty: 1 year further to warranty card recording
- Availability of spare parts: 10 years
- In compliance with standards ISO 7218, ISO 4833-2, FDA BAM (Bacteriological Analytical Manual), AOAC 977.27
- Manufactured under CE, UKCA, RoHS, WEEE
- Designed and made in France

DELIVERED WITH

- 1 turntable for Ø 90 mm Petri dishes (by default on the unit)
- 1 box with:
 - 1 blue dye for testing (Ref 415 530)
 - 1 syringe with silicone tubing
 - 1 Spiral® counting grid for Ø 90 mm Petri dishes (Ref 413 014)
 - 1 container stand for testing
 - 1 set of 2 filters (Ref 413 009)
- 1000 DB50 beakers (Ref 415 100)
- 1 EnzyClear® detergent (Ref 413 029)
- 4 connection kit for GL45 bottles (Ref 413 003)
- 1 double connection kit for GL45 bottles (Ref 414 004)
- 1 stylus (Ref 414 001)
- 1 power cable
- 1 user's manual
- 1 EnzyClear® quick user guide

OPTIONS/ACCESSORIES

- Ref 413 001 Housing protection for platers: Transparent housing for working under laminar flow
- Ref 413 002 2 L Bottle (GL45): 2 L GL45 glass bottle
- Ref 413 003 Connection kit for GL45 bottles: Autoclavable connection set for GL45 bottles
- Ref 413 004 Syringe (1 mL): Syringe (1000 µL capacity)
- Ref 413 014 Spiral® (90 mm) counting grid: Counting grid for Ø 90 mm Petri dishes
- Ref 413 016 Support for Eppendorf™ 0.5 mL: Support for 0.5 mL Eppendorf™ sample intake
- Ref 413 017 Support for Eppendorf™ 1.5 and 2 mL: Support for 1.5 and 2 mL Eppendorf™ sample intake
- Ref 413 018 Support for Dilucup™: Support sample intake for Dilucup™
- Ref 413 019 Adjustable plate for Petri dish (90-150 mm): Multiadaptable support for Petri dishes of 90 to 150 mm
- Ref 413 029 EnzyClear®: Spiral® plater recommended detergent for cleaning and pre-disinfection
- Ref 414 001 Stylus (1 mL): 1 mL stylus
- Ref 414 003 Maintenance kit for platers (1 mL): Maintenance kit (1 stylus + 1 syringe + 1 x tubings (Ref 413 006), 1 grease tube, 10 pump filters, 2 female quick connectors)
- Ref 415 100 Beakers (DB50): 1000 beakers for Spiral® plater (5 mL)
- Ref 513 023 PTFE filters 0.2 µm (Ø 25 mm): Set of 5 autoclavable filters (0.2 µm)

TUBINGS

- Ref 413 006 Tubings: Silicone tubings with filter for cleaning pump
- Ref 413 009 Filters for rinsing bath: Set of 2 filters for rinsing bath
- Ref 413 028 Tube weights (Ø 3 mm): Set of 5 tube weights (inner Ø 3 mm)
- Ref 414 004 Double connection kit for GL45 bottles: Double connection kit for GL45 bottle to be used on the same bottle
- Ref 415 530 Blue dye for platers: Blue dye for Spiral® tests (125 mL)