


Enzystel is a triple-enzymatic detergent for the manual cleaning of endoscopes and surgical instruments including soaking in ultrasonic baths. It removes and dissolves proteinaceous and organic matter such as blood, saliva and mucous from instruments quickly and effectively.



### HOW TO USE FOR MANUAL CLEANING


### HOW TO USE FOR AUTOMATED SYSTEMS



**1** Use one sachet or one 20ml bottle dose in 1 litre of water for instrument soaking (1 part Enzystel in 50 parts water).



**2** Soak the instrument for 10 minutes and rinse thoroughly.  
*Note: Solution should be changed weekly or when it becomes cloudy/soiled.*



**1** Use one sachet in 2 litres of water or one 10ml bottle dose in 1 litre of water for a flushing solution (1 part Enzystel in 100 parts water).



**2** Use for 5 minutes or as per manufacturer's recommended machine times.  
*Note: Solution should be changed weekly or when it becomes cloudy/soiled.*

### PRODUCT OPTIONS

# Enzystel™

Enzystel 1L  
MD315

Enzystel 50x20ml Sachets  
MD319



Enzystel is CE marked as a Class I Medical Device in accordance with the European Medical Devices Directive 93/42/EEC and the 2007/47/EC amendments thereto.

**Tristel™**

Created by: Tristel Solutions Limited, Lynx Business Park, Cambs, UK, CB8 7NY  
T +44 (0) 1638 721500 - E mail@tristel.com - W www.tristel.com

Copyright © Tristel Solutions  
Mkt-Fac-406-2  
February 2018

Page 1 of 1