

Hamilton 700 Series Microlitre Syringes

HAMILTON

Precision glass barrel syringes and stainless steel needles. The general purpose 700 series syringes are liquid tight for precision dispensing. Ideal for sample injection with gas, liquid and thin layer chromatography.



Fixed Needle, Sharp

Code	Type	Capacity, µL	Price
SYR7000	75	5	£44.30
SYR7002	701	10	£32.60
SYR7006	705	50	£53.30
SYR7008	710	100	£63.40
SYR7010	725	250	£68.00
SYR7012	750	500	£77.00
SYR7090	1701	10	£46.60
SYR7092	1702	25	£55.10
SYR7050	1705	50	£54.60
SYR7200	1725	250	£63.40
SYR7202	1750	500	£63.60



Removable Needle, Sharp

Code	Type	Capacity, µL	Price
SYR7020	75	5	£61.00
SYR7022	701	10	£58.80
SYR7024	702	25	£65.00
SYR7026	705	50	£67.00
SYR7028	710	100	£73.40
SYR7030	725	250	£87.70
SYR7032	750	500	£98.80
SYR7094	1701	10	£62.10
SYR7095	1702	25	£67.00
SYR7060	1705	50	£67.00
SYR7204	1725	250	£71.60
SYR7206	1750	500	£71.60

Hamilton 1700 Series Microlitre Syringes

HAMILTON

The 1700 series syringes are the gastight version of the original 700 series syringes. Other 1700 series syringes and accessories are available, please enquire for further information.



Fixed Needle, Sharp

Code	Type	Capacity, µL	Price
SYR7090	1701	10	£46.60
SYR7092	1702	25	£55.10
SYR7050	1705	50	£54.60
SYR7200	1725	250	£63.40
SYR7202	1750	500	£63.60