

Hamilton 1700 Series Microlitre Syringes Continued



Removable Needle, Sharp

Code	Type	Capacity, µL	Price
SYR7094	1701	10	£62.10
SYR7095	1702	25	£67.00
SYR7060	1705	50	£67.00
SYR7204	1725	250	£71.60
SYR7206	1750	500	£71.60



Luer Tip

Code	Type	Capacity, µL	Price
SYR7096	1701	10	£77.50
SYR7098	1702	25	£60.50
SYR7070	1705	50	£62.90
SYR7072	1710	100	£60.50
SYR7208	1725	250	£65.00
SYR7074	1750	500	£62.10



Hamilton 7000 Series Microlitre Syringes



The 7000 series syringes are designed on the dead-volume-free principal of positive displacement. A fine tungsten plunger extends through the whole length of the syringe to the needle end.

Fixed Needle, Sharp

Code	Type	Capacity, µL	Price
SYR7220	7101	1	£81.70
SYR7222	7102	2	£112.00
SYR7224	7105	5	£95.90



Fixed Needle, Blunt

Code	Type	Capacity, µL	Price
SYR7226	7101	1	£83.50
SYR7228	7102	2	£81.70
SYR7230	7105	5	£104.00

