



Hanna Checktemp® 1 Pocket Thermometer



A compact food probe with a stainless steel sensor on a 1m flexible cable. Highly accurate it is designed to take fast response temperature measurements in semi-solids and liquids. The flexible silicone cable also makes this meter ideal for taking readings in hard to reach locations. Featuring a water resistant casing and battery compartment both test to IP65, the latest model is ideal for use in food preparation and production including steamy kitchen environments and chilled storage.

- Auto shut-off: prolongs the battery life; also selectable from 8min, 60 min or off
- Tested to and compliant with the European standard for food thermometers EN13485
- Automatic CAL Check™ provides automatic calibration check every time Checktemp® 1 is switched on
- Easy read digital display, whatever the viewing angle
- Monitors the temperature sensor and alerts the user if it is faulty by displaying a clear ERROR message on the display
- Supplied with desktop stand, batteries and instructions

Code	Alt Ref	Range, °C	Accuracy, °C	Resolution	Dims, w x d x h, mm	Price
THE4502	HI-98509	-50 to 150	±0.2 (-30 to 120)	0.1	59 x 17 x 107	£71.00



Traceable® Fridge/Freezer Thermometer

Accurately monitor temperatures in freezers, water baths, heating blocks, incubators, and refrigerators with this enclosed temperature-buffered sensor. Featuring a triple display to simultaneously show min, max and current temperatures.

- Bottle insulates sensor from rapid temperature changes when refrigerator door is opened
- Ten-foot micro-cable permits refrigerator doors to close
- Min/max monitors high and low readings for any time period
- Programmable alarm sounds continuously when temperature exceeds set points
- High-impact, chemical-resistant ABS plastic case
- Supplied with individually-numbered Traceable® Certificate which assures accuracy from an ISO/IEC 17025:2005 (1750.01) calibration laboratory accredited by A2LA

Code	Temp Range, °C	Accuracy, °C	Resolution, °C	Dims, w x d x h, mm	Price
THE8100	-50 to 70	±1	1	70 x 19 x 108	£102.00



Alarm Thermometer



Modern, compact design combined with ease of use making this a cost-effective thermometer for liquid and semi-solid sampling. The stainless steel probe is connected to the unit with a 1 metre cable and can be clipped to the instrument body for single-handed use.

- Sampling rate of one reading per second
- Data hold function retains displayed value
- Min/max memory
- Programmable alarm for min/max values
- Auto power off after 1 hour, low battery indicator
- 100 x 4mm Ø stainless steel pointed tip probe with handle and 1m cable
- 30mm LCD display with function indicators
- Impact resistant ABS case incorporates a desk stand and compartment for the sensor cable when not in use
- Supplied with a probe protector cap & battery

Code	Temp Range, °C	Accuracy, °C	Resolution, °C	Dims, w x d x h, mm	Price
THE1502	-50 to 150	±1 (-30 to 150), otherwise ±2	0.1	60 x 20 x 120	£27.50