



Xylem Ebro EBI 20 Temperature Loggers

-ebro-
a xylem brand

With a memory capacity of 40,000 measurements the easy to use EBI 20 data loggers are suitable for the continuous documentation and monitoring of temperature. All EBI 20 data loggers are delivered with a factory calibration certificate and a user replaceable battery. The data loggers are particularly attractive because of their excellent price-performance ratio.

- Automatic recording of temperature data
- No network connection required
- Programming and evaluation with PC
- Waterproof
- LED alarm
- Min/max values on display
- Battery life approximately 1 year

Code	Description	Measuring Range, °C	Accuracy, °C	Price
THE1838	EBI 20-T1 Temperature logger with internal sensor	-30 to 70	±0.5 (-20 to 40), ±0.8 (40.1 to 70)	£58.00
THE1836	EBI 20-T1-Set Temperature logger set with interface and software	-30 to 70	±0.5 (-20 to 40), ±0.8 (40.1 to 70)	£158.00
THE1840	EBI 20-TE1 Temperature logger with external sensor	-30 to 70	±0.5 (-20 to 40), ±0.8 (40.1 to 70)	£120.00
THE1068	EBI 20-TE1-Set Temperature logger set with interface and software	-30 to 70	±0.5 (-20 to 40), ±0.8 (40.1 to 70)	£219.00
THE1844	EBI 20-TH1 Humidity and temperature logger	-30 to 70, 0 to 100% RH	±0.5 (-20 to 40), ±0.8 remaining, ±0.3% RH	£95.00
THE1074	EBI 20-TH1-Set Humidity and temperature logger set with interface and software	-30 to 70, 0 to 100% RH	±0.5 (-20 to 40), ±0.8 remaining, ±0.3% RH	£204.00
THE1842	EBI 20-TF Food logger with external probe	0 to 100	±0.5 (50 to 100), ±1 remaining	£136.00



SLS Select Stopwatch



- Impact-resistant ABS case
- Standard quartz time
- Split-time facility
- Full clock functions with alarm
- Clear display with large digit readout
- Supplied with neck cord and battery (LR43)

Code	Timer	Resolution	Price
TIM1030	23 hrs, 59 mins, 59 secs (up)	1/100 secs for 30 mins, 1 sec thereafter	£8.00



SLS Select Count Up/Down Timer and Clock



- Digital timer with 3 functions: countdown stopwatch and clock
- Countdown timer settable to nearest second up to 24 hours, with count-up from zero
- Alarm sounds for 1 minute at zero
- 80dB sound
- With pocket clip, bench stand and mounting magnet

Code	Timer	Price
TIM0250	23 hrs, 59 mins, 59 secs (up/down)	£24.30

