

Corning® Falcon® Multiwell Plates Continued



Code	No. of Wells	Bottom	Inner Pack	Pack	Price
353046	6	Flat	1	50	£139.00
353934	6	Flat	10	60	£165.00
353224	6	U	6	36	£100.00
353043	12	Flat	1	50	£181.00
353225	12	Flat	6	36	£114.00
353047	24	Flat	1	50	£146.00
353226	24	Flat	6	36	£103.00
353935	24	Flat	10	60	£171.00
353078	48	Flat	1	50	£200.00
353230	48	Flat	6	36	£145.00
353072	96	Flat	1	50	£165.00
353075	96	Flat	5	50	£149.00
353376	96	Flat	8	32	£172.00
353936	96	Flat	14	84	£300.00
353916	96	Flat	25	100	£356.00
353077	96	U	1	50	£160.00
353227	96	U	5	50	£160.00

Corning® Falcon® Dishes



- Flat, optically clear polystyrene surfaces for distortion-free microscopic visualisation of cells
- Uniform surface chemistry created by vacuum-gas plasma treatment promotes cell attachment
- Improved handling of small dishes with the unique Easy-Grip design
- Nonpyrogenic
- Sterilised by gamma irradiation
- Packaged in peel-open, medical-style bags
- Crystal-grade virgin polystyrene
- Lids designed for optimal gas exchange
- Stacking rings allow for easier stacking and handling
- Vacuum-gas plasma treatment permanently and consistently modifies the cell growth surface



Code	Description	Dims, Ø x h, mm	Inner Pack	Pack	Price
353001	Easy-Grip	35 x 10	20	500	£195.00
353002	Standard	60 x 15	20	500	£228.00
353003	Standard	100 x 20	20	200	£163.00
353025	20mm grid	150 x 25	10	100	£204.00
353037	Centre-well organ	60 x 15	20	500	£513.00

Liquid handling... Down to the last drop!

From pipettes to dispensers, tips to pipetting stations,
SLS is your No.1 source for liquid handling.

See pages 745–815



**SCIENTIFIC
LABORATORY
SUPPLIES**

