




Corning® Falcon® *in vitro* Fertilisation Plasticware



- Virgin crystalline polystyrene tested for USP Class VI cytotoxicity
- Ideal for applications like oocyte preparation as well as fertilisation, culture, and transfer
- Sterile, non-pyrogenic and non-embryotoxic (except lids)
- Comply with the College of American Pathologists (CAP) and American Fertility Society (AFS) recommended standards for IVF labs
- Flat dish bottoms for distortion-free optics and temperature stability and control
- Dish bottoms raised slightly to reduce abrasions and ensure optical clarity
- Bevel around interior perimeter to prevent settling in the corner (four well plate/dish)
- Lids designed for aseptic manipulation and consistent venting to maintain humidification
- Lids designed for consistent gas exchange except low-wall, tight fit lid which limits gas exchange
- Peel-open, medical-style packaging (four well plate/dish, low-wall dish)
- Multi-unit bags have reseal tabs (60mm dishes)

Code	Description	Tissue Culture Treated	Well Area, cm ²	Well Volume, mL	Inner Pack	Pack	Price
353652	60 x 15mm round IVF dish	Yes	21.29	23.0	20	500	£403.00
353653	60 x 15mm centre well IVF dish	Yes	2.45	2.5	20	500	£632.00
353655	50 x 9mm low wall IVF dish	No	18.3	—	4	100	£94.00
353654	17 x 11mm 4-well IVF plate	Yes	1.39	1.8	1	100	£500.00
353671	17 x 11mm 4-well IVF dish	Yes	1.39	1.8	1	100	£501.00



Corning® Falcon® Culture Slides



- Culture cells and then analyse them on a glass microscope slide
- Cells are grown in a plastic chamber affixed to a specially prepared glass microscope slide
- Cells can be fixed and stained in place without disruption of the cell monolayer
- Chamber is easily and safely removed with a supplied disposable safety removal tool
- Innovative sealing design eliminates leakage
- 1.2mm bevelled-edge slide, 25 x 75mm, soda-lime glass
- Pressure sensitive, biocompatible, non-migrating, acrylic-adhesive gasket
- Blue hydrophobic border defines cell culture areas
- Polystyrene vessel, lid and tool
- Wells numbered for easy identification
- Sterile
- Trays designed for use in incubator
- Supplied in inner packs of 12

Code	No. of Wells	Growth Area, cm ²	Well Volume, mL	Working Volume, mL	Pack	Price
354104	4	1.7	1.5 - 1.7	0.7 - 1.25	96	£571.00
354114	4	1.7	1.5 - 1.7	0.7 - 1.25	24	£170.00
354108	8	0.7	0.7 - 0.75	0.3 - 0.5	96	£571.00
354118	8	0.7	0.7 - 0.75	0.3 - 0.5	24	£170.00