

Corning® Microplate Surface Treatments Continued

CORNING




96 Well Non-Binding Surface (NBS™) Plates

- Non sterile
- Without lids

Supplied in packs of 100 (inner packs of 25)

Code	Colour	Well Volume, µL	Bottom	Price
3642	White	190	Flat	£979.00
3686	Black	190	Flat	£1130.00
3881	Black/clear	205	Flat	£1175.00
3884	White/clear	205	Flat	£1175.00
3605	White	330	Round	£908.00
3641	Clear	360	Flat	£656.00
3600	White	360	Flat	£1111.00
3650	Black	360	Flat	£1015.00
3604	White/clear	360	Flat	£1339.00
3651	Black/clear	360	Flat	£1015.00



Corning® 96 Well Polypropylene Storage Blocks

CORNING

- Uniform skirt heights for greater robotic gripping surface
- Solvent resistant polypropylene provides compatibility with many common organic solvents (e.g., DMSO, ethanol, methanol)
- RNase/DNase free
- Clear

Code	Shape	Well Volume, mL	Bottom	Sterility	Inner Pack	Pack	Price
3958	Round well	1	Round	Sterile	5	25	£176.00
3956	Round well	0.5	V	Sterile	10	50	£362.00
3957	Round well	0.5	V	Non sterile	10	100	£683.00
3960	Square well	2	V	Sterile	5	25	£293.00
3961	Square well	2	V	Non sterile	5	100	£999.00



Corning® 384 Well Medium Bind (Non-Treated) Plates

CORNING

- Flat bottomed, without lid
- Non sterile

Code	Colour	Description	Well Volume, µL	Inner Pack	Pack	Price
3821	Black	Low volume	50	10	50	£448.00
3540	Black/clear	Low volume	50	10	50	£552.00
3702	Clear	—	112	25	100	£636.00
3572	White	Low flange	112	10	50	£396.00
3573	Black	Low flange	112	10	50	£356.00
3762	Black/clear	Low flange	112	25	100	£1323.00