


Corning® 384 Well High Binding Surface Plates



- Flat bottomed, without lid
- Non sterile

Code	Plate	Well Volume, μL	Inner Pack	Pack	Price
3700	Clear	112	25	100	£695.00




Corning® 384 Well Non-Binding Surface (NBS™) Plates



- Flat bottomed, without lids
- Sterile

Code	Colour	Description	Well Volume, μL	Inner Pack	Pack	Price
3640	Clear	—	112	25	100	£1007.00
3544	Black/clear	Low volume	50	10	50	£619.00



Corning® 384 Well Polypropylene Storage Microplates



- Polypropylene microplates offer both small volume and large volume (blocks) well formats to meet assay and storage requirements
- Resistant to many common organic solvents (e.g. DMSO, ethanol, methanol)
- Certified DNase and RNase free
- Clear



Code	Well Volume, μL	Bottom	Sterility	Inner Pack	Pack	Price
3964	180	Round	Sterile	5	25	£276.00
3965	180	Round	Non sterile	5	100	£864.00
3342	240	V	Sterile	5	50	£371.00
3347	240	V	Non sterile	5	50	£292.00