

## Eppendorf ThermoStat™ C Block Heater/Cooler Continued

eppendorf

## SmartBlocks™



Exchangeable in just seconds, these individually sensor-controlled SmartBlocks™ are the best solution for maintaining accurate, uniform temperatures. A wide variety for all common lab vessels support a broad range of applications - for use with Eppendorf ThermoMixer® C (heating/cooling and mixing) and ThermoStat™ C (heating/cooling).

- Fast and simple block exchange due to Eppendorf QuickRelease™ technology
- Optimised block design guarantees maximum temperature transfer to the sample
- Individually sensor-controlled Eppendorf SmartBlocks™ allow for block-specific calibration to provide maximum temperature accuracy and homogeneity
- Insulated thermoblocks ensure safe and ergonomic use
- Automatic block recognition for safe handling

Code	Description	Price
<b>E5360000038</b>	For 24 x 1.5mL reaction vessels, incl. transfer rack 1.5/2mL	£517.00
<b>E5364000024</b>	For 24 x reaction vessels with Ø up to 12mm, height 25-76mm	£642.00
<b>E5366000021</b>	For 8 x 15mL conical tubes	£528.00
<b>E5362000035</b>	For 24 x 2.0mL reaction vessels, incl. transfer rack 1.5/2mL	£528.00
<b>E5309000007</b>	For 8 x 5.0mL tubes	£517.00
<b>E5365000028</b>	For 4 x 50mL conical tubes	£517.00
<b>E5307000000</b>	For 384 well PCR plates, incl. lid	£642.00
<b>E5306000006</b>	For 96 well PCR plates, incl. lid	£642.00
<b>E5318000001</b>	For 24 x 1.5–2mL cryo tubes with max. external Ø 12.5mm, all base designs	£517.00
<b>E5363000039</b>	For microplates and deepwell plates, incl. lid	£642.00
<b>E5361000031</b>	For 24 x 0.5mL reaction vessels, incl. transfer rack 0.5mL	£517.00
<b>E5310000002</b>	For Eppendorf deepwell plates 96/1000µL, incl. lid	£642.00
<b>E5316000004</b>	For Eppendorf deepwell plates 96/500µL, incl. lid	£642.00
<b>E3880000151</b>	Transfer rack for 24 x 1.5/2.0 mL microtubes	£41.00
<b>E3880000305</b>	Transfer rack for 24 x 0.5mL microtubes	£40.50
<b>E5363000233</b>	Lid for SmartBlock™ 0.50–2.0mL, plates, PCR 96, PCR 384, and DWP	£116.00



A

## SLS Select Aspirators



- Autoclavable polypropylene
- Wide shoulder handle for easy carrying and pouring
- Ideal for storage and transporting large batches of culture media, distilled water and other solutions

Code	Capacity, L	Dims, w x h, mm	Price
<b>ASP1206</b>	10	245 x 388	£70.00
<b>ASP1208</b>	20	270 x 540	£83.00