

## Azlon® Heavy Duty Aspirators



- Squat form for stable location on shelves, bench top or in transit
- Moulded handle for safe and secure carrying
- Clear moulded graduations allowing quick and easy volume assessment
- Leakproof heavy duty stopcock opens and closes in a quarter turn
- Autoclavable polypropylene for easy cleaning and sterilisation
- Caps have moulded drill guides for easy location of probes and feed lines
- Autoclavable at 121°C

Code	Alt Ref	Capacity, L	Thread Size	Dims, w x h, mm	Price
ASP0280	BNP05A	5	45	220 x 270	£38.00
ASP0282	BNP10A	10	80	280 x 360	£47.80
ASP0283	BNP20A	20	80	350 x 420	£64.00



A

## Thermo Scientific™ Nalgene™ Round Carboys with Spigot



- Carboys have a spigot boss molded in one piece with the carboy for convenience, strength and leak resistance
- Suitable for collecting and dispensing distilled water, acids and reagents
- Made only from high quality, laboratory grade plastic materials for dependably low leachables and extractables
- Molded in graduations in 1 gallon and 5L increments for ease in dispensing
- Integral shoulder handles make lifting comfortable for four largest sizes
- Leakproof
- LDPE carboys are not suitable for autoclaving; polypropylene carboys can be autoclaved before refilling to eliminate bacterial growth when storing sterile solutions



A

Code	Alt Ref	Material	Capacity, L	Dims, h x Ø, mm	Neck Ø, mm	Price
ASP1170	2318-0010	LDPE	4	335 x 154	26	£51.00
ASP1172	2318-0020	LDPE	10	389 x 250	64	£84.00
ASP1174	2318-0050	LDPE	20	528 x 286	64	£119.00
ASP1176	2318-0065	LDPE	25	594 x 287	64	£140.00
ASP1178	2318-0130	LDPE	50	678 x 379	64	£211.00
ASP2020	2319-0020	Polypropylene	10	389 x 250	64	£134.00
ASP2022	2319-0050	Polypropylene	20	528 x 286	64	£168.00
ASP2024	2319-0130	Polypropylene	50	678 x 379	64	£296.00

## SAMCO Rectangular Aspirator



- High density polythene
- Integrally moulded handle
- Supplied with screw in bung and tap
- Fitted with polythene leakproof cap

Code	Alt Ref	Capacity, L	Dims, w x d x h, mm	Price
ASP1140	P025/03	25	330 x 220 x 460	£52.00



## Spare Tap

Code	Alt Ref	Price
ASP1142	P025/02	£5.50

