


## Thermo Scientific™ Nunc™ Cell Culture Treated TripleFlasks™



- Experience 500cm<sup>2</sup> total growth area within the footprint of a standard 175cm<sup>2</sup> flask
- Total culture area is achieved by way of the three parallel layers and a straight neck design that works effectively within automation systems
- Nunclon™ Delta treatment offers maximum adhesion for a broad range of cell types
- Recommended working volume is 200mL
- Promotes monolayer adherent cell culture
- Also available barcoded with large code 128
- Nunclon™ Delta certificate is included in each case for monolayer formation, cloning efficiency, non cytotoxic, non pyrogenic, and sterility



Supplied in packs of 32

Code	Closure Type	Barcode	Price
132867	Solid	No	£500.00
132913	Filter	No	£506.00
132920	Filter	Yes	£609.00

## Thermo Scientific™ Nunc™ Cell Culture Multidishes



- Nunclon™ Delta surface treatment promotes maximum adhesion for a broad range of cell types
- Open format with the raised well rim enables convenient and secure manipulation of the cells
- Flat well bottom and stackable
- Raised well rims to lower the risk of cross contamination
- Sterilised by gamma irradiation and certified to be non pyrogenic
- Excellent optical quality
- Alphanumerical codes for well identification
- Round well styles feature vented lids for gas exchange, recessed areas around wells act as a water reservoir to reduce evaporation from wells, and lids with notched corners for easy orientation
- Rectangular well styles have a large culture area and fit standard glass cover slips: 24 x 60mm in 4 well plate; and 22 x 22mm in 8 well plates



Code	No. of Wells	Growth Area, cm <sup>2</sup>	Well Volume, mL	Lid	Bottom	Inner Pack	Pack	Price
176740	4	1.9	1	Yes	Flat	4	120	£314.00
167063	4	21.8	5	Yes	Flat	10	100	£992.00
140675	6	9.6	3	Yes	Flat	1	75	£267.00
140685	6	9.6	3	Yes	Flat	5	85	£290.00
167064	8	10.5	3	Yes	Flat	10	100	£1053.00
150628	12	3.5	2	Yes	Flat	1	75	£278.00
142475	24	1.9	1	Yes	Flat	1	75	£280.00
142485	24	1.9	1	Yes	Flat	5	85	£308.00
150687	48	1.1	0.5	Yes	Flat	1	75	£319.00
152640	48	1.1	0.5	Yes	Flat	5	85	£349.00
156545	96	—	0.4	No	Flat	10	180	£406.00
161093	96	—	0.4	Yes	Flat	10	160	£438.00
167008	96	0.32	0.4	Yes	Flat	1	50	£159.00
168055	96	—	0.4	No	Flat	1	50	£126.00
143761	96	—	0.3	No	U	1	50	£127.00
163320	96	—	0.3	Yes	U	1	50	£156.00
168136	96	—	0.3	Yes	U	10	160	£530.00