


Thermo Scientific™ Nunc™ *in vitro* Fertilisation Plasticware Continued

CE Marked Tubes for IVF

Code	Description	Inner Pack	Pack	Price
137860	11mL centrifuge tube	5	300	£237.00



Thermo Scientific™ Nunc™ UpCell™



- Nunc™ UpCell™ surface enables harvesting of cells with high visibility and intact surface proteins for culture passing, single cell analysis and cell transplantation research
- Temperature responsive cell culture surface
- Release single cells and cell sheets/monolayers by lowering the temperature to 32°C
- Harvesting of cells with intact surface proteins
- No trypsinisation - preserve cell surface proteins
- No cell scraping - get high cell viability
- Releases adherent cells by reduction of temperature of the cell culture
- For culture passing, single cell analyses and cell transplantation research
- Enables harvesting of cell sheets and creations of 3D tissue models held together by normal cell junctions and extracellular matrix
- Minimal hands-on time



Multiwell plates

- Sterile, clear and flat bottomed

Supplied in packs of 6 (individually wrapped)

Code	No. of Wells	Price
174901	6	£238.00
174900	12	£238.00
174899	24	£241.00
174898	48	£205.00

Dishes

- Supplied with lids with air vents, sterile

Code	Dims, Ø x h, mm	Description	Inner Pack	Pack	Price
174903	60 x 15	—	5	30	£967.00
174906	60 x 15	With grid	5	30	£965.00
174902	100 x 17	—	1	6	£387.00
174905	100 x 17	With grid	1	6	£388.00