


## Thermo Scientific™ Nunc™ Plates with Covalent or Affinity Binding Surfaces



- Increase the sensitivity of your ELISA techniques by using either covalent or affinity capture binding to immobilise your molecule to the microtitre plate
- Covalent binding surfaces chemically attach the molecule to the plate from the NH<sub>2</sub> end of the molecule, this means all molecules are oriented the same way and are impossible to wash away
- Affinity binding surfaces using streptavidin allow biotinylated molecules to bind strongly to the surface also in the same orientation
- Both methods do not require a blocking step for the ELISA
- No blocking required speeding up assay time
- Significantly reduces the amount of molecule required for coating
- Well volume: 400µL
- 5 year shelf life of a pre-coated plate



Code	Surface Treatment	Colour	Inner Pack	Pack	Price
<b>436006</b>	Covalent binding	Clear	5	30	£683.00
<b>436007</b>	Covalent binding	White	5	30	£985.00
<b>436008</b>	Covalent binding	Black	5	30	£980.00
<b>436014</b>	Affinity binding (streptavidin)	Clear	1	15	£508.00
<b>436015</b>	Affinity binding (streptavidin)	White	1	15	£661.00
<b>436016</b>	Affinity binding (streptavidin)	Black	1	15	£661.00

## Thermo Scientific™ Nunc™ Immuno Standard Modules



- Developed for ELISA techniques
- Alphanumeric orientation systems
- Fit standard equipment
- High optical quality
- With Nunc™ StarWell™ modules, incubation time can be reduced by more than 50% or the sensitivity similarly increased
- Frames and modules available separately
- Modules remain in frame when inverted
- For fluorescence/luminescence based immunoassays and binding assays



Code	Description	Colour	No. of Wells	Surface Treatment	Well Volume, µL	Working Volume, µL	Bottom	Inner Pack	Pack	Price
<b>475078</b>	Non breakable	Clear	8	MaxiSorp	300	250	Round	10	60	£562.00
<b>468667</b>	Non breakable	Clear	8	MaxiSorp	400	350	Flat	5	60	£558.00
<b>434797</b>	Non breakable	Clear	8	MaxiSorp	400	350	Flat	20	120	£1794.00
<b>463201</b>	Breakable	White	8	MaxiSorp	350	250	Flat, curved edge	10	60	£641.00
<b>446471</b>	Breakable	Black	8	MaxiSorp	350	250	Flat, curved edge	10	60	£641.00
<b>475086</b>	Non breakable	Clear	8	PolySorp	300	250	Round	10	60	£562.00
<b>469078</b>	Non breakable	Clear	8	PolySorp	400	350	Flat	10	60	£558.00
<b>463200</b>	Breakable	White	8	PolySorp	350	250	Flat, curved edge	10	60	£608.00
<b>446473</b>	Breakable	Black	8	PolySorp	350	250	Flat, curved edge	10	60	£511.00
<b>444865</b>	Non breakable	Clear	8	PolySorp	350	250	Flat, curved edge	10	60	£558.00
<b>473717</b>	Non breakable	Clear	12	PolySorp	350	250	Flat, curved edges	10	60	£472.00
<b>473709</b>	Non breakable	Clear	12	MaxiSorp	350	250	Flat, curved edge	10	60	£472.00
<b>467679</b>	Non breakable	Clear	16	PolySorp	400	350	Flat	10	60	£399.00
<b>467466</b>	Non breakable	Clear	16	MaxiSorp	400	350	Flat	5	60	£393.00
<b>446469</b>	Breakable	Clear	96	MaxiSorp	350	250	U	10	60	£579.00