

Hanna Potentiometric Titrators Continued



| Code | Model | Description | Price |
|---------|-----------|--|----------|
| TIT1030 | HI-931 | Potentiometric titrator with 1 analogue board | £4954.00 |
| TIT1016 | HI-932-C1 | Potentiometric titrator with 1 analogue board | £8255.00 |
| TIT1018 | HI-932-C2 | Potentiometric titrator with 2 analogue boards | £8751.00 |

Accessories

| Code | Alt Ref | Description | Pack | Price |
|---------|-----------|-------------------------|------|----------|
| TIT1020 | HI-930100 | Dosing pump | Each | £2374.00 |
| TIT1022 | HI-930101 | Peristaltic dosing pump | Each | £2698.00 |
| TIT1024 | HI-930190 | Blank burette holder | Each | £60.00 |
| TIT1026 | HI-930302 | Propellers | 3 | £22.20 |
| TIT1028 | HI-930303 | PVDF propellers | 3 | £51.00 |
| TIT1040 | HI-900205 | 5mL burette syringe | Each | £324.00 |
| TIT1042 | HI-900210 | 10mL burette syringe | Each | £324.00 |
| TIT1044 | HI-900225 | 25mL burette syringe | Each | £324.00 |
| TIT1046 | HI-900250 | 50mL burette syringe | Each | £283.00 |

2
YEARS

Hanna HI-903 Karl Fischer Volumetric Titrator



The HI-903 combines an ultra-high precision titrant delivery system with optically regulated magnetic stirring, sophisticated endpoint determination, dynamic dosing and background drift correction algorithms. The result is an extremely adaptable titrator capable of titrating with superior accuracy and precision even for samples with low moisture content. The HI-903 dispenses the titrant, detects the endpoint and performs all necessary calculations automatically.

- Powerful software and intuitive menus
- Large, colour LCD display making it simple to view results
- Equipped with a solvent handling system to reduce cell conditioning time
- Can be connected directly to your laboratory balance via serial interface
- Choose from included methods or develop a custom method for almost any application or sample type
- Using a USB flash drive or connecting the titrator to the HI 900PC application, methods (standard and user) can be upgraded, stored or deleted
- Clip-Lock™ exchangeable burette system enables users to exchange reagent burettes in a matter of seconds
- Supplied with dual platinum pin electrode, dosing pump, 5mL burette assembly with tubing, air pump assembly with tubing, beaker and bottle top assemblies and all fittings, desiccant cartridges x4 with indicating desiccant, stir bar, waste bottle, calibration key, USB cable, power cable, HI 900PC application, USB flash drive, quality certificate, ISO 8655 burette compliance report and instruction manual binder

| Code | Alt Ref | Measuring Range | Resolution, ppm | Dims, w x d x h, mm | Price |
|---------|-----------|-----------------|-----------------|---------------------|----------|
| TIT1032 | HI-903-02 | <100ppm to 100% | 1 (0.0001%) | 390 x 380 x 350 | £5608.00 |