



A

## Azlon® Instrument Tray



- Rectangular fabricated polypropylene tray
- Ideal for sterilising equipment
- Rigid and robust
- 22mm deep vertical sides
- Can be autoclaved at 121°C

Code	Alt Ref	Dims, l x w x d, mm	Price
<b>TRA1054</b>	TWR022	610 x 350 x 22	<b>£160.00</b>



## Kartell Photographic Trays



- Deep, PVC trays with ribbed bases
- Ideal for photographic developing
- White

Code	Alt Ref	Dims, l x w x d, mm	Price
<b>TRA3000</b>	280	150 x 200 x 45	<b>£14.63</b>
<b>TRA3002</b>	281	200 x 250 x 60	<b>£23.00</b>
<b>TRA3006</b>	283	300 x 350 x 80	<b>£40.80</b>
<b>TRA3008</b>	284	330 x 430 x 90	<b>£57.00</b>
<b>TRA3010</b>	285	420 x 520 x 90	<b>£81.00</b>
<b>TRA3014</b>	288	420 x 540 x 180	<b>£169.00</b>
<b>TRA3012</b>	286	530 x 675 x 100	<b>£140.00</b>



## Metal Tray Trolleys

- Constructed of tubular metal with trays and four castor tyre wheels
- Trays supplied equals the number of tiers stated
- Maximum recommended loading: 120kg
- 560mm (2 tier) or 260mm (3 tier) clearance between shelves
- Dims, h x l x w: 915 x 800 x 550mm

Code	Description	Tiers	Price
<b>TR01000</b>	Aluminium trays	2	<b>£351.00</b>
<b>TR01002</b>	Aluminium trays	3	<b>£535.00</b>
<b>TR01004</b>	Stainless steel trays	2	<b>£418.00</b>



## Metal Tray Trolleys, Blue

- All steel general purpose trolleys with tubular end frames
- Fully welded
- 100mm swivel nylon castors
- Carrying capacity of 200kg in a uniformly distributed load
- Shelves have 30mm retaining lip to all 4 sides
- Supplied in sky blue (RAL 5015)
- Dimensions, l x w x h: 1060 x 460 x 950mm

Code	Tiers	Price
<b>TR04000</b>	2	<b>£299.00</b>
<b>TR04002</b>	3	<b>£382.00</b>