


SLS Select Centrifuge Tubes



- Certified RNase and DNase free
- Sterile
- Manufactured from virgin ultra clear polypropylene
- Leak proof flat HDPE cap
- Printed graduations and a large writing area for sample identifications
- Racked tubes are packaged in recyclable card racks
- Lot traceability
- RCF: 15000xg (15mL) or 12000xg (50mL)

Supplied in packs of 500




Code	Capacity, mL	Bottom	Packaging	Price
SLS8102	15	Conical	Bagged	£178.00
SLS8103	15	Conical	Racked	£199.00
SLS8106	50	Conical	Bagged	£220.00
SLS8110	50	Conical	Racked	£206.00
SLS8114	50	Skirted	Bagged	£217.00

Eppendorf Centrifuge Tubes

eppendorf

- Ergonomic operation with easy open and close, and superior comfort for single-handed handling
- Stable, upright position to minimise contamination
- Superior materials and precision molding for maximal tightness and minimal sample loss
- Reliable, safe labelling due to large writing area and flat, light coloured cap
- Sterile
- Pyrogen, DNase, RNase, and DNA free (human & bacterial)
- Each individual lot is checked and certified for quality and purity
- For use up to a maximum RCF of 19500xg



Code	Description	Capacity, mL	Ø, mm	Inner Pack	Pack	Price
E0030122151	Standard	15	17	50	500	£143.00
E0030122178	Standard	50	30	25	500	£172.00
E0030122259	Forensic DNA grade	15	17	1	100	£71.00
E0030122267	Forensic DNA grade	50	30	1	48	£40.20

Eppendorf 25mL Conical Tubes

eppendorf

- Reduced height gives easy sample access and requires less storage space
- Available with SnapTec® caps for single-handed opening and closing for fast liquid extraction and addition
- Screw caps have flattened and grooved sides for user-friendly ergonomic handling, stable upright positioning, and reduced contamination risk
- Centrifugation stability up to 17000xg
- Tightness is ensured from -86 up to 100°C
- Premium USP class 6 raw material and no use of slip agents, plasticizers, and biocides during manufacturing ensure highest sample integrity and optimal sample and pellet visibility
- SnapTec® cap tube can be autoclaved up to 121°C (20 min)



Code	Description	Cap	Inner Pack	Pack	Price
E0030118405	Eppendorf Quality™	SnapTec®	40	200	£53.00
E0030122410	Eppendorf Quality™	Screw	50	200	£57.00
E0030118413	PCR clean	SnapTec®	40	200	£57.00
E0030122429	PCR clean	Screw	50	200	£61.00
E0030118421	Sterile	SnapTec®	25	150	£46.60
E0030122437	Sterile	Screw	25	200	£64.00