

DURAN® GLS 80® Laboratory Bottles Continued



Code	Alt Ref	Capacity, mL	Dims, w x h, mm	Pack	Price
BOT2117	218603656	250	95 x 110	10	£106.00
BOT2110	1112627	500	101 x 153	10	£141.00
BOT2112	1112713	1000	101 x 223	10	£181.00
BOT2114	1112715	2000	136 x 253	10	£358.00
BOT2111	218606953	3500	160 x 276	Each	£108.00
BOT2116	1113949	5000	182 x 315	Each	£149.00
BOT2118	1113950	10000	227 x 390	Each	£261.00
BOT2120	1113951	20000	288 x 485	Each	£566.00
BOT2124	1200265	30000	340 x 548	Each	£768.00
BOT2126	1200154	50000	400 x 590	Each	£1355.00

Accessories

Code	Alt Ref	Description	Pack	Price
BOT2119	1160166	Spare GLS 80® pouring ring	10	£17.20
BOT2122	1112716	Spare GLS 80® cap	10	£26.90

DURAN® GLS 80® Baffled Laboratory Bottles



- More effective mixing due to three moulded in, vertical baffles
- Manufactured from borosilicate glass 3.3
- Flat inner base
- Baffles form a continuous, smooth inner surface making clean up easy
- Easy access to the bottle contents due to the wide diameter of the bottle neck
- Special thread design for rapid opening and closing in less than a turn
- Stability due to the large base area
- Easy to read scale with graduation marks to estimate the content
- Fired on, highly durable white ceramic
- Special GLS 80® thread design for rapid opening and closing in less than a turn



Code	Alt Ref	Capacity, mL	Dims, w x h, mm	Pack	Price
BOT1132	212863658	250	95 x 110	Each	£20.60
BOT1134	212864457	500	101 x 153	Each	£27.60
BOT1136	212865453	1000	101 x 223	Each	£35.90

DURAN® protect Laboratory Bottles



- Safety coated version of the standard laboratory bottles, combining the best of both worlds: virtually inert USP/EP Type I borosilicate glass 3.3 glass on the inside and a protective safety plastic layer on the outside
- High performance materials and robust ISO 4796-1 design provides a long working life
- Easy to read graduated volume scale, and marking field in permanent white enamel
- Permanent retrace code on each bottle for manufacturing lot traceability
- Safety coating meets the requirements of ASTM E1522
- Fully autoclavable
- External polymer coating provides scratch, leak* and splinter* protection (*2L bottles and less only)
- Highly resistance, long service life coating
- UV protection up to approx. 380nm due to external coating
- Complete with screw cap (PP, integral lip seal) and pouring ring (PP) for drip free pouring and clean, safe working
- Size 10 has specially shaped glass edge for improved pouring out (so that an additional plastic pouring ring is not required)



A