

## Eppendorf 25mL Conical Tubes Continued

eppendorf

### Starter Packs

Code	Description	Price
<b>E0030118502</b>	200 x PCR clean SnapTec® cap tubes, 4 x single tube stands, 6 x adapters for rotors with bore for 50mL conical tubes	<b>£93.00</b>
<b>E0030122500</b>	200 x PCR clean screw cap tubes, 4 x single tube stands, 6 x adapters for rotors with bore for 50mL conical tubes	<b>£96.00</b>



## Eppendorf Centrifuge Bottles

eppendorf

- A reusable solution for large-volume, low-speed applications
- Ideal for cell harvesting applications
- Bottles and lids autoclavable at 121°C, UV-resistant and dishwasher-proof
- Low opacity of polypropylene gives high visibility of pellets
- Tight seal and easy to open and close
- Compatible with optional aerosol-tight caps for the safe handling of hazardous substances
- For use with the large benchtop centrifuge 5910 R: 1000mL bottle designed for swing-bucket rotor S-4xUniversal, 400mL bottle fits into swing-bucket rotor S-4x400

Supplied in packs of 2

Code	Capacity, mL	Dims, Ø x h, mm	Max RCF, xg	Price
<b>E5910707002</b>	400	81 x 117	5030	<b>£175.00</b>
<b>E5910770006</b>	1000	118 x 133	4256	<b>£175.00</b>



## Corning® Centrifuge Tubes

CORNING

- Sterile, certified non-pyrogenic
- Polypropylene
- RNase/DNase free
- Feature black printed graduations with a large white marking spot
- Choice of CentriStar™ cap or plug seal caps
- 15mL and 50mL tubes available in racks or bulk packed in resealable sleeves
- 50mL self standing centrifuge tubes are bulk packed in resealable sleeves

Code	Capacity, mL	Cap Type	Max RCF, xg	Packaging	Pack	Price
<b>430052</b>	15	Plug seal	12000	Racked	500 (10 racks of 50)	<b>£177.00</b>
<b>430766</b>	15	Plug seal	12000	Sleeved	500 (20 racks of 25)	<b>£178.00</b>
<b>430790</b>	15	CentriStar™	12000	Racked	500 (20 sleeves of 25)	<b>£177.00</b>
<b>430791</b>	15	CentriStar™	12000	Sleeved	500 (10 sleeves of 50)	<b>£175.00</b>
<b>430290</b>	50	Plug seal	15500	Racked	500 (10 racks of 50)	<b>£203.00</b>
<b>430291</b>	50	Plug seal	15500	Sleeved	500 (10 sleeves of 50)	<b>£182.00</b>
<b>430897</b>	50, self standing	Plug seal	3000	Sleeved	500 (20 racks of 25)	<b>£219.00</b>
<b>430828</b>	50	CentriStar™	15500	Racked	500 (20 sleeves of 25)	<b>£203.00</b>
<b>430829</b>	50	CentriStar™	15500	Sleeved	500 (20 sleeves of 25)	<b>£182.00</b>
<b>430921</b>	50, self standing	Flat	3000	Sleeved	500 (20 sleeves of 25)	<b>£219.00</b>