



Simport® Scientific Micrewtubes® with O-Ring Seal



- Polypropylene tubes available either plain or with printed graduations and white marking area for sample identification
- Washer seal ensures a positive leakproof seal
- Caps do not have the attachment loops for users who prefer to remove the caps completely from the tubes when filling or sampling
- Perfect for cryogenic work
- Non skirted tubes can be centrifuged up to 20000xg, skirted tubes can be centrifuged up to 17000xg
- Will withstand temperatures from -196 to 121°C

Code	Alt Ref	Capacity, mL	Bottom	Graduations	Sterile	Cap	Pack	Price
TUB2200	T334-2	0.5	Skirted	No	No	Loose	1000	£242.00
TUB2202	T334-2S	0.5	Skirted	No	Yes	Screwed on	500	£103.00
TUB2204	T334-2SPR	0.5	Skirted	Yes	Yes	Screwed on	500	£164.00
TUB2206	T334-4	1.5	Skirted	No	No	Loose	1000	£201.00
TUB2208	T334-4S	1.5	Skirted	No	Yes	Screwed on	500	£114.00
TUB2210	T334-4SPR	1.5	Skirted	Yes	Yes	Screwed on	500	£153.00
TUB2212	T334-5	1.5	Conical	No	No	Loose	1000	£201.00
TUB2214	T334-5S	1.5	Conical	No	Yes	Screwed on	500	£117.00
TUB2216	T334-5SPR	1.5	Conical	Yes	Yes	Screwed on	500	£164.00
TUB2218	T334-6	2.0	Skirted	No	No	Loose	1000	£179.00
TUB2220	T334-6S	2.0	Skirted	No	Yes	Screwed on	500	£114.00
TUB2222	T334-6SPR	2.0	Skirted	Yes	Yes	Screwed on	500	£156.00
TUB2224	T334-7	2.0	Conical	No	No	Loose	1000	£201.00
TUB2226	T334-7S	2.0	Conical	No	Yes	Screwed on	500	£165.00
TUB2228	T334-7SPR	2.0	Conical	Yes	Yes	Screwed on	500	£164.00



Eppendorf Safe-Lock Tubes



- The hinged lid on Eppendorf Safe-Lock microcentrifuge tubes eliminates any danger of leaks
- Ensures absolute safety when working with expensive, toxic samples, radioactive substances or DNA
- Volume markings and labelling marking surfaces
- The lid hook prevents the tube from opening accidentally (e.g. during heating)
- Tubes can be autoclaved when open (12°C, 20mins)
- Can be easily and single-handedly opened
- Excellent mechanical stability during centrifugation, up to 30000xg for 0.5, 1.5 and 2mL tubes
- Safe-Lock Tubes are available in Biopur®, PCR clean and Eppendorf Quality™ purity levels

Biopur®

Code	Capacity, mL	Colour	Pack	Price
E0030121570	0.5	Colourless	50	£27.50
E0030121589	1.5	Colourless	100	£28.10
E0030121597	2.0	Colourless	100	£37.40

PCR Clean

Code	Capacity, mL	Colour	Pack	Price
E0030123301	0.5	Colourless	500	£44.40
E0030123328	1.5	Colourless	1000	£44.40
E0030123344	2.0	Colourless	1000	£61.00