

Eppendorf Tubes® 5.0mL Continued

eppendorf

Code	Purity Type	Inner Pack	Pack	Price
E0030119487	Sterile	20	200	£45.00
E0030119401	Eppendorf Quality™	100	200	£36.40
E0030122305	Eppendorf Quality™	100	200	£46.00
E0030119460	PCR clean	100	200	£40.80
E0030119452	PCR clean	100	200	£46.50
E0030119479	Biopur®	1	50	£45.00
E0030122402	Forensic DNA grade	100	200	£61.00
E0030119606	Forensic DNA grade	50	200	£46.30

Starter Packs

Code	Description	Price
E0030119380	Starter pack: 400 PCR clean tubes, 2 tube racks, 8 adapters for 15mL conical tube rotors	£130.00

Accessories

Code	Description	Pack	Price
E0030119509	Tube clip	10	£21.50
E0030119827	16 well tube rack	2	£40.50
E0030140613	Storage box for 5 x 5 screw cap tubes (76mm height)	2	£17.75
E0030140532	Storage box for 5 x 5 snap cap tubes (64mm height)	4	£35.10


CLEARLine® Microtubes

ClearLine

- Flat cap adapted for easy writing
- Centrifugable up to 25000xg
- Made from transparent polypropylene
- Graduated
- Autoclavable at 121°C (for 20 minutes)
- DNase/RNase, human DNA and pyrogen free

Supplied in packs of 1000 (2 bags of 500)

Code	Alt Ref	Volume, mL	Price
TUB1030	CL022	1.5	£23.70
TUB1032	CL0269	2.0	£25.20



KEY

Throughout the catalogue the following symbols are used:



New product



Service and repairs available



Calibration service available



Made in Britain



Sustainable features



Autoclavable



No. years warranty



No. years warranty can be extended



Lifetime warranty