

Corning® Falcon® Round Bottom Test Tubes



- Provide reliable containment of laboratory fluid samples
- Polystyrene tubes are best suited for procedures requiring high optical clarity, polypropylene tubes are best suited for applications requiring greater thermal and chemical stability
- Dual-position snap caps offer both vented and fully closed positions, preventing sample loss with ease of handling
- Medical-style packaging materials
- Tube dimensions and volumes are approximate
- 5mL tubes are compatible with flow cytometers



Code	Capacity, mL	Material	Max Speed, x g	Cap	Sterility	Dims, Ø x h, mm	Inner Pack	Pack	Price
352003	5	Polystyrene	1400	Snap	Sterile	12 x 75	1	500	£126.00
352058	5	Polystyrene	1400	Snap	Sterile	12 x 75	25	500	£96.00
352054	5	Polystyrene	1400	Snap	Sterile	12 x 75	125	1000	£188.00
352052	5	Polystyrene	1400	No	Sterile	12 x 75	125	1000	£132.00
352063	5	Polypropylene	3000	Snap	Sterile	12 x 75	25	500	£101.00
352053	5	Polypropylene	3000	No	Sterile	12 x 75	125	1000	£154.00
352002	5	Polypropylene	3000	No	Non sterile	12 x 75	—	1000	£92.00
352027	8	Polystyrene	1400	Screw	Sterile	13 x 100	25	1000	£312.00
352001	14	Polystyrene	1400	Snap	Sterile	17 x 100	1	500	£138.00
352057	14	Polystyrene	1400	Snap	Sterile	17 x 100	25	500	£115.00
352051	14	Polystyrene	1400	Snap	Sterile	17 x 100	125	1000	£199.00
352017	14	Polystyrene	1400	No	Non sterile	17 x 100	125	1000	£143.00
352037	16	Polystyrene	1400	Screw	Sterile	16 x 125	1	500	£191.00
352025	16	Polystyrene	1400	Screw	Sterile	16 x 125	125	1000	£322.00
352006	14	Polypropylene	3000	Snap	Sterile	17 x 100	1	500	£143.00
352059	14	Polypropylene	3000	Snap	Sterile	17 x 100	25	500	£121.00
352018	14	Polypropylene	3000	No	Sterile	17 x 100	125	1000	£172.00
352045	19	Polystyrene	1400	Screw	Sterile	16 x 150	1	500	£199.00

Snap Caps

Polyethylene dual-position snap caps for 12 x 75mm tubes offering both vented and fully closed positions. Sterilised by gamma irradiation. RNase/DNase-free.


Code	Material	Pack	Price
352032	Polyethylene	2000	£115.00

SLS Select Disposable Test Tubes



- Round-bottomed disposable test tubes
- Caps are available and must be ordered separately
- Available in polystyrene or polypropylene
- Non sterile

Code	Nominal Capacity, mL	Material	Dims, Ø x h, mm	Pack	Price
TUB0096	3	Polystyrene	11 x 55	4000	£146.00
TUB0098	3	Polypropylene	11 x 55	4000	£147.00
TUB0105	3.5	Polystyrene	11.8 x 75	3000	£98.23
TUB0106	3.5	Polypropylene	11.8 x 75	2500	£86.86
TUB0108	12	Polystyrene	16 x 95	750	£69.00



Polyethylene Caps

Code	Description	Pack	Price
TUB0121	11mm	15000	£321.58
TUB0122	12mm	2000	£57.00
TUB0124	16mm	1000	£40.60

