



PYREX® Medium Wall Test Tubes

PYREX®

- Conforms to ISO 4142
- High quality test tubes designed and constructed for optimum performance
- Manufactured from PYREX® borosilicate glass which gives superior performance to neutral borosilicate and soda-lime glass
- Allows use at temperatures up to 500°C for short periods
- Excellent resistance to chemical attack

Rimless

Code	Alt Ref	Dims, Ø x h, mm	Wall Thickness, mm	Approx. Capacity, mL	Stopper Size	Pack	Price
TES2030	1622/02M	10 x 75	1.0	4	8	100	£47.80
TES2032	1622/04M	12 x 75	1.0	6	10	100	£50.00
TES2034	1622/06M	12 x 100	1.0	8	10	100	£54.50
TES2036	1622/09M	16 x 100	1.2	14	13	100	£71.00
TES2038	1622/10M	16 x 125	1.2	17	13	100	£71.00
TES2040	1622/12M	16 x 150	1.2	22	13	100	£71.00
TES2042	1622/14M	16 x 160	1.2	23	13	100	£85.00
TES2044	1622/15M	18 x 180	1.2	34	15	100	£89.00
TES2046	1622/16M	18 x 150	1.2	28	15	100	£86.00
TES2048	1622/20M	24 x 150	1.2	50	21	100	£127.00
TES2050	1622/22M	24 x 200	1.2	66	21	50	£91.00

Rimmed

Code	Alt Ref	Dims, Ø x h, mm	Wall Thickness, mm	Approx. Capacity, mL	Stopper Size	Pack	Price
TES2000	1620/02M	10 x 75	1.0	4	8	100	£47.80
TES2002	1620/04M	12 x 75	1.0	6	10	100	£50.20
TES2004	1620/06M	12 x 100	1.0	8	10	100	£54.70
TES2006	1620/09M	16 x 100	1.2	14	13	100	£71.00
TES2008	1620/10M	16 x 125	1.2	18	13	100	£71.00
TES2010	1620/12M	16 x 150	1.2	22	13	100	£71.00
TES2014	1620/15M	18 x 180	1.2	34	15	100	£89.00
TES2016	1620/16M	18 x 150	1.2	28	15	100	£86.00
TES2018	1620/20M	24 x 150	1.2	50	21	100	£127.00
TES2020	1620/22M	24 x 200	1.2	66	21	50	£91.00

SAMCO Test Tubes, Borosilicate Glass



- Suitable for all usual laboratory applications
- Withstand temperature change up to and including boiling of samples
- Highly chemically resistant



Heavy Wall, Rimmed

Code	Alt Ref	Dims, Ø x h, mm	Wall Thickness, mm	Pack	Price
TES1340	G032/01	10 x 75	1.5	100	£34.80
TES1342	G032/02	10 x 75	1.5	100	£37.00
TES1344	G032/03	12 x 100	1.5	100	£41.80
TES1346	G032/05	16 x 125	1.8	100	£65.00
TES1348	G032/06	16 x 150	1.8	100	£77.00
TES1350	G032/08	18 x 150	1.8	100	£92.00
TES1352	G032/10	24 x 150	1.8	50	£57.00