

Thermo Scientific™ Sterilin™ Silicone Tubing Continued

Thermo
SCIENTIFIC

Code	Alt Ref	Bore, mm	Wall Thickness, mm	Ext Ø, mm	Price
TUB7076	TSR0800240P	8	2.4	12.8	£262.00
TUB7078	TSR0800320P	8	3.2	14.4	£335.00
TUB7080	TSR0800400P	8	4	16	£497.00
TUB7084	TSR0850400P	8.5	4	16.5	£524.00
TUB7086	TSR0900250P	9	2.5	14	£320.00
TUB7088	TSR0900320P	9	3.2	15.4	£365.00
TUB7090	TSR0900400P	9	4	17	£453.00
TUB7094	TSR1000200P	10	2	14	£273.00
TUB7098	TSR1270320P	12.7	3.2	19.1	£470.00
TUB7102	TSR1900320P	19	3.2	25.4	£693.00



Thermo Scientific™ Nalgene™ Non-Phthalate PVC Tubing

Thermo
SCIENTIFIC

- Manufactured from the highest-grade resin materials, free from ortho-phthalate DEHP plasticisers
- Resists hardening and discoloration from a wide range of chemicals and does not impart odour or taste
- Ideal for a wide range of non-medical, general lab, research and pharmaceutical applications
- Equivalent to Tygon™ E-3603 tubing
- Marked with I.D. and regulatory compliance as well as at 10cm intervals for easy identification and measurement
- Durometer hardness is Shore A, 60
- Slips easily over glass and plastic fittings; grips securely
- Autoclavable
- Transparent

Supplied in 25m lengths

Code	Alt Ref	Bore, mm	Wall Thickness, mm	Ext Ø, mm	Price
TUB6268	8703-0406	4	1	6	£42.60
TUB6270	8703-0407	4	1.5	7	£54.00
TUB6272	8703-0508	5	1.5	8	£73.00
TUB6274	8703-0609	6	1.5	9	£97.00
TUB6276	8703-0610	6	2	10	£111.00
TUB6278	8703-0710	7	1.5	10	£116.00
TUB6280	8703-0811	8	1.5	11	£121.00
TUB6282	8703-0812	8	2	12	£127.00
TUB6284	8703-0913	9	2	13	£134.00
TUB6286	8703-1014	10	2	14	£155.00
TUB6288	8703-1216	12	2	16	£171.00
TUB6290	8703-1520	15	2.5	20	£76.00



Gas Burner Tubing

- Neoprene rubber tubing with reinforced ends
- Ideal for use with bunsen burners

Code	Length, mm	Price
TUB2100	450	£11.82
TUB2102	600	£16.37
TUB2104	900	£14.97