


SAMCO Screw Tubing Clip



- Closes rubber and flexible PVC tubing of up to 14mm external diameter
- Hinged side allows it to be clamped around the tubing, already connected to the apparatus
- Screw adjustment ensures fine flow control
- Moulded in acetal copolymer which is easy to clean, can be sterilised and is resistant to most acids and alkalis

Code	Alt Ref	Tube Ø, mm	Price
CLI2060	L145/01	14	£2.20



Unex® H Type Hose Clips

- Requires no tools - the hexagonal cap is intended for hand tightening
- Large enough to be manipulated by gloved hands
- Constructed from a stainless steel band, brass screw and bush
- Durable and rustproof
- Light, compact and powerful

Supplied in packs of 10

Code	Tube Ø, mm	Price
CLI1000	8 - 10	£39.30
CLI1002	9 - 12	£39.30
CLI1004	10 - 14	£39.30
CLI1006	13 - 17	£39.30
CLI1008	15 - 19	£39.30
CLI1010	18 - 22	£39.30
CLI1012	21 - 25	£39.30
CLI1014	23 - 27	£46.10
CLI1016	26 - 30	£46.10
CLI1018	29 - 33	£46.10
CLI1020	31 - 35	£46.10
CLI1024	37 - 41	£46.10



Kartell Tube Connectors, Straight



- Tapered stepped arms
- Autoclavable polypropylene
- Dimensions are nominal, and elasticity of hoses will influence choice of size

Supplied in packs of 10

Code	Alt Ref	Nominal OD, mm	Bore, mm	Price
CON7060	510	4 - 6	2	£7.00
CON7062	511	6 - 8	3.5	£25.00
CON7064	512	8 - 10	4.5	£34.00
CON7066	513	10 - 12	6.5	£35.00
CON7068	514	12 - 14	8.5	£15.00

