





Kartell Unequal Tube Connectors, Straight

Kartell
LABWARE

- Inlet and outlet have different diameters
- Autoclavable polypropylene
- Tapered arms
- Dimensions are nominal, and elasticity of hoses will influence choice of size

Supplied in packs of 10

Code	Alt Ref	Nominal OD, mm	Bore, mm	Price
CON7080	877	4 - 8 & 10 - 12	1.6 & 4.6	£15.00
CON7082	878	4 - 8 & 12 - 16	1.6 & 7.5	£16.00
CON7084	879	8 - 12 & 14 - 16	4.6 & 7.5	£20.00



Kartell Disconnectors, Straight

Kartell
LABWARE

- The two pieces fit tightly yet pull apart easily
- Useful for joining hoses of glass, rubber or plastic, replacing expensive adapter
- Supplied only with the same size at each end
- With tapered and rigid arms
- Polyethylene
- Dimensions are nominal, and elasticity of hoses will influence choice of size

Supplied in packs of 10

Code	Alt Ref	Nominal OD, mm	Bore, mm	Price
CON7100	434	4 - 6	1.6	£23.00
CON7102	435	6 - 8	2.7	£27.00
CON7104	436	8 - 10	3.6	£32.00
CON7106	535	10 - 12	5.5	£31.00
CON7110	537	14 - 16	8.8	£37.00



Portex® Tube Connectors, T Shape

- Rigid, straight arms
- Polypropylene

Supplied in packs of 10

Code	Ø, mm	Price
CON7202	9	£30.00
CON7204	11	£38.00

View our Safety & PPE catalogue

To browse our new Safety and PPE catalogue, visit the SLS website:

scientificlabs.co.uk/literature (UK & Northern Ireland)

scientificlabs.ie/literature (Republic of Ireland)

scientificlabs.com/literature (International)

