

## GETINGE Lancer ULTIMA Series Labware Washing Systems Continued

**ULTIMA Mid Range Freestanding Washers**

The GETINGE Lancer ULTIMA series mid range freestanding washers include 3 high-performance models that offer a range of labware loading options using baskets, injector racks or combination racks for the best high capacity processing of multiple labware types in the industry.

	910LX	1300LX	1400LX
Dims, w x d x h, mm	609 x 752 x 1483	609 x 752 x 1643	789 x 808 x 1731
Effective Chamber Volume, L	160	180	304
Wash Temperature, °C	Up to 95		
Drying	HEPA filtered direct injection		
View In Progress (VIP)	Yes with chamber light		
Touchscreen Controller	3.5" colour stores up to 40 programs		
Low Pressure Valve & Pump Kit	Yes - for trouble free DI water connection		
Steam Condenser	Included	Option	Option
Water Softener	Included	Option	Option
Chemical Storage	Yes	Yes	Yes
Chemical Dosing Pump	2 pumps to dispense detergent & acid rinse		
Chemical Level Sensors	2 sensors for notification when chemicals are low		
Number of Wash Levels	2	3	3
Number of Runner Positions	2	4	4
Noise Level, dB(A)	<63		

Code	Alt Ref	Model	Price
<b>WAS2030</b>	1260317	910LX	<b>£16472.00</b>
<b>WAS2006</b>	01260339	1300LX	<b>£19972.00</b>
<b>WAS2008</b>	1260340	1400LX	<b>£28885.00</b>

**ULTIMA Premium Range Freestanding Washers**

The GETINGE Lancer ULTIMA series premium range includes all the features of the mid-range however the range includes for a larger touchscreen at 7", a more powerful PLC which allows the addition of options for greater process control, independent validation monitoring, conductivity reading and on-board printers.

	1400LXP	1600LXP	1800LXA
Dimensions, w x d x h, mm	789 x 808 x 1731	866 x 858 x 1814	1398 x 841 x 1999
Effective Chamber Volume, L	304	408	458
Wash Temperature, °C	Up to 95		
Drying	HEPA filtered direct injection - programmable to 1°C increments to 110°C		
View In Progress (VIP)	Yes with chamber light	Glass door	
Touchscreen Control Panel	7" colour stores up to 40 programs		
Low Pressure Valve & Pump Kit	Yes - for trouble free DI water connection		
Steam Condenser	Option	Option	Option
Water Softener	Option	Option	Option
Chemical Storage	Yes	Yes	Yes
Chemical Dosing Pumps	2 pumps to dispense detergent & acid rinse		
Chemical Level Sensors	2 sensors for notification when chemicals are low		
Number of Wash Levels	3	4	6
Number of Runner Positions	4	5	6
Noise Level, dB(A)	<63	<69	<65

Code	Alt Ref	Model	Price
<b>WAS2022</b>	01260341	1400LXP	<b>£30490.00</b>
<b>1600LXP</b>	01260342	1600LXP	<b>£32530.00</b>
<b>WAS2024</b>	01260018	1800LXA	<b>£40295.00</b>