

## Miele PG 8536 CD Washer Disinfectors

**Miele**

A glassware washer and disinfector with integrated pumps and process chemical drawer. The unit features the freely programmable Profitronic Control System and PerfectTouch control which allows easy programming by the end user, and up to 64 programme slots.

- 18 standard cleaning programmes
- DryPlus hot air drying with H13 HEPA filtration
- Spray arm sensing
- 2 integrated dispenser pumps for detergent/neutralising agent
- Drawer for 4 x 5L canisters
- Ultrasound flow control
- Optional: conductivity meter
- 3 phase 400V
- Network interface for process documentation
- Dims, w x d x h, mm: 900 x 700 x 1175



Note: services (water, drainage and 3-phase power) need to be prepared prior to delivery, installation and commissioning.

Code	Alt Ref	Description	Price
<b>WAS1076</b>	7831480	PG 8536 CD (for pressurised water supply)	£20653.00
<b>WAS1078</b>	7831490	PG 8536 CD (for unpressurised water supply)	£21188.00

## Miele Baskets and Insert Accessories

**Miele**

Code	Alt Ref	Description	Price
<b>WAS7200</b>	A100	Upper basket for the connection of up to two injector modules	£520.00
<b>WAS7202</b>	A101	Upper basket for various inserts, height adjustable with open front vertical clearance 160mm	£331.00
<b>WAS7204</b>	A102	Upper basket for various inserts, height adjustable with open front vertical clearance 205mm	£383.00
<b>WAS7206</b>	A103	Upper basket spray arm for mesh trays and other shallow loads	£383.00
<b>WAS7208</b>	A150	Lower basket for the connection of up to two injector modules, automatically self-sealing docking valves	£452.00
<b>WAS7210</b>	A151	Open front lower basket for various DIN mesh trays and various inserts	£210.00
<b>WAS7212</b>	A200	Load carrier mobile unit for pipettes	£1817.00
<b>WAS7214</b>	A202	Two level load carrier for mesh trays	£744.00
<b>WAS7232</b>	A303	Module for pipettes for the optimum loading of up to 98 pipettes	£1019.00
<b>WAS7234</b>	A304	Module for phials for the optimum loading of up to 98 vials	£1458.00
<b>WAS7236</b>	A306	Module for measuring cylinders	£1790.00
<b>A3003</b>	A300/3	Module for optimum loading of laboratory glassware – can be equipped with 8 nozzles	£257.00
<b>A3015</b>	A301/5	Module for optimum loading of laboratory glassware – can be equipped with 18 nozzles	£399.00
<b>A860</b>	A860	Holder rack for injector modules A 301/4 and A 301/5	£142.00
<b>A861</b>	A861	Holder rack for injector modules A 301/4 and A 301/5	£194.00

## A world of products in one place...

- Over 500,000 products
- Easy ordering capability
- Detailed product information
- Promotions and offers

