



SLS Lab Pro PurA-Q Choice Type 2 or 3 Water Purification Systems



Available in two water quality versions providing general laboratory water, suitable for glass washers and autoclave feed, stability chambers, or Type 2 for buffer, reagent & media preparation, sample dilution etc. With many unique features including:

- Low cost, one box solution
- Wall or bench mounted
- Suits all general applications
- 50% recovery rate keeps water costs to a minimum
- Energy consumption <60W
- Unique disposable RO maintains peak efficiency
- Consumables changed in seconds
- Built in 35L reservoir
- Easy to install, low maintenance
- Minimal water wastage
- Mercury free UV

| Specifications | |
|--|-----------------|
| Type 2 Water Quality, MΩ-cm | 1 to 15 |
| Type 3 Water Quality, µS/cm | <40 |
| TOC (Type 2 only), ppb | <5 |
| Organics Rejection Rate (Type 3 only), % | <99 |
| Feed Water Pressure, bar | 0.1 to 6 |
| Dims, w x d x h, mm | 485 x 330 x 500 |
| Weight, kg | 20 |

| Code | Model | Water Type | Make Up Rate, L/hr | Daily Usage (max), L/day | Tank UV Kit | Price |
|---------|------------------------|-------------|--------------------|--------------------------|-------------|----------|
| WAT2228 | PurA-Q Choice 10 T2 | Type 2 (DI) | 10 | 100 | No | £4340.00 |
| WAT2230 | PurA-Q Choice 10 T2 UV | Type 2 (DI) | 10 | 100 | Yes | £4333.00 |
| WAT2220 | PurA-Q Choice 10 T3 | Type 3 (RO) | 10 | 100 | No | £3757.00 |
| WAT2222 | PurA-Q Choice 10 T3 UV | Type 3 (RO) | 10 | 100 | Yes | £4004.00 |
| WAT2408 | PurA-Q Choice 20 T2 | Type 2 (DI) | 20 | 200 | No | £4339.00 |
| WAT2418 | PurA-Q Choice 20 T2 UV | Type 2 (DI) | 20 | 200 | Yes | £4435.00 |
| WAT2406 | PurA-Q Choice 20 T3 | Type 3 (RO) | 20 | 200 | No | £4250.00 |
| WAT2419 | PurA-Q Choice 20 T3 UV | Type 3 (RO) | 20 | 200 | Yes | £4333.00 |

Accessories and Consumables

| Code | Description | Price |
|---------|---|----------|
| WAT2456 | Wall mounting kit | £174.00 |
| WAT2454 | UV upgrade kit - tank mounted | £709.00 |
| WAT2234 | Distribution pump kit | £2214.00 |
| WAT2172 | Valve kit for washing machine feed (single POU) | £163.00 |
| WAT2436 | Pre-treatment module (RO) 20L/hr | £782.00 |
| WAT2438 | DI module | £243.00 |
| WAT2440 | Air vent filter | £225.00 |
| WAT2442 | UV lamp | £362.00 |
| WAT1136 | Annual consumable pack for WAT2228 - pre-treatment module (RO) 10L/hr, DI cartridge, air vent filter | £1091.00 |
| WAT1138 | Annual consumable pack for WAT2230 - pre-treatment module (RO) 10L/hr, 2 x DI cartridge, air vent filter, UV lamp | £1379.00 |
| WAT1140 | Annual consumable pack for WAT2220 - pre-treatment module (RO) 10L/hr, air vent filter | £704.00 |
| WAT1142 | Annual consumable pack for WAT2222 - pre-treatment module (RO) 10L/hr, air vent filter, UV lamp | £995.00 |
| WAT1144 | Annual consumable pack for WAT2408 - pre-treatment module (RO) 20L/hr, 2 x DI cartridge, air vent filter | £1133.00 |
| WAT1146 | Annual consumable pack for WAT2418 - pre-treatment module (RO) 20L/hr, 2 x DI cartridge, air vent filter, UV lamp | £1422.00 |
| WAT1148 | Annual consumable pack for WAT2406 - pre-treatment module (RO) 20L/hr, air vent filter | £746.00 |
| WAT1150 | Annual consumable pack for WAT2419 - pre-treatment module (RO) 20L/hr, air vent filter, UV lamp | £1034.00 |