



## Sartorius arium® mini Water Purification Systems

SARTORIUS

A compact ultrapure water series that is specially designed for Type 1 water requirements of about 10L per day. With the capacity to deliver ultrapure water at a rate of 1L/min and a width of only 28cm, this space-saving unit easily integrates into nearly any lab environment. Regardless of the type of feed water you have, use arium® mini or arium® mini plus with our unique bagtank technology or connect arium® mini essential directly to your deionised water supply line as needed.

- **arium® mini:** the standard version, which is independent of a permanently installed water tap, water is supplied to this system via the integrated bag, which is filled by a built in pump
- **arium® mini plus:** the version with expanded features, it can be directly connected to a feed water tap, in this system, feed water is treated first by a pre-treatment cartridge
- **arium® mini essential:** a space saving system which connects to a pre-treated water supply line (RO/DI/EDI)

Code	Alt Ref	Description	Price
<b>WAT4284</b>	H20-MM-T	arium® mini	<b>£2973.00</b>
<b>WAT4286</b>	H20-MM-UV-T	arium® mini with UV lamp	<b>£3524.00</b>
<b>WAT4280</b>	H20-MA-T	arium® mini plus	<b>£4011.00</b>
<b>WAT4282</b>	H20-MA-UV-T	arium® mini plus with UV lamp	<b>£4054.00</b>
<b>WAT1084</b>	H20-MU-T	arium® mini essential	<b>£3361.00</b>
<b>WAT1086</b>	H20-MU-UV-T	arium® mini essential with UV lamp	<b>£4334.00</b>

### Accessories

Code	Alt Ref	Description	Price
<b>WAT4288</b>	H20-CPR	Pre-treatment cartridge	<b>£307.00</b>
<b>WAT4290</b>	H20-S-Pack	Scientific kit ultrapure cartridge	<b>£249.00</b>
<b>WAT4292</b>	H20-CEL1	UV lamp (185/254nm)	<b>£221.00</b>
<b>WAT4294</b>	H20-CBS-5-S	5L bag	<b>£40.00</b>
<b>WAT4232</b>	5441307H4--CE	Sterile final filter 0.2µm [Sartopore® 2 150 capsule, 0.2µm]	<b>£88.00</b>
<b>WAT4300</b>	H20-CUF	Cell plus ultra filter	<b>£261.00</b>

## Sartorius arium® Water Purification Systems

SARTORIUS

	comfort	advance	pro
iJust Software	✓	✓	N
Bagtank System	✓	✓	N
Display with Touch-Activated Functions	✓	✓	✓
High Quality Features	✓	✓	✓
Maximum System Flexibility	✓	✓	✓
Highest Water Quality	✓	N	✓
SOP Monitoring	✓	✓	✓
Favorites Function	✓	N	✓
EDI Technology	comfort II	advance EDI	N
ASTM Type 1 Water	comfort I, comfort II	N	✓
Type 2 Water	comfort II	advance EDI	N
Type 3 Water	comfort I	advance RO	N