

Sartorius arium® Water Purification Systems Continued

SARTORIUS



arium® pro

The arium® pro ultrapure water systems provide the ultimate flexibility and an excellent cost-benefit ratio because their device configurations are specially tailored to your applications. All systems meet and exceed the ASTM Type 1 water quality standards and ensure the best reproducible results. Up to 2L/min of ultrapure water with a conductivity of 0.055µS/cm (~18.2 MQ+cm) can be dispensed in consistently high quality.

- Options of UV and Ultrafiltration
- Cartridge packs available for analytical, biological, standard and feedwater applications
- Option of a TOC monitor

Code	Alt Ref	Description	Mount	Price
WAT4242	H2ObasicT	Basic with reduced operating features type 1 ion exchange water system	Bench mounted	£4602.00
WAT4243	H2ObasicB	Basic with reduced operating features type 1 ion exchange water system	Wall mounted	£4602.00
WAT4244	H2OproDIT	DI type 1 ion exchange	Bench mounted	£6144.00
WAT4245	H2OproDIB	DI type 1 ion exchange	Wall mounted	£6144.00
WAT4247	H2OproUFT	UF type 1 ion exchange and UF	Bench mounted	£7335.00
WAT4248	H2OproUFB	UF type 1 ion exchange and UF	Wall mounted	£7595.00
WAT4250	H2OproUVT	UV type 1 ion exchange and UV	Bench mounted	£7335.00
WAT4251	H2OproUVB	UV type 1 ion exchange and UV	Wall mounted	£7335.00
WAT4253	H2OproUVTTOC	UV type 1 ion exchange and UV, with TOC monitor	Bench mounted	£9000.00
WAT4254	H2OproUVBTOC	UV type 1 ion exchange and UV, with TOC monitor	Wall mounted	£9000.00
WAT4256	H2OproVFT	VF type 1 ion exchange and UV and UF	Bench mounted	£8506.00
WAT4257	H2OproVFB	VF type 1 ion exchange and UV and UF	Wall mounted	£8506.00
WAT4259	H2OproVFTTOC	VF type 1 ion exchange and UV and UF, with TOC monitor	Bench mounted	£10283.00
WAT4260	H2OproVFBTOC	VF type 1 ion exchange and UV and UF, with TOC monitor	Wall mounted	£10283.00



arium® bagtank

The innovative arium® bagtank system replaces conventional storage tanks and, unlike conventional systems, utilises a replaceable bag for storage. Easy to exchange, the arium® bagtank protects pure water against secondary airborne contamination eliminating the need for time consuming and labour intensive cleaning procedures. Available with or without an integrated pump, the addition of rollers on the arium® bagtank let you move it easily to wherever you need to use it.

Code	Alt Ref	Capacity, L	Description	Price
WAT4268	H2OAOV20	20	Without pump	£841.00
WAT4269	H2OAOV50	50	With pump	£1496.00
WAT4271	H2OAOV50W	50	Without pump	£1171.00
WAT4270	H2OAOV100	100	With pump	£1854.00
WAT4272	H2OAOV100W	100	Without pump	£1528.00

Accessories

Code	Alt Ref	Description	Pack	Price
WAT4220	H2OCPACK	Comfort ion exchange cartridge	Each	£348.00
WAT4222	611CEL1	Dual wavelength UV lamp	Each	£336.00
WAT4223	613CPM4V	RO module	Each	£822.00
WAT4225	H2OCCS	RO module cleaning set	Each	£64.00
WAT4230	H2OAFS1	Foot switch	Each	£427.00
WAT4231	611APR1	Printer	Each	£845.00
WAT4232	5441307H4--CE	Sterile final filter 0.2µm (Sartopore® 2 150 capsule, 0.2µm)	Each	£88.00
WAT4296	H2O-CPFCO-1	Comfort prefilter	Each	£309.00
WAT4298	H2O-CSO-1	Softener cartridge	Each	£84.00