

## Xylem WTW photoLab® 7000 Spectrophotometers Continued



### PC Software

Code	Alt Ref	Description	Price
<b>WAT5111</b>	902763	PhotoLab® Data spectral PC software	<b>£304.00</b>

## Xylem WTW Water Test Kits



Depending on the application, there are a variety of test kits available for routine investigations. Photometers and test kits together form a system in which each is coordinated with the other depending on optics and the wavelength used.

- Quick, reliable, cost effective
- The right test for every application
- Guaranteed results through AQA/IQC
- Suitable for laboratory, infrequent work, or ease with very large sample sizes
- Round cuvette tests: reaction cuvette with barcode or method selection, 16 mm: sample addition, insertion, measurement and reading with minimum work effort, AQA support for stored results
- Reagent tests: large measurement range for 10, 20 and 50mm rectangular cuvettes, recording of the smallest concentrations in rectangular cuvettes up to 50mm, AQA support for stored results

Code	Alt Ref	Parameter	Measuring Range	Type	Pack	Price
<b>WAT4302</b>	250495	Ammonium (NH <sub>4</sub> )	0.010 to 2000mg/L NH <sub>4</sub> -N, 0.01 - 2.58mg/L NH <sub>4</sub> <sup>+</sup>	Round cuvette test	25	<b>£94.00</b>
<b>WAT4306</b>	250424	Ammonium (NH <sub>4</sub> )	4.0 - 80.0mg/L NH <sub>4</sub> -N, 5.2 - 103.0mg/L NH <sub>4</sub> <sup>+</sup>	Round cuvette test	25	<b>£94.00</b>
<b>WAT4304</b>	250329	Ammonia (NH <sub>3</sub> )/ Ammonium (NH <sub>4</sub> )	0.09 - 3.00mg/L NH <sub>3</sub> (pH 8.5/25°C), 0.5 - 16.0mg/L NH <sub>4</sub> -N, 0.6 - 20.6mg/L NH <sub>4</sub> <sup>+</sup>	Round cuvette test	25	<b>£95.00</b>
<b>WAT4308</b>	252081	Ammonia (NH <sub>3</sub> )/ Ammonium (NH <sub>4</sub> )	0.002 - 0.730mg/L NH <sub>3</sub> (pH 8.5/25°C), 0.010 - 3.00mg/L NH <sub>4</sub> -N, 0.013 - 3.86mg/L NH <sub>4</sub> <sup>+</sup>	Reagent test	250	<b>£165.00</b>
<b>WAT4310</b>	250420	Chlorine (Cl <sub>2</sub> )	0.03 - 6.00 Cl <sub>2</sub> free & total	Round cuvette test	200*	<b>£76.00</b>
<b>WAT4312</b>	252070	Chemical Oxygen Demand (COD)	10 - 150mg/L COD (148°C, 2h)	Round cuvette test	25	<b>£95.00</b>
<b>WAT4314</b>	252071	Chemical Oxygen Demand (COD)	25 - 1500mg/L COD (148°C, 2h)	Round cuvette test	25	<b>£95.00</b>
<b>WAT4316</b>	250309	Chemical Oxygen Demand (COD)	500 - 10000mg/L COD (148°C, 2h)	Round cuvette test	25	<b>£95.00</b>
<b>WAT4318</b>	250410	Nitrate (NO <sub>3</sub> )	0.5 - 18.0mg/L NO <sub>3</sub> -N, 2.2 - 79.7mg/L NO <sub>3</sub>	Round cuvette test	25	<b>£117.00</b>
<b>WAT4320</b>	252019	Nitrate (NO <sub>3</sub> )	23 - 225mg/L NO <sub>3</sub> -N, 102 - 996mg/L NO <sub>3</sub>	Round cuvette test	25	<b>£169.00</b>
<b>WAT4322</b>	251993	Nitrate (NO <sub>3</sub> )	0.2 - 30.0mg/L NO <sub>3</sub> -N, 1 - 133.0mg/L NO <sub>3</sub>	Round cuvette test	50	<b>£78.00</b>
<b>WAT4324</b>	252074	Nitrite (NO <sub>2</sub> )	0.010 - 0.700mg/L NO <sub>2</sub> -N, 0.03 - 2.30mg/L NO <sub>2</sub>	Round cuvette test	25	<b>£89.00</b>
<b>WAT4326</b>	252096	Organic acids (volatile)	50 - 3000mg/L	Round cuvette test	25	<b>£192.00</b>
<b>WAT4328</b>	252075	Phosphate (PO <sub>4</sub> )	0.05 - 5.00mg/L PO <sub>4</sub> -P, 0.05 - 5.0mg/L P <sub>ges</sub> , 0.2 - 15.3mg/L PO <sub>4</sub>	Round cuvette test	25	<b>£106.00</b>
<b>WAT4330</b>	252076	Phosphate (PO <sub>4</sub> )	0.5 - 25.0mg/L PO <sub>4</sub> -P, 0.5 - 25.0mg/L P <sub>ges</sub> , 1.5 - 76.7mg/L PO <sub>4</sub>	Round cuvette test	25	<b>£106.00</b>
<b>WAT4332</b>	252045	Phosphate (PO <sub>4</sub> )	1.0 - 100.0mg/L PO <sub>4</sub> -P, 3.0 - 307.0mg/L PO <sub>4</sub>	Reagent test	100	<b>£151.00</b>
<b>WAT4334</b>	252036	Total Organically-bound Carbon (TOC)	5.0 - 80.0mg/L TOC	Round cuvette test	25	<b>£519.00</b>
<b>WAT4338</b>	252038	Screw caps	—	—	25	<b>£147.00</b>

\* 200 = 100 Cl<sub>2</sub> free + 100 Cl<sub>2</sub> total