



### Water Test Strips

Test strips to detect residual levels of peracetic acid using a visual qualitative method. Available to test high or low levels, and gives results in just seconds.

- Desiccant vial with flip top lid
- Water resistant label
- 2 year shelf life

Supplied as vials of 50 strips

Code	Description	Range, ppm	Price
<b>WAT4800</b>	Low level test strips	0 - 5 - 10 - 20 - 30 - 50	<b>£45.80</b>
<b>WAT4802</b>	High level test strips	0 - 100 - 250 - 500	<b>£48.00</b>

## Looking for Turbidity Standards?

Visit [scientificlabs.co.uk/turbidity](https://scientificlabs.co.uk/turbidity) (UK & Northern Ireland), [scientificlabs.ie/turbidity](https://scientificlabs.ie/turbidity) (Republic of Ireland) or [scientificlabs.com/turbidity](https://scientificlabs.com/turbidity) (International).



Our service options:

- Project Design & Management
- Installation
- Qualification
- Repair
- Maintenance
- Calibration
- Validation
- Spare Parts
- Training

## Clarity in Servicing



**Our Knowledge + Your Equipment = Service Tailored to You**

**Email:** [servicecentre@scientific-labs.com](mailto:servicecentre@scientific-labs.com)

**Telephone:** 0115 977 5075