



A

SLS Select Wide Mouth Bottles, PP



- Autoclavable
- Polypropylene
- With cap

Code	Capacity, mL	Dims, h x Ø, mm	Pack	Price
BOT2220	30	62 x 35	12	£8.60
BOT2222	60	85 x 40	12	£11.20
BOT2224	125	100 x 51	12	£13.40
BOT2226	250	132 x 61	12	£20.80
BOT2230	1000	217 x 93	6	£27.60



SLS Select Wide Mouth Bottles, LDPE



- Screw thread with cap attached
- Low density polyethylene

Supplied in packs of 6

Code	Capacity, mL	Dims, h x Ø, mm	Price
SLS3548	500	169 x 73	£14.90
SLS3550	1000	217 x 93	£23.70



SLS Select Amber Wide Mouth Bottles, HDPE



- Amber
- High density polyethylene

Code	Capacity, mL	Dims, h x Ø, mm	Pack	Price
BOT2260	30	62 x 35	12	£10.20
BOT2262	60	85 x 40	12	£22.70
BOT2264	125	100 x 51	12	£16.00
BOT2266	250	132 x 61	12	£24.30
BOT2268	500	169 x 73	6	£16.10
BOT2270	1000	217 x 93	6	£35.50



Thermo Scientific™ Nalgene™ Wide Mouth Lab Quality Bottles, HDPE



- Excellent chemical resistance to most acids, bases and alcohols
- Ideal for storing, shipping and packaging liquids
- High quality, laboratory grade plastic materials for dependably low leachables and extractables
- Uniform walls provide durability and resistance to splitting or punctures
- Suitable for freezer use to -100°C
- Linerless caps provide the ultimate in leakproof protection without the use of a liner that can wrinkle, cause leaks, or contaminate your reagents
- Leakproof

Code	Alt Ref	Capacity, mL	Dims, Ø x h, mm	Closure Size, mm	Pack	Price
BOT0080	2104-0001	30	34 x 63	28	12	£12.80
BOT0082	2104-0002	60	39 x 86	28	12	£15.20
BOT0084	2104-0004	125	50 x 99	38	12	£21.00
BOT0086	2104-0008	250	61 x 131	43	12	£27.10
BOT0088	2104-0016	500	73 x 168	53	12	£44.30
BOT0090	2104-0032	1000	91 x 199	63	6	£29.10
BOT0093	2104-0048	1500	91 x 284	63	6	£40.60