

Thermo Scientific™ Nalgene™ Wide Mouth Economy Bottles, HDPE

Thermo
SCIENTIFIC

- A lighter-weight and more flexible version of Nalgene™ lab quality bottles
- Recommended for sample collection and storage of daily working solutions
- Wide mouth makes these bottles easy to fill with solids or liquids
- Linerless caps ensure leakproof performance without a liner that can wrinkle, cause leaks and contaminate your solutions
- Leakproof



Code	Alt Ref	Capacity, mL	Dims, Ø x h, mm	Closure Size, mm	Pack	Price
BOT5050	2189-0001	30	34 x 63	28	72	£57.00
BOT5052	2189-0002	60	39 x 86	28	72	£61.00
BOT5054	2189-0004	125	50 x 99	38	72	£65.00
BOT5056	2189-0008	250	61 x 131	43	72	£92.00
BOT5058	2189-0016	500	73 x 168	53	48	£102.00
BOT5059	2189-0032	1000	91 x 199	63	24	£78.00

Thermo Scientific™ Nalgene™ Wide Mouth Bottles, PPCO

Thermo
SCIENTIFIC

- Store specimens where autoclaving is required
- High quality, laboratory grade plastic materials for reliably low leachables and extractables.
- Excellent chemical resistance to most acids, bases and alcohol
- Translucent for easy viewing of liquid levels
- Linerless caps for leakproof performance without a liner that can wrinkle, cause leaks and contaminate your solutions
- Autoclavable
- Leakproof
- Manufactured from non toxic materials



Code	Alt Ref	Capacity, mL	Dims, Ø x h, mm	Closure Size, mm	Pack	Price
BOT0100	2105-0001	30	34 x 63	28	12	£13.60
BOT0102	2105-0002	60	39 x 86	28	12	£15.80
BOT0104	2105-0004	125	50 x 99	38	12	£24.60
BOT0106	2105-0008	250	61 x 131	43	12	£32.20
BOT0108	2105-0016	500	73 x 168	53	12	£54.00
BOT0110	2105-0032	1000	91 x 199	63	6	£31.80

Thermo Scientific™ Nalgene™ Wide Mouth Bottles, LDPE

Thermo
SCIENTIFIC

- Translucent low density polyethylene
- Ideal for storing chemicals and standards, as well as collecting and storing distilled water
- High quality, laboratory grade plastic materials for dependably low leachables and extractables
- Extremely low trace metal content in LDPE makes it ideal for trace metal analysis
- Excellent chemical resistance to most acids, bases and alcohols
- Linerless caps ensure leakproof performance without a liner that can wrinkle, cause leaks and contaminate solutions
- Leakproof
- Can be used in freezer to -100°C for long term storage
- Also available sterile and in bulk packs - please enquire for details



Code	Alt Ref	Capacity, mL	Dims, Ø x h, mm	Closure Size, mm	Pack	Price
BOT3542	2103-0002	60	39 x 86	28	12	£15.20
BOT3544	2103-0004	125	50 x 99	38	12	£21.00
BOT3546	2103-0008	250	61 x 131	43	12	£26.90
BOT3548	2103-0016	500	73 x 168	53	12	£43.20
BOT3550	2103-0032	1000	91 x 199	63	6	£29.10