



Thermo Scientific™ Nalgene™ Large Wide Mouth Bottles, HDPE

Thermo
SCIENTIFIC

- Excellent chemical resistance to most acids, bases and alcohols make this a versatile bottle
- Compatible with a wide variety of lab reagents
- Can be used in freezer to –100°C for long term storage
- Linerless cap closes the bottle without a liner that can wrinkle, cause leaks and contaminate your solutions
- Not liquid leak proof

Code	Alt Ref	Capacity, L	Dims, Ø x h, mm	Price
BOT0120	2120-0005	2	119 x 244	£17.10
BOT0122	2120-0010	4	153 x 288	£24.10



Thermo Scientific™ Nalgene™ Large Wide Mouth Bottles, PPCO

Thermo
SCIENTIFIC

- Specifically designed for applications where large volumes require autoclaving
- High quality, laboratory grade plastic materials for reliably low leachables and extractables
- Excellent chemical resistance makes this bottle suitable for a wide variety of uses in the lab or field
- Linerless cap closes the bottle without a liner that can wrinkle, cause leaks and contaminate your solutions
- Withstands temperatures up to 121°C
- Not liquid leak proof

Code	Alt Ref	Capacity, L	Dims, Ø x h, mm	Price
BOT0124	2121-0005	2	120 x 244	£20.00
BOT0126	2121-0010	4	153 x 288	£26.70



Thermo Scientific™ Nalgene™ Heavy Duty Wide Mouth Bottles, HDPE

Thermo
SCIENTIFIC

- HDPE heavy duty bottles with closure
- Especially designed for extra rugged use
- Store and ship liquids or dry chemicals
- Thick walled and wide mouthed
- Can be used as reservoirs and waste traps in automated instrumentation
- White polypropylene linerless closures are leakproof

Code	Alt Ref	Capacity, L	Dims, Ø x h, mm	Closure Size	Pack	Price
BOT0276	2125-1000	1	93 x 226	53B	6	£82.00
BOT0277	2125-2000	2	119 x 260	53B	2	£36.00
BOT0278	2125-4000	4	155 x 338	83B	Each	£42.10