

Thermo Scientific™ Nalgene™ Narrow Mouth Lab Quality Bottles, HDPE Continued



Code	Alt Ref	Capacity, mL	Dims, Ø x h, mm	Closure Size, mm	Pack	Price
BOT0040	2002-9125	4	16 x 41	13	12	£9.60
BOT0042	2002-9025	8	25 x 45	20	12	£9.60
BOT0044	2002-9050	15	25 x 58	20	12	£10.20
BOT0046	2002-0001	30	34 x 61	20	12	£11.40
BOT0048	2002-0002	60	85 x 37	20	12	£13.80
BOT0050	2002-0004	125	50 x 101	24	12	£16.11
BOT0054	2002-0008	250	61 x 133	24	12	£23.70
BOT0056	2002-0016	500	73 x 170	28	12	£34.10
BOT0060	2002-0032	1000	91 x 216	38	6	£28.00

Thermo Scientific™ Nalgene™ Narrow Mouth Economy Bottles, HDPE



- Contain samples and daily working solutions like buffers and bench reagents
- Lightweight versions of the Nalgene™ lab quality bottles
- Suitable for lighter duty, general purpose applications
- Leakproof
- High quality, laboratory grade plastic materials for dependably low leachables and extractables
- Linerless caps provide the ultimate in leakproof protection without the use of liner that can wrinkle, cause leaks or contaminate your reagents



Code	Alt Ref	Capacity, mL	Dims, Ø x h, mm	Closure Size, mm	Pack	Price
BOT5040	2089-0001	30	34 x 61	20	72	£46.90
BOT5042	2089-0002	60	37 x 85	20	72	£57.00
BOT5044	2089-0004	125	50 x 101	24	72	£62.00
BOT5046	2089-0008	250	61 x 133	24	72	£80.00
BOT5048	2089-0016	500	73 x 170	28	48	£79.00
BOT5049	2089-0032	1000	91 x 216	38	24	£63.00

Thermo Scientific™ Nalgene™ Narrow Mouth Bottles, LDPE



- Suitable for applications such as trace metal analysis
- Flexible
- Impact resistant
- Excellent chemical resistance to most acids, bases and alcohols
- Extremely low in trace metal content
- Leakproof
- Provide excellent visibility of contents
- High quality, laboratory grade plastic materials for dependably low leachables and extractables
- Can be used in freezer to -100°C for long term storage
- Translucence allows easy viewing of liquid levels
- Linerless caps provide the ultimate in leakproof protection without use of liner that can wrinkle, cause leaks or contaminate your reagents



Code	Alt Ref	Capacity, mL	Dims, Ø x h, mm	Closure Size, mm	Pack	Price
BOT3220	2003-9025	8	25 x 45	20	12	£9.60
BOT3222	2003-9050	15	25 x 58	20	12	£10.20
BOT3224	2003-0001	30	34 x 61	20	12	£11.40
BOT3226	2003-0002	60	37 x 85	20	12	£13.80
BOT3228	2003-0004	125	50 x 101	24	12	£16.20
BOT3230	2003-0008	250	61 x 133	24	12	£23.10
BOT3232	2003-0016	500	73 x 170	28	12	£34.10
BOT3234	2003-0032	1000	91 x 216	38	6	£28.00