

Azlon® Narrow Mouth Bottles, PP Continued



Code	Alt Ref	Capacity, mL	Dims, Ø x h, mm	Pack	Price
BOT3060	BNP0030PN	30	32 x 70	10	£7.30
BOT3064	BNP0125PN	125	51 x 106	10	£12.10
BOT3066	BNP0150PN	150	51 x 118	10	£13.20
BOT3068	BNP0250P	250	60 x 142	10	£19.70
BOT3070	BNP0500P	500	75 x 163	10	£25.90
BOT3072	BNP1000P	1000	90 x 203	5	£22.50
BOT3074	BNP2000P	2000	120 x 234	5	£45.90



Thermo Scientific™ Nalgene™ Rectangular Bottles, HDPE



- Transport, sample or store liquids and dry materials
- Space efficient rectangular shape takes up less shelf space than traditional round bottles
- Excellent chemical resistance to most acids, bases and alcohols
- Sturdy wall construction with straight sides and a wide mouth for easy filling
- High quality, laboratory grade, plastic materials for dependably low leachables and extractables
- Suitable for freezer for storage up to -100°C
- Linerless polypropylene screw closures provide the ultimate in leakproof protection without the use of a liner that can wrinkle, cause leaks, or contaminate your reagents
- Leakproof

Code	Alt Ref	Capacity, mL	Closure Size, mm	Dims, l x w x h, mm	Pack	Price
BOT0150	2007-0004	125	28	61 x 38 x 102	12	£22.80
BOT0152	2007-0008	250	38	76 x 51 x 117	12	£33.30
BOT0154	2007-0016	500	48	94 x 60 x 147	12	£53.00
BOT0156	2007-0032	1000	53	125 x 71 x 53	6	£61.00
BOT0158	2007-0064	2000	63	152 x 84 x 242	6	£67.00



Thermo Scientific™ Nalgene™ Square Bottles, PETG



- Ideal for storage and sampling of active pharmaceutical ingredients and bulk intermediates, and also well suited to the preparation and containment of buffers, cultures or prolonged storage of pH sensitive liquids such as culture media
- Durable, break resistant polyethylene terephthalate copolyester, glycol modified (PETG)
- Suitable for frozen storage down to -40°C
- Certified sterile and non pyrogenic
- Heat shrink band around closure and neck for tamper evident seal
- Leakproof bottles and HDPE closures provided
- Supported by extensive validation documentation
- 2L size has molded in handgrips

Code	Alt Ref	Capacity, mL	Dims, w x h, mm	Closure Size, mm	Pack	Price
BOT5330	2019-0030	30	38 x 64	20	96	£215.00
BOT5332	2019-0060	60	41 x 82	24	96	£279.00
BOT5334	2019-0125	125	54 x 110	38	48	£227.00
BOT5336	2019-0250	250	61 x 146	38	48	£274.00
BOT5338	2019-0500	500	74 x 177	38	24	£200.00
BOT5340	2019-1000	1000	94 x 220	38	24	£270.00
BOT5342	2019-2000	2000	116 x 271	53B	12	£224.00