



Corning® Octagonal Bottles, PET



- PET resin makes the bottle lightweight, break resistant and as transparent as glass
- Ideal for storage and sampling of media, biologicals, and other aqueous solutions
- Packaged in two convenient, shrink wrapped trays within an outer bag to assure cleanliness
- Break resistant bottles are an ergonomic and easy to handle alternative to glass
- Bottles are sterilised to 10⁻⁶ SAL and non pyrogenic
- Leakproof, tamper evident HDPE screw cap
- Molded graduations for accurate measurements
- USP class VI compliant bottles are validated non-cytotoxic and non-hemolytic

Code	Capacity, mL	Dims, w x h, mm	Graduations, mL	Pack	Price
431729	30	38 x 61	10	120	£251.00
431730	60	40 x 81	10	96	£263.00
431731	125	53 x 115	25	48	£206.00
431732	250	64 x 142	50	48	£257.00
431733	500	74.2 x 180	50	24	£199.00
431734	1000	96 x 220	50	24	£255.00



Thermo Scientific™ Nalgene™ Heavy Duty Vacuum Bottles, PPCO



- Hold full vacuum for up to 24 hours
- Ideal for heavy-duty biotechnology lab use
- Thick walled PPCO bottles can withstand repeated room temperature vacuum applications
- Use as a waste aspirator bottle or as autoclavable vessel for scaleup activities in cell culture applications
- Large size with wide mouth
- Chemical resistant
- Rugged 53B and 83B caps with TPE gaskets

Code	Alt Ref	Capacity, L	Dims, Ø x h, mm	Closure Size	Pack	Price
BOT0322	2126-1000	1	93 x 226	53B	6	£188.00
BOT0324	2126-2000	2	119 x 260	53B	2	£76.00
BOT0326	2126-4000	4	155 x 338	83B	Each	£62.00
BOT0328	2126-5000	5	156 x 414	83B	Each	£94.00



SLS Select Carboys



- Autoclavable polypropylene
- Wide shoulder handle for easy carrying and pouring
- Ideal for storage and transporting large batches of culture media, distilled water and other solutions

Code	Capacity, L	Colour	Dims, w x h, mm	Price
ASP1200	10	Translucent	249 x 390	£46.10
ASP1204	10	Amber	249 x 390	£50.10
ASP1202	20	Translucent	285 x 535	£67.00