

## Azlon® Integral Wash Bottles



- Integral spout to minimise the risk of contamination
- Wide neck enables ease of filling
- Colour coded vented cap with PTFE membrane to prevent pressure build up and dripping when using solvents
- Unprinted to allow custom labelling
- LDPE

Code	Alt Ref	Capacity, mL	Cap Colours	Pack	Price
<b>BOT0034</b>	WGW800P	500	Red, yellow, blue, green, orange	5	£21.70



## Azlon® Integral Multilingual Wash Bottles



- Integral spout to minimise contamination
- Indelibly printed in English, French, German and Spanish
- Choose from five solvent labels all including Globally Harmonised Symbols, NFPA diamond, risk phrases and safety information
- Colour coded vented cap with PTFE membrane to prevent pressure build up and dripping when used with solvents
- LDPE
- Capacity: 500mL

Supplied in packs of 5

Code	Alt Ref	Label	Cap Colour	Price
<b>BOT0036</b>	WGW801P	Distilled water	White	£25.40
<b>BOT0038</b>	WGW802P	Acetone	Red	£28.50
<b>BOT0302</b>	WGW803P	Methanol	Green	£28.50
<b>BOT0304</b>	WGW804P	Isopropanol	Blue	£28.50
<b>BOT0306</b>	WGW805P	Ethanol	Orange	£28.50



## Azlon® Narrow Neck Wash Bottles



- Narrow neck
- Traditional and popular square shoulder
- One piece, white PP cap
- Closure features a precision jet tip aperture for fine control
- LDPE

Supplied in packs of 5

Code	Alt Ref	Capacity, mL	Price
<b>BOT0210</b>	WGF230P	250	£13.40
<b>BOT0220</b>	WGF236P	500	£16.10

