


Pencil Form/Artist Brushes


Pure soft hair brushes in metal mount on a wooden handle. Overall length: 170mm.

Code	Head Width, mm	Size	Price
BRU2050	1	00	£4.10
BRU2052	1.5	0	£4.80
BRU2054	2.5	3	£5.43
BRU2056	3.5	4	£5.51
BRU2058	4	5	£7.66



Sieve Brushes

Code	Description	Price
SIE1110	Double ended nylon	£16.04
SIE1111	Small brass/nylon	£28.00
SIE1113	Large bristle/hair	£18.00



Specimen Buckets

- Complete with airtight lids
- White polypropylene

Code	Capacity, L	Description	Pack	Price
BUC2002	1	Without handle	10	£25.00
BUC2012	5	With handle	10	£60.00
BUC2016	10	With handle	Each	£14.46
BUC2020	25	With handle	Each	£18.60
BUC2024	60	Without handle	Each	£27.80



Sampling Systems Buckets



- Swing handle for easy lifting
- 10, 12 and 15L buckets have a skirt band around the base for stability and extra protection to the base
- Very durable
- Graduated
- Crevice free interior is easy to clean after use
- 304 stainless steel and finished to a quality satin finish



Code	Alt Ref	Capacity, L	Dims, Ø x h, mm	Lid	Price
SAM1336	8095A-250	6	257 x 240	No	£96.00
SAM1337	8095A-250S	10	300 x 258	No	£88.00
SAM1340	A264-10L	10	300 x 295	Yes	£30.00
SAM1338	A264-12	12	311 x 290	No	£96.00
SAM1341	A264-12L	12	330 x 315	Yes	£31.00
SAM1339	A264-15	15	323 x 330	No	£100.00
SAM1342	A264-15L	15	350 x 360	Yes	£38.00