



Rodwell Crystal Autoclave



The Crystal is a powerful workhorse, combining large capacity with clever design features to save you time and money in the running of your laboratory.

- MP25 controller
- Square section chamber
- Cost-effective front loading hinge / bolt door
- Fixed seal
- Space-saving design
- Integral drip tray under door
- Efficient fan assisted cooling
- Loading trolley
- Scale prevention systems
- Waste discard system
- Range of cycle recording systems
- 3ph/N 50/60Hz

Code	Alt Ref	Capacity, L	Chamber Dims, w x d x h, mm	Dims, w x d x h, mm	Price
AUT0503	CRYS3	300	510 x 1030 x 510	830 x 1450 x 1800	£27949.00



Rodwell Gemini Autoclave



The Gemini uniquely combines two popular autoclave vessels into one robust frame, ideal for any busy laboratory.

- 2 x 158L chambers in one frame (316L total capacity)
- MP25 controller
- Independently controlled chambers
- Robust space saving design
- Integral drip tray
- Grade 316 stainless steel cylindrical chambers
- Anti door bolt over tight mechanism
- Complete versatility for almost any application
- Enables separation of clean and dirty loads

Code	Alt Ref	Capacity, L	Chamber Dims, Ø x d, mm	Dims, w x d x h, mm	Price
AUT1293	GEM21	2 x 158	460 x 960	650 x 1310 x 1920	£32500.00



Rodwell Monarch Autoclaves



The Monarch 50 is ideal as an under-bench solution, and the Monarch 75 is perfect for large vessels such as fermenters.

- MP25 controller
- Top loading, cylindrical chamber
- Angled control panel
- Low loading height
- Load capacity for up to 3 baskets
- Counterbalanced lid for easy opening
- Monarch 50 fits underneath a standard bench
- Power supply: single phase 3kW "plug & play"

Code	Alt Ref	Capacity, L	Chamber Dims, Ø x d, mm	Dims, w x d x h, mm	Price
AUT1300	MON7	50	350 x 520	470 x 690 x 840	£9900.00
AUT1302	MON8	75	350 x 780	470 x 690 x 1110	£10725.00

Options and Accessories

Code	Alt Ref	Description	Price
AUT2008	AC28	Waste discard system	£709.00
AUT1357	AC4	Fan assisted cooling	£464.00
AUT2019	AC164	Condense collector	£277.00
AUT2017	AC70	Additional loading basket	£273.00